SERVOTOME® CLASSIC

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ENGLISH

C O N T E N T

I - PRESENTATION AND INSTALLATION	
Presentation	6
Technical description	6-7
Installation	7-8
II - RECOMMENDATIONS AND CONTRA INDICATIONS	
Recommendations	8-9
Contra indications	9
III - MAINTENANCE AND STERILISATION	
Maintenance of the housing	10
Maintenance and sterilisation of the electrodes, electrode holder and cord	10
Troubleshooting	11
IV - LIABILITY	12
V - NOTES	12
VI - SYMBOLS	13
POWER SETTING CHART	14-15

I - PRESENTATION AND INSTALLATION

1.1 PRESENTATION

SERVOTOME CLASSIC is a generator of high frequency current, which with a large range of electrodes can perform a high quality cutting together with coagulation control during the surgical operation.

SERVOTOME CLASSIC is delivered with (Fig. 1) :

- 1 control pedal and cord (2 m long) (Fig. 1-1).
- 1 disconnectable earthed power cord (2 m long) (Fig. 1-2).
- 1 electrode holder (Fig. 1-3) and cord (Fig. 1-4).
- 1 neutral plate and cord (Fig. 1-5)
- 1 set of 10 electrodes (Fig. 1-6)

1.2 TECHNICAL DESCRIPTION

a) Front panel (Fig. 2)

- On / Off switch (**I/O**) (Fig. 2-1)
- "Switched On" light (Fig. 2-2).
- Power adjustment of the depth of the cut (Fig 2-3).
- Orange operating light (Fig. 2-4).
- Cutting current adjustment (depth of coagulation) (Fig. 2-5)
- Plug for the electrode holder cord (Fig. 2-6).

b) Rear panel (Fig. 3)

- Neutral plate socket (Fig. 3-1)
- Pedal socket (Fig. 3-2)
- Power socket (Fig. 3-3).
- Fuse housing (Fig. 3-4).

c) **Technical data** 

Power supply : 230 V~ 50/60 Hz (110 V~ optional).

Class I equipment, type BF IP 21

Intermittent operation : 10s / 30s

Power consumption : 220 VA max.

Output power (Fig. 4) : 50 W / 700 Ω

Output frequency : 1,8 MHz \pm 0,7

Overall dimensions:

Depth : 170 mm

Width : 140 mm,

Height : 65 mm

Weight : 1.5 Kg

Manufactured in compliance with the regulations and standards in force (Directive 93/42/EEC).

Note : Upon specific request any technical information required to repair elements of the unit identified as repairable by SATELEC is available to the technical staff of SATELEC approved dealers.

1.3 INSTALLATION

a) Recommendations

All electrical connections to the mains to be performed by an approved dental technician

SERVOTOME CLASSIC to only be used by qualified and trained practitioners (dental surgeon)

b) Connection

Warning : All accessories (Electrode holder, cord...) must be cleaned, decontaminated and sterilised before being used for the first time. (See chapter III - Maintenance and Sterilisation).

SERVOTOME CLASSIC front panel

- Plug the electrode holder cord into the socket (Fig. 2-6)

SERVOTOME CLASSIC back panel

Plug:

- The neutral plate cord into its socket (Fig. 3-1).
- The pedal cord DIN pin into its socket (Fig. 3-2).
- Plug the mains cord into a wall socket (Fig. 3-3).

c) Starting the unit

- Plug the unit into the mains, ensuring that the wall socket is equipped with an earth connection.
- Place the neutral plate on the seat level with the patient's shoulder-blade.

Respect the indication "patient side" as printed on the neutral plate.

- Turn the On/Off (**I/O**) switch to "I". The "switched on" green indicator lights (Fig. 2-2).

- Adjust the power (Fig. 2-3) and the coagulation (Fig. 2-5). (See Fig. 5).
- Press the pedal. The cutting effect is available straight away. The orange indicator lights (Fig. 2-3) and a sound signal is activated. This stops as soon as the pressure on the pedal is released.

Note : The unit must be left to rest for a time equal to at least 3 times the utilisation time.

II - RECOMMENDATIONS AND CONTRA INDICATIONS

2.1 RECOMMENDATIONS

In order to lower the risk of accidental burns, care must be taken.

- The entire surface of the neutral plate must be in contact with the clothed patient's body.
- The patient must not be in contact with any metal parts linked to earth (for example: the operating table, holder etc....).
- Ensure that the unit is well connected to the mains, well insulated. The use of an earth socket is compulsory.
- Skin to skin contact is inadvisable (practitioner/patient).
- When SERVOTOME CLASSIC is used simultaneously with physiological monitoring equipment on a patient, if they are not equipped with either smoothing resistances or self HF, their sensors must be placed as far away as possible from the surgical electrodes. Monitoring needle electrodes are not recommended.

- Surgical electrode holder cords must be placed away from the patient or any conductive materials.
 - The output power must be as low as possible for the purpose required.
 - Some materials such as cotton, wool and surgical gauze saturated with oxygen can catch fire when in contact with sparkles generated during the normal utilisation of SERVOTOME CLASSIC.
 - Remove active electrodes that are not being used away from the patient.
 - Accessories (electrode holder, electrodes, neutral plate) must be replaced if damaged or worn.
 - SERVOTOME CLASSIC housing together with the accessories do not require any specific maintenance other than normal maintenance.
 - The electrode holder, its cord and the electrodes must be cleaned, decontaminated and sterilised before being used.
 - The user must check the accessories and cords on a regular basis in order to detect any insulation problem and to have them replaced if necessary.
 - A very low power can originate from a poor connection of the neutral plate.
 - Because of the "high frequency" radiations, other electrical equipment may be disturbed when SERVOTOME CLASSIC is operated.
- Note :** If there is a fault, it is recommended to contact your SERVOTOME CLASSIC local supplier rather than any other supplier.

2. 2 CONTRA INDICATIONS

- SERVOTOME CLASSIC must not be used on patients or by practitioners wearing a pacemaker.
- SERVOTOME CLASSIC must not be used in explosion areas or close to gaseous anaesthetic.
- SERVOTOME CLASSIC must only be used with accessories supplied by SATELEC.

Other misuses listed below are to inform the operator about electrosurgery limitations.

- Agitated, emotive or extremely nervous patients.
- Incomplete anaesthetic.
- Delicate surgery (mucoperiosteal surgery, grafts, etc...).
- Extremely fragile tissues.
- Lack of knowledge or training in electrosurgery theories.
- Lack of knowledge of the patient and his/her health .

Warning : The electrode must be pushed firmly in the electrode holder.

III - MAINTENANCE AND STERILISATION

3.1 MAINTENANCE OF THE HOUSING

The SERVOTOME CLASSIC housing can be cleaned with alcohol or disinfectant wipes for dental use.

3.2 MAINTENANCE AND STERILISATION OF THE ELECTRODES, ELECTRODE HOLDER AND CORD

a) Cleaning

Electrodes :

They can be cleaned with alcohol or with disinfectant wipes for dental use.

Electrode holder and cord :

They can be cleaned with alcohol or with disinfectant wipes for dental use. The screw thread can be cleaned with a piece of cotton saturated with alcohol.

Do not clean the electrode holder in a ultrasonic tank.

b) Decontamination

Decontaminate the electrodes by immersing them into an alkaline solution such as glutaraldehyde or quaternary ammonium for approximately 20 min, depending on concentration. (Follow supplier's recommendations).

c) Sterilisation

It is absolutely necessary to sterilise the electrodes every time before use.

Sterilisation conditions and cycles according to equipment:

Equipment: autoclave

Wet heat sterilisation:

134° C under 2 bar (200 KPa) for 20 min

Important :

Different metallic parts should not be in contact during sterilisation. Any contact would lead to the creation of electrolytic effects that would result in local damage. In order to avoid this phenomenon, place every instrument in a cloth or a sterilisation bag.

3.3 TROUBLESHOOTING

FAULT	POSSIBLE CAUSE	REMEDY
No operation	<ul style="list-style-type: none"> • Faulty power cord connection • Fuse(s) blown • Faulty pedal cord connection 	<ul style="list-style-type: none"> • Check mains plug • Replace fuse(s) • Check the pedal : Plug incorrectly inserted
The electrode does not cut well or does not cut at all	<ul style="list-style-type: none"> • Neutral plate incorrectly positioned • Dirty electrode • Electrode movement too rapid 	<ul style="list-style-type: none"> • Respect the indication "patient side" as printed on the neutral plate. Check that the entire plate is in contact with the patient (wearing light clothes) and at the level of his shoulder-blades for instance. • Turn the On/Off switch to "O" and clean the electrode (see chapter 3) • Slow the operative act down.
Electrode bonding to the tissues	<ul style="list-style-type: none"> • Power adjustment too low 	<ul style="list-style-type: none"> • Increase power until cutting mode is reached. It is not necessary to further increase.
The electrode sparkling while cutting	<ul style="list-style-type: none"> • Power adjustment too high 	<ul style="list-style-type: none"> • Decrease power until cutting recommended setting is reached. It is not necessary to increase.

IV - LIABILITY

V - NOTES

The manufacturer's liability shall not be pledged in the following cases:

- Failure to respect the manufacturer's instructions for installation (Network voltage).
- Maintenance or repairs by persons not authorised by the manufacturer.
- Use on an electrical installation that does not comply with the regulations in force
- Other uses as those specified in this manual.
- The manufacturer reserves the right to modify the device and/or the instruction manual without prior notice.

SERVOTOME CLASSIC must not be used on patients or by practitioners wearing a pacemaker

VI - SYMBOLS

5.1 RATING PLATE

 Alternating current


 BF type unit

 Caution, refer to accompanying documents


5.2 UNIT SYMBOLS

 Power level adjustment

 Coagulation level adjustment

 Neutral plate linked to earth for high frequency currents

 Neutral plate connector

 Pedal connector

CE
0459 EEC labelling conforming to (93/42/EEC) directive

SERVOTOME Classic System
Appareil d'Electrochirurgie H.F
High Frequency Surgical Equipment

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Appareil de classe : I  
(Class equipment)

Service intermittent : 10s / 30s
(Intermittent Operation)

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
Made in France by S A T E L E C
Z.I du Phare 33708 MERIGNAC

230V : 2 x 1,25 AT

115V : 2 x 2 AT

Primary fuse values according to the mains.

Procedure	Section	Color	Power	Coagulation	Electrodes
Drainage of abscess	Electrosection	Blue	5/6	6/7	
		Yellow	6/7	0	
Frenectomy	Electrosection Electrosection with coagulation	Blue	6/7	4/5	
	Electrosection Electrosection with coagulation	Yellow	6/7	0	
Impacted teeth	Electrosection	Blue	6	7	
		White	7/8	4	
Gingivectomy plastic surgery	Electrosection	Blue	5/6	6/7	
		White	7/8	4	
		Yellow	6/7	0	
Stripping	Electrosection	Blue	5/6	6/7	
		White	7/8	4	
Elongations coronaires	Electrosection Electrosection with coagulation	Blue	5/6	6/7	
	Electrosection Electrosection with coagulation	Yellow	6/7	0	
Disengagement of necks	Slightly coagulating electrosection	Blue	5	6	

Procedure	Section	Color	Power	Coagulation	Electrodes
Marginal gingiva, vestibular or lingual	Electrosection	Blue	5/6	6/7	
		Yellow	6/7	0	
Gingival evictions (interproximal gingiva)	Electrosection	Yellow	6/7	0	
	Slightly coagulating electrosection				
Fractured tooth	Electrosection	Blue	5/6	6/7	
	Slightly coagulating electrosection				
	Electrosection	White	7/8	4	
	Slightly coagulating electrosection	Yellow	6/7	0	
Insertion of dental dam	Electrosection	Blue	5	6	
		Yellow	6/7	0	
Hyperplasia, Mucous, Proliferation, Biopsy, Benign tumors	Electrosection + Coagulation	Blue	5/6	6/7	
		Yellow	6/7	0	
Coagulation	Coagulation	Red	5/6	6/7	

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 50/60Hz
(Mains)

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Appareil de classe : I  
(*Class equipment*)

Service intermittent : 10s / 30s
(*Intermittent Operation*)

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










Procedure	Section	Color	Power	Coagulation	Electrodes
Drainage of abscess	Electrosection	Blue	5/6	6/7	 
		Yellow	6/7	0	 
Frenectomy	Electrosection Electrosection with coagulation	Blue	6/7	4/5	 
	Electrosection Electrosection with coagulation	Yellow	6/7	0	 
Impacted teeth	Electrosection	Blue	6	7	 
		White	7/8	4	 
Gingivectomy plastic surgery	Electrosection	Blue	5/6	6/7	 
		White	7/8	4	 
		Yellow	6/7	0	 
Stripping	Electrosection	Blue	5/6	6/7	 
		White	7/8	4	 
Elongations coronaires	Electrosection Electrosection with coagulation	Blue	5/6	6/7	 
	Electrosection Electrosection with coagulation	Yellow	6/7	0	 
Disengagement of necks	Slightly coagulating electrosection	Blue	5	6	 

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