# IPS e.max® CAD

**[en] Instructions for Use**

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## 1 Intended use

### Intended purpose

Anterior and posterior single-tooth restorations, restoration of teeth with 3-unit bridges up to the second premolar as the terminal abutment

### Patient target groups

Patients with permanent teeth

### Intended users / Special training

- Dentists (chairside fabrication of restorations, clinical procedure)
  - Dental laboratory technicians (fabrication of restorations in the laboratory)
- No special training required.

### Use

For dental use only.

### Description

IPS e.max CAD is a tried-and-tested lithium disilicate glass-ceramic block (LS<sub>2</sub>) for the fabrication of fixed anterior and posterior restorations.

IPS e.max CAD can be processed in an authorized CAD/CAM machine<sup>[1]</sup> in the intermediate crystalline state ( $\geq 130$  MPa). After wet-processing of the block, the restoration is crystallized in a ceramic furnace.<sup>[2]</sup>

### Technical data

| Property                                          | Specification     | Typical mean value |
|---------------------------------------------------|-------------------|--------------------|
| CTE (25–500 °C) [ $10^{-6}$ K <sup>-1</sup> ]     | 10.1 ± 0.5        | –                  |
| Flexural strength (biaxial) [MPa]                 | $\geq 360$        | 530 <sup>[1]</sup> |
| Chemical solubility [ $\mu\text{g}/\text{cm}^2$ ] | < 100             | –                  |
| Type/Class                                        | Type II / Class 3 | –                  |

According to ISO 6872:2015

<sup>[1]</sup> mean biaxial flexural strength evaluated over 10 years of quality measurements

### Indications

- Missing tooth structure in anterior and posterior teeth
- Partial edentulism in the anterior and posterior region

#### Types of restorations:

- Veneers
- Inlays
- Onlays (e.g. occlusal veneers, partial crowns)
- Crowns
- Three-unit bridges up to the second premolar as the terminal abutment

### Contraindications

- Patients with substantially reduced residual dentition
- The use of the product is contraindicated if the patient is known to be allergic to any of its ingredients.

### Limitations of use

- Inlay, cantilever and Maryland bridges
- Hybrid abutments and hybrid abutment crowns (use IPS e.max CAD Abutment Solutions for these types of restorations)
- Width of the pontic: anterior region > 11 mm, premolar region > 9 mm
- Temporary cementation
- Complete veneering of molar crowns
- Very deep sub-gingival preparations
- Untreated bruxism (the use of a splint is indicated after incorporation)
-  Do not re-use

Do not re-use

Additional limitations of use for minimally invasive crowns:

- Layer thicknesses below 1 mm
- Preparations with sharp edges
- Preparations that are not anatomically supported and feature varying layer thicknesses
- Conventional and self-adhesive cementation
- Build-up materials other than composite resin
- Absence of canine guidance
- Crowns on implants

### Processing restrictions

In the following situations, a successful procedure cannot be ensured:

- Falling short of the required minimum thicknesses
- Milling of the blocks in a non-compatible CAD/CAM system
- Crystallization in a non-approved and non-calibrated ceramic furnace
- Mixing of IPS e.max CAD Crystall./Glaze, Shades and Stains with other dental ceramics (e.g. IPS Ivocolor® Glaze, Shades and Essences).
- Layering with a veneering ceramic other than IPS e.max Ceram

IPS e.max CAD Crystall./Glaze Spray:

- Intraoral application
- Application on IPS e.max CAD restorations, if they are to be veneered with IPS e.max Ceram
- Restorations whose inner surfaces cannot be reliably and precisely covered with IPS Objekt Fix Putty/Flow (e.g. inlays)

### System requirements

IPS e.max CAD must be processed with an authorized CAD/CAM system.<sup>[1]</sup>

### Side effects

There are no known side effects to date.

### Interactions

There are no known interactions to date.

### Clinical benefit

- Reconstruction of chewing function
- Restoration of esthetics

### Composition

#### Lithium disilicate glass-ceramic

After the manufacturing process of the glass-ceramic, a stable and inert network is formed in which the different elements are incorporated via oxygen bridges. The composition is determined as oxides.

| Oxide                                                                                                                                                                                      | in wt. %    |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------|
| SiO <sub>2</sub>                                                                                                                                                                           | 57.0 – 80.0 |
| Li <sub>2</sub> O                                                                                                                                                                          | 11.0 – 19.0 |
| K <sub>2</sub> O                                                                                                                                                                           | 0.0 – 13.0  |
| P <sub>2</sub> O <sub>5</sub>                                                                                                                                                              | 0.0 – 11.0  |
| ZrO <sub>2</sub>                                                                                                                                                                           | 0.0 – 8.0   |
| ZnO                                                                                                                                                                                        | 0.0 – 8.0   |
| Al <sub>2</sub> O <sub>3</sub>                                                                                                                                                             | 0.0 – 5.0   |
| MgO                                                                                                                                                                                        | 0.0 – 5.0   |
| Colouring oxides (CeO <sub>2</sub> , MnO <sub>2</sub> , V <sub>2</sub> O <sub>5</sub> , Tb <sub>4</sub> O <sub>7</sub> , Er <sub>2</sub> O <sub>3</sub> , La <sub>2</sub> O <sub>3</sub> ) | 0.0 – 8.0   |

## 2 Application

### Processing techniques and compatible materials

#### – On the blue restoration

- *Polishing technique (self glaze):* e.g. with OptraGloss® or dental lab polishers, followed by crystallization firing without individual characterization and glazing
- *Staining technique*
  - a) Glazing with IPS e.max CAD Crystall./Glaze Spray, followed by speed crystallization (speed crystallization and glaze firing in one step)
  - b) Staining and glazing with IPS e.max CAD Crystall./Glaze Paste/Fluo, followed by crystallization and stains/glaze firing in one step
  - c) Staining and glazing with IPS e.max CAD Crystall./Glaze Spray, followed by crystallization and stains/glaze firing in one step

#### – On the tooth-coloured restoration

- a) *Staining technique:* Characterization/glaze firing of the tooth-coloured restorations with either IPS e.max CAD Crystall./ or IPS Ivocolor materials
- b) *Cut-back technique:* Incisal firing using IPS e.max Ceram materials. Stain/glaze firing using IPS Ivocolor materials
- c) *Layering technique:* Dentin/incisal firing using IPS e.max Ceram materials. Stain/glaze firing using IPS Ivocolor materials

See Instructions

The respective instructions for use must be observed.

## Block concept

| Degree of translucency              | Block sizes/<br>shades                                              | Processing technique                                                                                        |                                                                                                            |                                                                                                            |                                                                                                            | Types of restorations             |                            |                  |                  |                  |                  |                  |
|-------------------------------------|---------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------|-----------------------------------|----------------------------|------------------|------------------|------------------|------------------|------------------|
|                                     |                                                                     | Polishing<br>technique<br> | Staining<br>technique<br> | Cut-back<br>technique<br> | Layering<br>technique<br> | Occlusal veneer<br><sup>[a]</sup> | Thin veneer <sup>[a]</sup> | Veneer           | Inlay, onlay     | Partial crown    | Crown            | 3-unit bridge    |
| <b>HT</b><br>High<br>Translucency   | I12, C14,<br>B40, B40 L<br>available in<br>20 shades <sup>[f]</sup> | ✓                                                                                                           | ✓                                                                                                          | ✓                                                                                                          |                                                                                                            | ✓                                 | ✓ <sup>[d]</sup>           | ✓ <sup>[d]</sup> | ✓ <sup>[d]</sup> | ✓ <sup>[d]</sup> | ✓ <sup>[d]</sup> |                  |
| <b>MT</b><br>Medium<br>Translucency | C14,<br>available in<br>7 shades                                    | ✓                                                                                                           | ✓                                                                                                          | ✓                                                                                                          |                                                                                                            | ✓                                 | ✓ <sup>[d]</sup>           | ✓ <sup>[d]</sup> | ✓ <sup>[d]</sup> | ✓ <sup>[d]</sup> | ✓ <sup>[d]</sup> |                  |
| <b>LT</b><br>Low<br>Translucency    | I12, C14,<br>C16, B32<br>in 20<br>shades <sup>[f]</sup>             | ✓                                                                                                           | ✓                                                                                                          | ✓                                                                                                          |                                                                                                            |                                   |                            | ✓ <sup>[d]</sup> |                  | ✓ <sup>[d]</sup> | ✓ <sup>[d]</sup> | ✓ <sup>[b]</sup> |
| <b>M0</b><br>Medium<br>Opacity      | C14,<br>in M0<br>0 – 4 <sup>[f]</sup>                               |                                                                                                             |                                                                                                            |                                                                                                            | ✓                                                                                                          |                                   |                            |                  |                  |                  | ✓ <sup>[c]</sup> |                  |
| <b>I</b><br>Impulse                 | C14<br>in 01 and<br>02                                              | ✓                                                                                                           | ✓                                                                                                          | ✓                                                                                                          |                                                                                                            | ✓                                 | ✓                          | ✓                |                  |                  |                  |                  |

<sup>[a]</sup> The cut-back technique must not be used when fabricating thin anterior veneers and occlusal veneers.

<sup>[b]</sup> Only up to the second premolar as the terminal abutment

<sup>[c]</sup> Up to the second premolar

<sup>[d]</sup> Max. 2 units if the IPS Speed Tray is used

<sup>[f]</sup> The shade range may vary depending on the translucency level/block size or the CAD/CAM machine.

## Compatible cementation methods

|                                                                          | Adhesive cementation<br>e.g. Variolink® Esthetic,<br>Multilink® Automix | Self-adhesive cementation<br>e.g. SpeedCEM® Plus | Conventional cementation<br>e.g. Vivaglass® CEM,<br>ZirCAD® Cement |
|--------------------------------------------------------------------------|-------------------------------------------------------------------------|--------------------------------------------------|--------------------------------------------------------------------|
| Preparation requirements                                                 | Non-retentive preparation                                               | Retentive preparation                            | Retentive preparation                                              |
| Veneers                                                                  | ✓                                                                       | –                                                | –                                                                  |
| Inlays, onlays (e.g. occlusal veneers,<br>partial crowns)                | ✓                                                                       | –                                                | –                                                                  |
| Minimally invasive crowns                                                | ✓                                                                       | –                                                | –                                                                  |
| Crowns                                                                   | ✓                                                                       | ✓                                                | ✓                                                                  |
| Three-unit bridges up to the second<br>premolar as the terminal abutment | ✓                                                                       | ✓                                                | ✓                                                                  |

For further information, please refer to the flow chart "Pretreatment and cementation".











## Shade selection

Clean the teeth prior to shade determination. Determine the shade based on the moist and unprepared tooth or the neighbouring teeth. Determine the shade of the prepared tooth based on the largest discoloured area of the preparation. The IPS e.max Shade Navigation App helps to select the most suitable block.

## Preparation

Prepare the tooth based on the guidelines for all-ceramic restorations and make sure to observe the minimum layer thicknesses:

- No angles or sharp edges
- Shoulder preparation with a rounded internal angle and/or pronounced chamfer
- The indicated dimensions reflect the minimum thickness for IPS e.max CAD restorations.
- The thickness of the incisal edge of the preparation, particularly in anterior teeth, must be at least 1.0 mm in order to ensure optimum processing in the CAD/CAM machine.

| Minimally invasive posterior crown (adhesive cementation mandatory)               | Minimally invasive anterior crown (adhesive cementation mandatory)                | Inlay                                                                             | Onlay                                                                             | Thin veneer                                                                        |
|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|------------------------------------------------------------------------------------|
|  |  |  |  |  |
| Veneer                                                                            | Posterior crown/bridge abutment in the premolar region                            | Anterior crown/bridge abutment in the anterior region                             | Occlusal veneer (table top)                                                       | Partial crown                                                                      |
|  |  |  |  |  |

Dimensions in mm

## Minimum layer thickness of restorations when using the staining technique

| Cementation                                                 | Mandatory adhesive cementation |                              |                                             |                                                               | Optional adhesive, self-adhesive or conventional cementation |                  |                                                  |                  |
|-------------------------------------------------------------|--------------------------------|------------------------------|---------------------------------------------|---------------------------------------------------------------|--------------------------------------------------------------|------------------|--------------------------------------------------|------------------|
|                                                             | Thin veneer                    | Inlay                        | Onlay (e.g. occlusal veneer, partial crown) | Minimally invasive crown in the anterior and posterior region | Crown                                                        |                  | Bridge                                           |                  |
|                                                             |                                |                              |                                             |                                                               | Anterior region                                              | Posterior region | Anterior region                                  | Posterior region |
| Types of restorations                                       |                                |                              |                                             |                                                               |                                                              |                  |                                                  |                  |
| Minimum layer thickness IPS e.max CAD – Polishing technique |                                |                              |                                             |                                                               |                                                              |                  |                                                  |                  |
| Minimum layer thickness IPS e.max CAD – Staining technique  |                                |                              |                                             |                                                               |                                                              |                  |                                                  |                  |
| Incisal/occlusal                                            | 0.5                            | 1.0<br>Depth of the fissures | 1.0                                         | 1.0                                                           | 1.5                                                          | 1.5              | 1.5                                              | 1.5              |
| Circular                                                    | 0.4                            | 1.0<br>Isthmus width         | 1.0                                         | 1.0                                                           | 1.2                                                          | 1.5              | 1.2                                              | 1.5              |
| Connector dimension                                         | –                              | –                            | –                                           | –                                                             | –                                                            | –                | 16 mm <sup>2</sup><br>In general: height ≥ width |                  |

The design of bridge connectors should be extended in a vertical direction rather than in a horizontal direction.

## Cut-back and layering technique

When using the cut-back or layering technique, a reduced, tooth shape-supporting framework is built up to full contour using the IPS e.max Ceram layering materials.

In the fabrication of veneered or partially veneered restorations, the available space in large preparations must be filled by appropriately dimensioning the high-strength IPS e.max CAD component and not by adding additional IPS e.max Ceram layering material (see below).

| Cementation           | Mandatory adhesive cementation |       |       |               | Optional adhesive, self-adhesive or conventional cementation |                 |              |                 |                 |
|-----------------------|--------------------------------|-------|-------|---------------|--------------------------------------------------------------|-----------------|--------------|-----------------|-----------------|
| Types of restorations | Veneer                         | Inlay | Onlay | Partial crown | Crown                                                        |                 |              | Bridge          |                 |
|                       |                                |       |       |               | Anterior region                                              | Premolar region | Molar region | Anterior region | Premolar region |

### Minimum layer thickness IPS e.max CAD – Cut-back technique

|                     |     |   |   |     |     |     |     |                                                  |     |
|---------------------|-----|---|---|-----|-----|-----|-----|--------------------------------------------------|-----|
| Incisal/occlusal    | 0.4 | – | – | 1.3 | 0.8 | 1.0 | 1.3 | 0.8                                              | 1.0 |
| Circular            | 0.6 | – | – | 1.5 | 1.2 | 1.5 | 1.5 | 1.2                                              | 1.5 |
| Connector dimension | –   | – | – | –   | –   | –   | –   | 16 mm <sup>2</sup><br>In general: height ≥ width |     |

### Minimum layer thickness IPS e.max CAD – Layering technique

|                  |   |   |   |   |                            |     |   |   |   |
|------------------|---|---|---|---|----------------------------|-----|---|---|---|
| Incisal/occlusal | – | – | – | – | 0.8                        | 1.0 | – | – | – |
| Circular         | – | – | – | – | 0.8                        | 0.8 | – | – | – |
| Design type      | – | – | – | – | Supporting the tooth shape |     | – | – | – |

**IPS e.max CAD is the high-strength component of the restoration and must, therefore, always make up at least 50% of the total layer thickness of the restoration. The total layer thickness of the restoration (depending on the type of restoration) consists of:**

|                                                       |     |     |     |     |     |     |     |     |
|-------------------------------------------------------|-----|-----|-----|-----|-----|-----|-----|-----|
| Total layer thickness of the restoration              | 0.8 | 1.0 | 1.2 | 1.5 | 1.8 | 2.0 | 2.5 | 3.0 |
| Minimum thickness of the IPS e.max CAD framework      | 0.4 | 0.5 | 0.6 | 0.8 | 1.0 | 1.1 | 1.3 | 1.6 |
| Maximum layer thickness of the IPS e.max Ceram veneer | 0.4 | 0.5 | 0.6 | 0.7 | 0.8 | 0.9 | 1.2 | 1.4 |

Dimensions in mm

## Finishing

Observe the grinding instrument recommendations<sup>[4]</sup> and the minimum layer thicknesses when finishing and contouring the ceramic structures. Adjustments by grinding should be carried out while the restoration is still in the pre-crystallized (blue) state, at low speed and exerting only light pressure to prevent delamination and chipping at the margins. Overheating of the ceramic must be avoided. Smooth out the attachment point of the block, paying particular attention to the proximal contacts. If required, carry out individual shape adjustments. Finish the functional areas (occlusal contact surfaces) of the restoration with a fine diamond to smooth out the surface structure created by the CAD/CAM process. Do not "post-separate" bridge framework connectors. This may result in undesired predetermined breaking points, which will subsequently compromise the stability of the all-ceramic restoration. If desired, the restoration can be carefully tried in intraorally in the blue state in order to adjust the occlusion/articulation. Always sonicate the restoration in an ultrasonic water bath or steam-clean them with a steam jet before crystallization. Make sure that the restoration is completely free of any milling additive residue of the CAD/CAM milling unit before further processing. Residue of the milling additive remaining on the surface may result in bonding problems and discolouration. The restoration must not be blasted with Al<sub>2</sub>O<sub>3</sub> or glass polishing beads.

## Completion

### – Polishing technique (self glaze) on the blue restoration

For polishing, observe the grinding instrument recommendations<sup>[4]</sup>. Avoid overheating of the restoration. Use diamond rubber polishers<sup>[5]</sup> for pre-polishing, and high-gloss rubber polishers<sup>[6]</sup> for high-gloss polishing. After polishing, sonicate the restoration in an ultrasonic water bath or steam-clean it with a steam jet. Subsequently, attach the restoration to the crystallization pin as described under "Fixation of restorations on the IPS e.max CAD Crystallization Pin". Place the ceramic structure either on the IPS Speed Tray or on the IPS e.max CAD Crystallization Tray and position the tray in the centre of the furnace. The firing program is selected based on the material and firing tray used (see "Crystallization and firing parameters").

### – Staining technique on the blue restoration

#### a) Glazing with IPS e.max CAD Crystall./Glaze Spray, followed by speed crystallization (speed crystallization and glaze firing in one step).

Position the restoration as described under "Fixation of restorations on the IPS e.max CAD Crystallization Pin". Shake the IPS e.max CAD Crystall./Glaze Spray vigorously immediately before use until the mixing ball in the container is moving freely (approximately 20 seconds). If the spray is not sufficiently shaken, mainly the propellant is discharged with a spraying burst. As a result, the restoration is insufficiently coated with glazing powder. Observe a distance of 10 cm between the nozzle and the surface to be sprayed. Hold the spray can as upright as possible during spraying. Spray the restoration from all sides with short bursts while simultaneously rotating the restoration so that an even, covering layer is created. Shake the can again between individual sprays. Subsequently, spray the restoration a second time from all sides with short bursts while simultaneously rotating the restoration so

that an even, covering layer is created. Shake the can again between individual sprays. Wait briefly until the glaze layer is dry and has assumed a whitish colour. Areas that do not show an even layer have to be sprayed again. Then position max. 2 restorations on the IPS Speed Tray and conduct the speed crystallization firing using the stipulated firing parameters (speed crystallization). Observe the directions under "How to proceed after firing". For adjustments, please observe the information under "Corrective firing".

- b) **Staining and glazing** using IPS e.max CAD Crystall./Glaze Paste/Fluo, followed by crystallization and stains/glaze firing in one step. Position the restoration as described under "Fixation of restorations on the IPS e.max CAD Crystallization Pin". Subsequently, apply IPS e.max CAD Crystall./Glaze Paste/Fluo evenly on the outer surfaces of the restoration using a brush. If the ready-to-use glaze needs to be diluted, it may be mixed with a small amount of IPS e.max CAD Crystall./Glaze Liquid. Do not apply the glaze too thickly. Avoid "pooling", especially on the occlusal surface. Too thin a glaze layer may lead to an unsatisfactory gloss. If characterizations are desired, the restoration can be individualized using IPS e.max CAD Crystall./Shades and/or IPS e.max CAD Crystall./Stains before crystallization firing. Extrude the ready-to-use shades and stains from the syringe and mix them thoroughly. The shades and stains may be slightly thinned using IPS e.max CAD Crystall./Glaze Liquid. However, the consistency should still remain pasty. Apply the mixed shades and stains directly to the unfired glaze layer using a fine brush. Then place the restoration in the centre of the IPS e.max CAD Crystallization Tray or place a maximum of 6 units on the tray and conduct the crystallization firing using the stipulated firing parameters for the respective translucency. Observe the directions under "How to proceed after firing". For adjustments, please observe the information under "Corrective firing".
- c) **Staining and glazing** using IPS e.max CAD Crystall./Glaze Spray, followed by crystallization and stain/glaze firing in one step. Position the restoration as described under "Fixation of restorations on the IPS e.max CAD Crystallization Pin". Extrude the ready-to-use shades or stains from the syringe and mix them thoroughly. The shades and stains may be slightly thinned using IPS e.max CAD Crystall./Glaze Liquid. However, the consistency should still remain pasty. Apply the mixed shades and stains directly to the blue restoration using a fine brush. Spray the restoration with IPS e.max CAD Crystall./Glaze Spray from all sides to coat it with an even layer. Proceed according to the instructions under a). Then place the restoration in the centre of the IPS e.max CAD Crystallization Tray or place a maximum of 6 units on the tray and conduct the crystallization firing using the stipulated firing parameters for the respective translucency. Observe the directions under "How to proceed after firing". For adjustments, please observe the information under "Corrective firing".

#### Corrective firing

If additional characterizations or adjustments are required after crystallization, a corrective firing using IPS e.max CAD Crystall./Shades and Stains and Glaze can be conducted. Also use the IPS e.max CAD Crystallization Tray for the corrective firing cycle. Use the IPS e.max CAD Crystall./Add-On including the respective mixing liquid to carry out minor shape adjustments (e.g. proximal contact points). The adjustments can be carried out in both the crystallization and corrective firing sequences.

#### – On the tooth-coloured restoration, crystallization without the application of materials

- a) **Staining technique:** Stain/Glaze firing of the tooth-coloured restoration with either **IPS e.max CAD Crystall./ or IPS Ivocolor materials**.

Position the restoration as described under "Fixation of restorations on the IPS e.max CAD Crystallization Pin". For characterization and glazing, use either IPS e.max CAD Crystall./Shades, Stains, Glaze (see b under "Corrective/stains/glaze firing") or IPS Ivocolor.

**When using IPS Ivocolor:** To ensure better wetting, a small quantity of IPS Ivocolor Mixing Liquid may be slightly rubbed into the area that needs to be characterized. Mix the IPS Ivocolor Shades and Essences to the desired consistency using the respective IPS Ivocolor Liquids. More intensive shades are achieved by repeating the staining procedure and firing, not by applying thicker layers. Use the IPS Ivocolor Shades Incisal to imitate the incisal area and create the effect of translucency on the crown in the incisal and occlusal third. The cusps and fissures can be customized using the Essences. Then place the restoration in the centre of the IPS e.max CAD Crystallization Tray or place a maximum of 6 units on the tray and conduct the crystallization firing using the stipulated firing parameters for the respective translucency.

IPS e.max CAD Crystall./Shades, Stains, Glaze and IPS Ivocolor Shades, Essence, Glaze must not be mixed with one another or applied one after the other. Observe the directions under "How to proceed after firing".

- b) and c) **Cut-back and layering technique:** using **IPS e.max Ceram materials**. Stain/glaze firing using IPS Ivocolor materials. When performing the cut-back, the minimum layer thicknesses of the framework must be observed. See points "Finishing" and "Minimum layer thicknesses" and as described under "Fixing the restoration on an IPS e.max CAD Crystallization Pin". Then place the restoration in the centre of the IPS e.max CAD Crystallization Tray or place a maximum of 6 units on the tray and conduct the crystallization firing using the stipulated firing parameters for the respective translucency. Observe the directions under "How to proceed after firing".

**Option A: Wash firing using IPS e.max Ceram:** If enough space is available, conduct the wash firing with the required IPS e.max Ceram Transpa Incisal and/or Impulse material. Use the IPS Build-Up Liquids allround or soft to mix the materials. Apply the wash in a thin coat on the entire framework.

**Option B: Wash firing using IPS Ivocolor:** If space is limited or to enhance the in-depth chroma effect, the wash firing can be conducted using IPS Ivocolor Shade, Essence and Glaze. Mix the paste or powder with the IPS Ivocolor Mixing Liquid allround or long-life to achieve the desired consistency and apply it in a thin coat on the entire framework.

**Option C: Wash firing using IPS Ivocolor and IPS e.max Ceram (sprinkle technique):** If space is limited or to enhance the in-depth chroma effect, the sprinkle technique can be used. Apply IPS Ivocolor Shade, Essence and Glaze as described in Option B. Then sprinkle the respective IPS e.max Ceram material (e.g. Dentin) on the restoration using a dry brush. Carefully remove excess with blown air and/or tap it off.

The following applies for all options: The wash (foundation) must be fired before the actual layering procedure is started. Place the honey-comb firing tray in the furnace and conduct the wash (foundation) firing using the respective parameters.

**1<sup>st</sup>/2<sup>nd</sup> Dentin / Incisal firing:** With the IPS e.max Ceram layering materials, both the anatomical shape is completed and a customized esthetic appearance is achieved. The IPS e.max Ceram layering materials can be mixed with either IPS Build-Up Liquid allround or soft. If required, a second firing is conducted.

### Stain/glaze firing using IPS Ivocolor materials

Preparing for stain and glaze firing: Finish the restoration using diamond grinding instruments and give it a true-to-nature shape and surface texture, such as growth lines and convex/concave areas. Areas which should exhibit a higher gloss after glaze firing can be smoothed out and pre-polished using silicone discs. If gold and/or silver dust was used to visualize the surface texture, the restoration has to be thoroughly steam-cleaned with a steam jet. Make sure to remove all gold or silver dust in order to avoid any discolouration.

Stain firing is conducted using IPS Ivocolor Shades and/or Essences, and glaze firing is conducted using IPS Ivocolor Glaze Powder/Fluo or Paste/Fluo. Depending on the situation, the firing sequences may be conducted together or separately. The firing parameters are identical.

Thoroughly steam-clean the restoration with a steam jet and dry with oil-free air. Subsequently, mix the IPS Ivocolor Shades and Essences to the desired consistency using the respective IPS Ivocolor Liquids. In order to promote wetting of the stains and glaze materials, lightly moisten the surface with IPS Ivocolor Mixing Liquid. Subsequently, apply the glazing material in an even layer on the entire restoration. Customize the cusps and fissures using IPS Ivocolor Essence. Apply minor shade adjustments on the applied glazing material using IPS Ivocolor Shades. More intensive shades are achieved by repeating the staining procedure and firing, not by applying thicker layers. The degree of gloss of the glazed surface is controlled via the consistency of the IPS Ivocolor Glaze and the applied quantity, not by means of the firing temperature. In order to achieve a higher degree of gloss, avoid excessively diluting the glaze and/or apply more glazing material. Conduct the stain and glaze firing in a ceramic furnace on the corresponding honey-comb tray using the stipulated firing parameters for the cut-back and layering technique (stain firing using IPS Ivocolor).

### Fixing the restoration on an IPS e.max CAD Crystallization Pin

1. Select the largest possible IPS e.max CAD Crystallization Pin (S, M, L) that best "fills" the inside of the restoration, but does not come into contact with the surrounding crown walls.
2. Subsequently, fill the inside of the restoration with IPS Object Fix Putty or Flow up to the restoration margin. Immediately re-seal the IPS Object Fix Putty/Flow syringe after extruding the material. Once removed from the aluminium bag, the syringe is ideally stored in a resealable plastic bag or a container with a humid atmosphere.
3. Press the selected IPS e.max CAD Crystallization Pin deeply into the IPS Object Fix Putty or Flow material so that it is adequately secured.
4. Smooth out displaced auxiliary firing paste using a plastic spatula, so that the pin is securely in place and the restoration margins are optimally supported.
5. Prevent any IPS Object Fix residue from adhering to the outer surface of the restoration. Remove any adhering residue with a brush dampened with water and then dry.

**Important:** For crystallization, IPS e.max CAD restorations must not be placed directly on the IPS e.max CAD Crystallization Tray/IPS Speed Tray and Pins, i.e. without auxiliary firing paste.

### How to proceed after firing

Remove the restoration from the furnace once the firing cycle is complete (wait for the acoustic signal of the furnace) and allow it to cool to room temperature in a place protected from draft. Hot objects must not be touched with metal tongs. Remove the restoration from the hardened IPS Object Fix Putty / Flow. Remove any residue in an ultrasonic water bath or with the steam jet. Residue must not be removed by blasting with Al<sub>2</sub>O<sub>3</sub> or glass polishing beads. If the restoration needs to be adjusted by grinding<sup>[4]</sup>, make sure that no overheating of the ceramic occurs. Finally, polish the ground areas to a high-gloss finish.

### Shade combination table

For the characterization and adjustment of the shade of IPS e.max CAD restorations, IPS e.max CAD Crystall./Shades, Stains or IPS Ivocolor Shades, Essences are used.

- IPS e.max CAD Crystall./Shades, Stains: For use on blue and tooth-coloured IPS e.max CAD restorations
- IPS Ivocolor Shades, Essences: For use on tooth-coloured IPS e.max CAD restorations

The combination table must be observed.

|                                          | BL1                                                  | BL2 | BL3              | BL4 | A1                | A2   | A3               | A3.5 | A4                   | B1   | B2               | B3   | B4                 | C1   | C2            | C3 | C4                 | D2   | D3               | D4   |  |  |  |  |
|------------------------------------------|------------------------------------------------------|-----|------------------|-----|-------------------|------|------------------|------|----------------------|------|------------------|------|--------------------|------|---------------|----|--------------------|------|------------------|------|--|--|--|--|
| IPS e.max CAD Crystall/<br>Shade         | 0                                                    |     |                  |     | 1                 |      |                  |      |                      | 2    |                  |      |                    | 3    |               |    |                    | 4    |                  |      |  |  |  |  |
| IPS e.max CAD Crystall/<br>Shade Incisal | I1                                                   |     |                  |     |                   |      |                  |      | I2                   |      | I1               |      |                    |      | I2            |    |                    |      |                  |      |  |  |  |  |
| IPS e.max CAD Crystall/<br>Stains        | white, cream, sunset, copper, olive, khaki, mahogany |     |                  |     |                   |      |                  |      |                      |      |                  |      |                    |      |               |    |                    |      |                  |      |  |  |  |  |
|                                          | BL1                                                  | BL2 | BL3              | BL4 | A1                | A2   | A3               | A3.5 | A4                   | B1   | B2               | B3   | B4                 | C1   | C2            | C3 | C4                 | D2   | D3               | D4   |  |  |  |  |
| IPS Ivocolor<br>Shade                    | SD 0                                                 |     |                  |     | SD 1              | SD 2 |                  |      | SD 3                 | SD 4 |                  | SD 5 |                    | SD 6 |               |    |                    | SD 7 |                  | SD 6 |  |  |  |  |
| IPS Ivocolor<br>Shade Incisal            | SI 1                                                 |     |                  |     |                   |      |                  |      |                      |      | SI 2             |      |                    |      | SI 3          |    |                    |      |                  |      |  |  |  |  |
| IPS Ivocolor<br>Essence                  | E 01<br>white                                        |     | E 02<br>creme    |     | E 03<br>lemon     |      | E 04<br>sunset   |      | E 05<br>copper       |      | E 06<br>hazel    |      | E 07<br>olive      |      | E 08<br>khaki |    | E 09<br>terracotta |      | E 10<br>mahogany |      |  |  |  |  |
|                                          | E 11<br>cappuccino                                   |     | E 12<br>espresso |     | E 13<br>terra     |      | E 14<br>profundo |      | E 15<br>ocean        |      | E 16<br>sapphire |      | E 17<br>anthracite |      | E 18<br>black |    | E 19<br>rose       |      | E 20<br>coral    |      |  |  |  |  |
|                                          |                                                      |     |                  |     | E 21<br>basic red |      |                  |      | E 22<br>basic yellow |      |                  |      | E 23<br>basic blue |      |               |    |                    |      |                  |      |  |  |  |  |

### Crystallization and firing parameters

Ceramic furnaces without controlled (long-term) cooling function cannot be used. The ceramic furnace must be calibrated before the first crystallization and regularly every six months thereafter. Depending on the mode of operation, more frequent calibration may be required. Observe the manufacturer's directions.

### Crystallization MO, Impulse, LT, MT, HT

with or without the application of **IPS e.max CAD Crystall./ materials**

| Furnaces<br>Programat      | Stand-by<br>temperature<br>B<br>[°C] | Closing time *<br>S<br>[min] | Heating rate<br>t ↗<br>[°C/min] | Firing<br>temperature<br>T1<br>[°C] | Holding<br>time<br>H1<br>[min] | Heating<br>rate<br>t ↗<br>[°C/min] | Firing<br>temperature<br>T2<br>[°C] | Holding<br>time<br>H2<br>[min] | Vacuum 1<br>11<br>12<br>[°C] | Vacuum 2<br>21<br>22<br>[°C] | Long-term<br>cooling<br>L<br>[°C] | Cooling rate<br>tl<br>[°C/min] |
|----------------------------|--------------------------------------|------------------------------|---------------------------------|-------------------------------------|--------------------------------|------------------------------------|-------------------------------------|--------------------------------|------------------------------|------------------------------|-----------------------------------|--------------------------------|
| P300<br>P500<br>P700       | 403                                  | 6:00                         | 60                              | 770                                 | 0:10                           | 30                                 | 850                                 | 10:00                          | 550/770                      | 770/850                      | 700                               | 0                              |
| P310<br>P510<br>P710       | 403                                  | 6:00                         | 60                              | 780                                 | 0:10                           | 30                                 | 860                                 | 10:00                          | 550/780                      | 780/860                      | 710                               | 0                              |
| CS/CS2/<br>CS3/CS4/<br>CS6 | Select the corresponding program     |                              |                                 |                                     |                                |                                    |                                     |                                |                              |                              |                                   |                                |

### Crystallization LT, MT, HT

with or without the application of **IPS e.max CAD Crystall./ materials**

| Furnaces<br>Programat      | Stand-by<br>temperature<br>B<br>[°C] | Closing time *<br>S<br>[min] | Heating rate<br>t ↗<br>[°C/min] | Firing<br>temperature<br>T1<br>[°C] | Holding<br>time<br>H1<br>[min] | Heating<br>rate<br>t ↗<br>[°C/min] | Firing<br>temperature<br>T2<br>[°C] | Holding<br>time<br>H2<br>[min] | Vacuum 1<br>11<br>12<br>[°C] | Vacuum 2<br>21<br>22<br>[°C] | Long-term<br>cooling<br>L<br>[°C] | Cooling rate<br>tl<br>[°C/min] |
|----------------------------|--------------------------------------|------------------------------|---------------------------------|-------------------------------------|--------------------------------|------------------------------------|-------------------------------------|--------------------------------|------------------------------|------------------------------|-----------------------------------|--------------------------------|
| P300<br>P500<br>P700       | 403                                  | 6:00                         | 90                              | 820                                 | 0:10                           | 30                                 | 840                                 | 7:00                           | 550/820                      | 820/840                      | 700                               | 0                              |
| P310<br>P510<br>P710       | 403                                  | 6:00                         | 90                              | 830                                 | 0:10                           | 30                                 | 850                                 | 7:00                           | 550/830                      | 830/850                      | 710                               | 0                              |
| CS/CS2/<br>CS3/CS4/<br>CS6 | Select the corresponding program     |                              |                                 |                                     |                                |                                    |                                     |                                |                              |                              |                                   |                                |

### Speed crystallization (observe the block concept)

Max. 2 units with or without the application of **IPS e.max CAD Crystall./Glaze Spray** on an IPS Speed Tray

| Furnaces<br>Programat      | Stand-by<br>temperature<br>B<br>[°C] | Closing time *<br>S<br>[min] | Heating rate<br>t ↗<br>[°C/min] | Firing<br>temperature<br>T1<br>[°C] | Holding<br>time<br>H1<br>[min] | Heating<br>rate<br>t ↗<br>[°C/min] | Firing<br>temperature<br>T2<br>[°C] | Holding<br>time<br>H2<br>[min] | Vacuum 1<br>11<br>12<br>[°C] | Vacuum 2<br>21<br>22<br>[°C] | Long-term<br>cooling<br>L<br>[°C] | Cooling rate<br>tl<br>[°C/min] |
|----------------------------|--------------------------------------|------------------------------|---------------------------------|-------------------------------------|--------------------------------|------------------------------------|-------------------------------------|--------------------------------|------------------------------|------------------------------|-----------------------------------|--------------------------------|
| P300<br>P500<br>P700       | 403                                  | 1:00                         | 110                             | 800                                 | 0:00                           | 50                                 | 850                                 | 3:00                           | 690/800                      | 800/850                      | 700                               | 40                             |
| P310<br>P510<br>P710       | 403                                  | 0:30                         | 120                             | 850                                 | 0:00                           | 70                                 | 870                                 | 3:30                           | 690/850                      | 850/870                      | 705                               | 0                              |
| CS/CS2/<br>CS3/CS4/<br>CS6 | Select the corresponding program     |                              |                                 |                                     |                                |                                    |                                     |                                |                              |                              |                                   |                                |

### Corrective firing/stain firing/glaze firing

with **IPS e.max CAD Crystall./ materials**

| Furnaces<br>Programat      | Stand-by<br>temperature<br>B<br>[°C] | Closing time *<br>S<br>[min] | Heating rate<br>t ↗<br>[°C/min] | Firing<br>temperature<br>T1<br>[°C] | Holding<br>time<br>H1<br>[min] | Heating<br>rate<br>t ↗<br>[°C/min] | Firing<br>temperature<br>T2<br>[°C] | Holding<br>time<br>H2<br>[min] | Vacuum 1<br>11<br>12<br>[°C] | Vacuum 2<br>21<br>22<br>[°C] | Long-term<br>cooling<br>L<br>[°C] | Cooling rate<br>tl<br>[°C/min] |
|----------------------------|--------------------------------------|------------------------------|---------------------------------|-------------------------------------|--------------------------------|------------------------------------|-------------------------------------|--------------------------------|------------------------------|------------------------------|-----------------------------------|--------------------------------|
| P300<br>P500<br>P700       | 403                                  | 6:00                         | 90                              | 820                                 | 0:10                           | 30                                 | 840                                 | 3:00                           | 550/820                      | 820/840                      | 700                               | 0                              |
| P310<br>P510<br>P710       | 403                                  | 6:00                         | 90                              | 830                                 | 0:10                           | 30                                 | 850                                 | 3:00                           | 550/830                      | 830/850                      | 710                               | 0                              |
| CS/CS2/<br>CS3/CS4/<br>CS6 | Select the corresponding program     |                              |                                 |                                     |                                |                                    |                                     |                                |                              |                              |                                   |                                |

**Firing parameters for the staining technique**  
with IPS Ivocolor Shade, Essence, Glaze

|                        | Stand-by temperature<br>B [°C] | Closing time *<br>S [min] | Heating rate<br>t↗ [°C/min] | Firing temperature<br>T [°C] | Holding time<br>H [min] | Vacuum 1<br>V1 [°C] | Vacuum 2<br>V2 [°C] | Long-term cooling **<br>L [°C] | Cooling rate<br>tl [°C/min] |
|------------------------|--------------------------------|---------------------------|-----------------------------|------------------------------|-------------------------|---------------------|---------------------|--------------------------------|-----------------------------|
| Stain and Glaze firing | 403                            | IRT/<br>6:00              | 60                          | 710                          | 1:00                    | 450                 | 709                 | 0                              | 0                           |

\* IRT normal mode

\*\* Note: If the layer thicknesses exceed 2 mm, long-term cooling to 500 °C is required.

**Note:** Due to their geometry, the restorations may feature varying layer thicknesses. When the objects cool after the firing cycle, the different cooling speeds in the areas with different thicknesses may result in a build-up of internal tension. In the worst case, these internal tensions may result in fracture of the restoration. By using slow cooling (long-term cooling L), these tensions can be minimized. For monolithic restorations (staining technique) featuring layer thicknesses of more than 2 mm, long-term cooling L must be used.

**Firing parameters for the cut-back and layering technique**  
with IPS e.max Ceram/IPS Ivocolor Shade, Essence, Glaze

|                                           | Stand-by temperature<br>B [°C] | Closing time *<br>S [min] | Heating rate<br>t↗ [°C/min] | Firing temperature<br>T1 [°C] | Holding time<br>H1 [min] | Heating rate<br>t↗ [°C/min] | Firing temperature<br>T2 [°C] | Holding time<br>H2 [min] | Vacuum 1<br>11<br>12 [°C] | Vacuum 2<br>21<br>22 [°C] | Long-term cooling<br>L [°C] | Cooling rate<br>tl [°C/min] |
|-------------------------------------------|--------------------------------|---------------------------|-----------------------------|-------------------------------|--------------------------|-----------------------------|-------------------------------|--------------------------|---------------------------|---------------------------|-----------------------------|-----------------------------|
| Wash firing (foundation)                  | 403                            | IRT/<br>04:00             | 90                          | 650                           | 00:00                    | 20                          | 730                           | 02:00                    | 400/650                   | 650/729                   | 0                           | 0                           |
| 1 <sup>st</sup> Dentin and Incisal firing | 403                            | IRT/<br>04:00             | 90                          | 650                           | 00:00                    | 20                          | 730                           | 02:00                    | 400/650                   | 650/729                   | 0                           | 0                           |
| 2 <sup>nd</sup> Dentin and Incisal firing | 403                            | IRT/<br>04:00             | 90                          | 650                           | 00:00                    | 20                          | 730                           | 02:00                    | 400/650                   | 650/729                   | 0                           | 0                           |
| Stain firing using IPS Ivocolor           | 403                            | IRT/<br>06:00             | 60                          | 710                           | 01:00                    | –                           | –                             | –                        | 450                       | 709                       | 0                           | 0                           |
| Glaze firing using IPS Ivocolor           | 403                            | IRT/<br>06:00             | 60                          | 710                           | 01:00                    | –                           | –                             | –                        | 450                       | 709                       | 0                           | 0                           |
| Add-On with Glaze firing                  | 403                            | IRT/<br>06:00             | 60                          | 710                           | 01:00                    | –                           | –                             | –                        | 450                       | 709                       | 0                           | 0                           |
| Add-On after Glaze firing                 | 403                            | IRT/<br>06:00             | 50                          | 700                           | 01:00                    | –                           | –                             | –                        | 450                       | 699                       | 0                           | 0                           |

\* IRT normal mode

### 3 Safety information

- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, website: [www.ivoclar.com](http://www.ivoclar.com), and your responsible competent authority.
- The current Instructions for Use are available in the download section of the Ivoclar Vivadent AG website ([www.ivoclar.com](http://www.ivoclar.com)).
- The Summary of Safety and Clinical Performance (SSCP) can be retrieved from the European Database on Medical Devices (EUDAMED) at <https://ec.europa.eu/tools/eudamed>. Basic UDI-DI: 76152082ACERA001EQ

#### Warnings

- IPS Natural Die Material Separator contains hexane. Hexane is highly flammable and detrimental to health. Avoid contact of the material with the skin and eyes. Do not inhale vapours and keep away from sources of ignition.
- Do not inhale ceramic dust during finishing. Use an extraction unit and wear a face mask.
- IPS e.max CAD Crystall/Glaze Spray: Do not inhale the spray mist. Container is pressurized. Protect from direct sunlight and temperatures over 50 °C. Do not forcibly open or burn the can, even after it has been completely emptied.
- If applied frequently or for a prolonged period of time, peroxide-based professional bleaching agents (carbamide peroxide; hydrogen peroxide) as well as acidulated phosphate fluorides used for caries prevention can make the surface of existing IPS e.max CAD restorations rough and matte.
- Observe the Safety Data Sheet (SDS) (available in the download section of the Ivoclar Vivadent AG website [www.ivoclar.com](http://www.ivoclar.com)).

#### Disposal information

Remaining stocks or removed restorations must be disposed of according to the corresponding national legal requirements.

#### Residual risks

Users should be aware that any dental intervention in the oral cavity involves certain risks. Some of these risks are listed below:

- Chipping / fracture / decementation of the restoration material may lead to ingestion or aspiration of material and additional dental treatment.
- Cement excess may lead to irritation of the soft tissue / gingiva. Bone loss may occur and periodontitis may develop as the inflammation progresses.

### 4 Shelf life and storage

This product does not require any special storage conditions.

## 5 Additional information

Keep out of the reach of children!

Not all products are available in all countries.

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of use. The user is responsible for testing the materials for their suitability and use for any purpose not explicitly stated in the Instructions.

<sup>[1]</sup> e.g. PrograMill, CEREC/inLab, PlanMill. The complete list is available at [www.ivoclar.com](http://www.ivoclar.com).

CEREC/inLab, PlanMill, CEREC SpeedFire are not registered trademarks of Ivoclar Vivadent AG.

<sup>[2]</sup> e.g. Programat CS6, CS4, CS3, CS2, CS, CEREC SpeedFire

<sup>[3]</sup> The availability of shades may vary depending on the degree of translucency/block size or CAD/CAM machine.

<sup>[4]</sup> Ivoclar Vivadent flow chart "Recommended grinding instruments for extraoral and intraoral use".

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| Transluzenz-<br>stufe               | Block-<br>größen/<br>Farben                            | Verarbeitungstechnik                                                                                    |                                                                                                      |                                                                                                           |                                                                                                          | Restaurationsarten               |                              |                  |                  |                  |                  |                    |
|-------------------------------------|--------------------------------------------------------|---------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------|----------------------------------|------------------------------|------------------|------------------|------------------|------------------|--------------------|
|                                     |                                                        | Polier-<br>technik<br> | Mal-<br>technik<br> | Cut-Back-<br>Technik<br> | Schicht-<br>Technik<br> | Okklusales Veneer <sup>[a]</sup> | Dünnes Veneer <sup>[a]</sup> | Veneer           | Inlay, Onlay     | Teilkronen       | Krone            | 3-gliedrige Brücke |
| <b>HT</b><br>High<br>Translucency   | I12, C14, B40,<br>B40 L<br>in 20 Farben <sup>[f]</sup> | ✓                                                                                                       | ✓                                                                                                    | ✓                                                                                                         |                                                                                                          | ✓                                | ✓ <sup>[d]</sup>             | ✓ <sup>[d]</sup> | ✓ <sup>[d]</sup> | ✓ <sup>[d]</sup> | ✓ <sup>[d]</sup> |                    |
| <b>MT</b><br>Medium<br>Translucency | C14<br>in 7 Farben                                     | ✓                                                                                                       | ✓                                                                                                    | ✓                                                                                                         |                                                                                                          | ✓                                | ✓ <sup>[d]</sup>             | ✓ <sup>[d]</sup> | ✓ <sup>[d]</sup> | ✓ <sup>[d]</sup> | ✓ <sup>[d]</sup> |                    |
| <b>LT</b><br>Low<br>Translucency    | I12, C14, C16,<br>B32 in<br>20 Farben <sup>[f]</sup>   | ✓                                                                                                       | ✓                                                                                                    | ✓                                                                                                         |                                                                                                          |                                  |                              | ✓ <sup>[d]</sup> |                  | ✓ <sup>[d]</sup> | ✓ <sup>[d]</sup> | ✓ <sup>[b]</sup>   |
| <b>MO</b><br>Medium<br>Opacity      | C14<br>in MO 0 – 4 <sup>[g]</sup>                      |                                                                                                         |                                                                                                      |                                                                                                           | ✓                                                                                                        |                                  |                              |                  |                  |                  | ✓ <sup>[c]</sup> |                    |
| <b>I</b><br>Impulse                 | C14<br>in 01 und 02                                    | ✓                                                                                                       | ✓                                                                                                    | ✓                                                                                                         |                                                                                                          | ✓                                | ✓                            | ✓                |                  |                  |                  |                    |

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**Ivoclar Vivadent LLC**
**Ivoclar Vivadent Marketing Ltd.**
**Ivoclar Vivadent AB**
**Ivoclar Vivadent Limited**
**Ivoclar Vivadent, Inc.**