Esthet•X[®] HD

High Definition Micro Matrix Restorative

INSTRUCTIONS FOR USE	ENGLISH	2
Matériau de restauration haute définition à micro-matrice MODE D'EMPLOI	FRANÇAIS	6
Materiale da restauro a micromatrice ad alta definizione ISTRUZIONI PER L'USO	ITALIANO	11
Hochauflösendes Mikromatrix-Füllungsmaterial GEBRAUCHSANWEISUNG	DEUTSCH	16
Material restaurador con micro matriz de alta definición INSTRUCCIONES DE USO	ESPAÑOL	21
Hogedefinitie-micromatrixrestauratiemateriaal GEBRUIKSAANWIJZING	NEDERLANDS	26
High Definition mikro-matrix tandersättningsmaterial BRUKSANVISNING	SVENSKA	31
Mikromatricerestaurationsmateriale i højdefinition BRUGSANVISNING	DANSK	35
Didelės raiškos mikromatricos restauravimo medžiaga NAUDOJIMO INSTRUKCIJOS	LIETUVIŲ K.	39
Augstas izšķirtspējas mikromatricas restaurācijas līdzeklis LIETOŠANAS INSTRUKCIJA	LATVIEŠU	44
Výplň Micro Matrix s vysokou rozlišovací schopností NÁVOD NA POUŽITÍ	SLOVENSKÝ	49
Микро матричен възстановъчен материал с висока прецизност ИНСТРУКЦИИ ЗА УПОТРЕБА	БЪЛГАРСКИ	54
Υλικό αποκαταστάσεων μικρομετρικής μήτρας υψηλής πιστότητας ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ	ΕΛΛΗΝΙΚΑ	59

Manufactured by
Dentsply Caulk
38 West Clarke Avenue
Milford, DE 19963, USA
Tel.: 1-302-422-4511
www.dentsplysirona.com

Dentsply DeTrey GmbH
De-Trey-Strasse 1
78467 Konstanz, Germany
Tel.: 49-7531-583-0
www.dentsplysirona.com

Esthet•X® HD

High Definition Micro Matrix Restorative

INSTRUCTIONS FOR USE - ENGLISH

Caution: This is a medical device. For dental use only by dental professionals.
U.S.A.: Rx only.

1. PRODUCT DESCRIPTION

Esthet•X® HD High Definition Micro Matrix Restorative is a visible light cured, radiopaque, composite restorative material for anterior and posterior restorations of primary and permanent teeth. It is to be used following the application of a suitable dentin/enamel adhesive.

Esthet•X® HD High Restorative Material complies with ISO 4049:2019 Polymer-based restorative materials.

1.1 Indications for use

- Esthet•X® HD Restorative material is indicated as a direct restorative for all cavity classes in anterior and posterior teeth.
- Esthet•X® HD Restorative material is indicated as a direct esthetic veneering restorative material and cosmetic reshaping (e.g., diastema closure, incisal lengthening).
- Esthet•X® HD Restorative material is indicated for the indirect fabrication of inlays and on-lays.

1.2 Contraindications

Esthet•X® HD Restorative material is contraindicated for use with patients who have a history of hypersensitivity to methacrylate resins.

1.3 Delivery Forms (some delivery forms may not be available in all countries)

- This one-component, visible light cured composite restorative is available in:
- Pre-dosed Compules® Tips
 - 3 g Easy•Twist syringes

The regular body shades are optimized to match the VITA^{®1} system. Additional available body, opacious dentin and translucent enamel shades are available for esthetically demanding applications existing beyond the traditional VITA® System. A custom blended shade guide is available to guide the clinician in shade selection and blending technique.

Esthet•X® HD Restorative material may be used with separately available BisGMA-based tints and opaques at the discretion and responsibility of the dental practitioner for individual custom characterization of the ultimate esthetic direct-placement restoration.

Available Shades

Opacious Dentin	Regular Body		Translucent Enamel
White-O = W-O	White = W	B5/Dark Yellow	Clear Enamel = CE
A2-O	Extra Light = XL	C1	White Enamel = WE
A4-O	A1	C2	Yellow Enamel = YE
B2-O	A2	C3	Amber Enamel = AE
C1-O	A3	C4	Grey Enamel = GE
C4-O	A3.5	C5/Extra Grey-Brown	
D3-O	A4	D2	
	B1	D3	
	B2	Universal	
	B3		

1.4 Composition

The resin matrix of Esthet•X® HD Restorative material consists of a Bis-GMA adduct, a Bis-EMA adduct, and triethylene glycol dimethacrylate, Camphorquinone (CQ), Photoinitiator, Stabilizer, Pigments.

The filler combination consists of barium fluoroborosilicate glass with a mean particle size below 1 µm and nanofiller silica (particle size 0.04 µm); 52-54 volume %.

1.5 Compatible Adhesives

Esthet•X® HD Restorative material is used following application of a suitable dentin/enamel adhesive and is chemically compatible with conventional (meth)acrylate-based dentin/enamel adhesives including Dentsply Sirona adhesives designed for use with visible light cured composite restoratives (see complete Instructions for Use of selected adhesive).

1.6 Compatible Cements

Inlays and Onlays fabricated from Esthet•X® HD Restorative are compatible with resin-based cements including Dentsply Sirona cements designed for use with composite inlays/onlays (see complete Instructions for Use of selected cement).

2. GENERAL SAFETY NOTES

Be aware of the following general safety notes and the special safety notes in other chapters of these Instructions for Use.

Safety alert symbol

This is the safety alert symbol. It is used to alert you to potential personal injury hazards. Obey all safety messages that follow this symbol to avoid possible injury.

2.1 Warnings

Esthet•X® HD Restorative material contains polymerizable methacrylate monomers which may be irritating to skin, eyes and oral mucosa and may cause sensitization by skin contact and allergic contact dermatitis in susceptible persons.

- **Avoid eye contact** to prevent irritation and possible corneal damage. In case of contact with eyes, rinse immediately with plenty of water and seek medical attention.
- **Avoid skin contact** to prevent irritation and possible allergic response. In case of contact, reddish rashes may be seen on the skin. If contact with skin occurs, immediately remove material with cotton and wash thoroughly with water and soap. In case of skin sensitization or rash, discontinue use and seek medical attention.
- **Avoid contact with oral soft tissues/mucosa** to prevent inflammation. If accidental contact occurs, immediately remove material from the tissues. Flush mucosa with plenty of water after the restoration is completed and expectorate/evacuate the water. If sensitization of mucosa persists, seek medical attention.

- Esthet•X® HD Restorative material curing reaction is exothermic. Increased temperature may irritate pulp or soft tissue. Cover dentin close to the pulp (less than 1 mm) with a hard setting calcium hydroxide liner (Dycal® Liner) and/or treat according to adhesive manufacturer's instructions.
- Esthet•X® HD Restorative material is light cured only. Incomplete light-curing may lead to compromised physical properties, residual uncured monomer, post-operative sensitivity or premature bond/restoration failure. See Step-by-Step Instructions for wavelength, output and maximum increment requirements. Check output and curing cycle before each use.
- When using a Compules® Tips Gun, dental professionals and patients should wear protective eyewear.
- Esthet•X® HD Restorative material should not be used with patients who have a known hypersensitivity to any of the components.

2.2 Precautions

- This product is intended to be used only as specifically outlined in these Instructions for use. Any use of this product inconsistent with the Instructions for use is at the discretion and is the sole responsibility of the practitioner.
- Devices marked "single use" on the labeling are intended for single use only. Discard after use. Do not reuse in other patients in order to prevent cross-contamination.
- Contact with saliva and blood during composite placement may cause failure of the restoration. Use of rubber dam or adequate isolation is recommended.
- Wear suitable protective eyewear, mask, clothing and gloves. Protective eyewear is recommended for patients.
- Use only clean and properly disinfected and sterilized, or ready to use accessories. Do not use dispensers which have not been properly pre-use processed or reprocessed. Before use and each re-use dispensers must be properly cleaned and autoclaved (see Hygiene, Processing and Disposal below).
 - As additional precautionary measure, dispensers may be protected from gross debris but not from all contamination by applying a protective barrier.
- Compules® Tips are provided ready for use and should not be processed. Use only new, clean unused Compules® Tips, stored in covered storage, such as drawers and cabinets, away from potential contamination, and removed from their original packaging using either clean gloved hands or a clean suitable instrument at time of use. To reduce risk of infection, do not use tips if packaging is damaged or compromised.
- The syringe cannot be reprocessed. To prevent the syringe and from exposure to spatter or spray of body fluids or contaminated hands it is mandatory that the syringe is handled with clean/disinfected gloves.
 - As additional precautionary measure, syringe may be protected from gross debris but not from all contamination by applying a protective barrier.
- Do not insert or use restorations which have not been properly processed. Before delivery, restorations must be properly cleaned and disinfected to reduce the risk of cross contamination (see Hygiene, Processing and Disposal).
- Material should extrude easily from the Compules® Tip. Use a gentle, even motion when exerting pressure on the Compules® Tips Gun. DO NOT USE EXCESSIVE FORCE. Excessive pressure may result in unanticipated extrusion of the material or cause the Compules® Tip to rupture or eject from the Compules® Tips Gun. If the material doesn't extrude easily remove from the operating field, discard the Compules® Tip, and install a new Compules® Tip. Use of Compules® Tips with the Dentsply Caulk Compules® Tips Gun is recommended. Use with another manufacturer's dispensing device is at the discretion and sole responsibility of the practitioner.
- Do not use the product if package is broken.
- Numerous factors can affect restoration longevity. Marginal leakage can result in restoration loss or recurrent caries. Check restoration marginal integrity at every periodic examination
- The syringe should be tightly closed immediately after use.
- **Interactions:**
 - Eugenol and hydrogen peroxide-containing materials should not be used in conjunction with this product because they may interfere with hardening and cause softening of the polymeric components of the material.
 - Esthet•X® HD Restorative material is light-cured material. Therefore, it should be protected from ambient light. Proceed immediately once the material has been placed.
 - If mineral-impregnated (e.g., ferric compounds) retraction cords and/or hemostatic solutions are used in conjunction with adhesive procedures, marginal seal may be adversely affected, allowing microleakage, subsurface staining and/or restoration failure. If gingival retraction is necessary, use of plain, non-impregnated cord is recommended.

2.3 Adverse Reactions

- Product may irritate the eyes and skin.
 - **Eye contact:** irritation and possible corneal damage.
 - **Skin contact:** irritation or possible allergic response. Reddish rashes may be seen on the skin.
 - **Mucous Membranes:** inflammation (see Warnings).
- Product may cause pulpal effect (see Warnings).

2.4 Storage

Inadequate storage conditions may shorten the shelf life and may lead to malfunction of the product.

- Store with original caps tightly closed.
- Store at temperatures between 2 °C/35 °F-24 °C/75 °F.
- Use the product at room temperature.
- Keep out of sunlight.
- Do not freeze.
- Excessive humidity can adversely affect the properties of Esthet•X® HD Restorative material. Protect from moisture.
- Do not use after expiration date.
- Do not use if the lot number and/or expiration date are missing or illegible.

The TruMatch® Restorative material shade guide and Compules® Tips Gun should be processed and stored according to the instructions in Hygiene below.

3. STEP-BY-STEP INSTRUCTIONS

Technique Tips

- Esthet•X® HD Restorative material is a "sculptable," non-sticky, material that resists slumping and allows the operator to create precise anatomical form.
 - When using Esthet•X® HD Restorative material, it is best not to "spread" the material, but to "push" or "sculpt" the material.
 - It is recommended to dispense/extrude material from a Compules® Tip with a slow, steady force.
 - Rapid and/or forceful extrusion of the material from a Compules® Tip may cause material to appear whitish and/or separated. Upon placement and manipulation, material will coalesce without compromise in strength or esthetics.
- Esthet•X® HD Restorative material was developed to provide the highest surface luster in order to replicate natural enamel. To achieve maximum polish and maintenance it is highly recommended to fully complete the polish of Esthet•X® HD Restorative material by finishing the restoration with Enhance® Finishing System and polishing with the Enhance® PoGo® One Step Diamond Micro-Polisher System and/or Prisma® Gloss™ Composite Polishing Pastes.

¹ VITA is not a registered trademark of Dentsply Sirona.

Optimizing Shade Selection

1. The opacity of Esthet•X® HD Restorative material Body Shades has been scientifically designed to allow the tooth structure enamel and dentin to blend together with the composite with a life-like translucent result. Additional Opaque and Enamel Shades allow replacement of each component of lost tooth structure. Layered opacities of material visually blend to form a restoration with natural vitality and translucence.
2. The final shade will ultimately result from an interaction of the thickness of the layer(s) of restorative composite and surrounding tooth structure (enamel and/or dentin). The intensity of the shade is related to the thickness or thinness of opacious, body and/or translucent shade layers of the Esthet•X® HD Restorative material used. Restoration shade approximation to the selected shade tab of the TruMatch® shade guide occurs with restorations allowing adequate thicknesses of each recipe component, e.g., Class IV, labial-lingual Class III, etc. For restorations less than 2 mm total thickness, e.g., conservative cosmetic veneering, and/or where significant shade lightening is desired, placement of a liquid opaquener may provide a more desirable foundation for subsequent shade layering.
3. In all cases, the use of a mock-up is an excellent method for shade verification. The technique includes the use of the shades and opacities as prescribed by the selected TruMatch® shade tab in the appropriate thicknesses to approximate that of the final Esthet•X® HD Restorative material restoration. The materials are placed on the clean, hydrated unetched tooth in the area of the cavity or restoration on the tooth. Each layer is applied without etch or bonding agent, then cured thoroughly. Layering and/or altering shades in varying thickness will allow trial customization of the final shade. The dentist, patient and ancillary dental support personnel can co-jointly view the shade selection result. After viewing, the composite material can be removed by the use of a dental explorer or scaler.
4. Be aware of room and ambient light effects on shade selection; incandescent versus fluorescent lighting versus natural daylight (ideal light is northern exposure/outdoors/overcast).
5. Room accessories can influence shade selection as with colored walls/wallpaper reflection. Also, the patient dental napkin should be removed. Observe tooth shades for short periods of time, using a blue background to "neutralize" the effect of extended shade viewing comparisons. Viewing a blue background will have a relaxing effect on the viewing operator's photo/color optic discrimination. Eyes should be rested. It is useful to have ancillary corroboration on shade selection by dental personnel as well as by the dental patient.

3.1 Direct Restoration

3.1.1 Shade Selection

Before selecting the shade, teeth should be clean, hydrated and free of extrinsic material or stain. The TruMatch® shade guide included in the Introductory Kit and the Complete Organizer System has been created to simulate the result of layering specific shades. Select the shade tab that simulates best the final desired result. Following the recipe on the shade tabs will create the desired shade. Mixing layered shades and various thicknesses will allow customized shade adjustability. Please refer to the Optimizing Shade Selection section above. Additionally, use of the Technique Guide will detail the step-by-step procedure for accurate shade reproduction.

3.1.2 Cavity Preparation

- 2.1 Anterior Restorations: Use conservative cavity preparations for all Class III, IV and Class V restorations. Refinement of the cavo-surface margin (beveling) for enhancement of acid-etching, enamel bonding, and finishing procedures is recommended. **Technique Tip:** Esthetic masking of underlying intrinsic staining may require deeper preparation, allowing for adequate restorative material blending.
- 2.2 Posterior Restorations: Cavity design requirements are essentially a conventional preparation with refinement of the cavo-surface margin for enhancement of acid etching. No residual amalgam or other base material should be left in the internal forms of the preparation which would interfere with light transmission and the hardening of the restorative.

3.1.3 Placement of Matrix

The use of a Mylar, Bimatrix, sectional (e.g. Palodent® Sectional Matrix System) or thin matrix band (e.g. AutoMatrix® Retainerless Matrix System) and subsequent burnishing of the matrix band will improve final interproximal contact and contour. PREWEDGING/BiTine® RING PLACEMENT IS ADVOCATED TO ACHIEVE SLIGHT SEPARATION AND FACILITATE ACCEPTABLE PROXIMAL CONTACT.

3.1.4 Pulp Protection, Tooth Conditioning/Dentin Pretreatment, Adhesive Application

Refer to adhesive manufacturer's instructions for pulp protection, tooth conditioning and/or adhesive application. Once the surfaces have been properly treated, they must be kept uncontaminated. Proceed immediately to placement of Esthet•X® HD Restorative material.

3.1.5 Placement of Esthet•X® HD Micro Matrix Restorative Material

- 5.1 Easy•Twist Syringe: Remove the cap from the Easy•Twist syringe by using lateral force. It is easier to remove the cap in a "snap-off" rather than a "pull-off" motion. Dispense the necessary amount of Esthet•X® HD Restorative material from the Easy•Twist syringe onto the mix pad by turning the handle slowly in a clockwise direction. To prevent oozing of the material when dispensing is completed, point the front tip of the Easy•Twist syringe upwards and turn the handle counter-clockwise and recap the Easy•Twist syringe. **Technique Tip:** Protect dispensed Esthet•X® HD Restorative material from premature polymerization via ambient light while current increment is adapted and light cured.
- 5.2 Pre-dosed Compules® Tips: Load Compules® Tips Gun with predosed Compules® Tip. Insert a Compules® Tip into the notched opening of the Compules® Tips Gun barrel. Be certain that the collar on the Compules® Tip is inserted first. Remove the colored cap from the Compules® Tip. The Compules® Tip may be rotated 360° to gain the proper angle of entrance into the cavity. To dispense the material into a cavity preparation, use a slow, steady pressure. Do not use excessive force. To remove the used Compules® Tip, be sure that the Compules® Tips Gun plunger is pulled back completely by allowing the handle to open to its widest position. Apply a downward motion to the front end of the Compules® Tip and remove. Sterilize Compules® Tips Gun for storage according to instructions (see Hygiene, Processing and Disposal section).

Material Extrusion – To reduce risk of injury due to excessive force.

- Apply slow and steady pressure on the Compules® Tips Gun.
- Do not use excessive force – Compules® Tip rupture or ejection may result.

5.3 Optional Flowable Liner: Placement of a compatible flowable liner (available separately) prior to placement of Esthet•X® HD Restorative material is optional. Follow manufacturer's instructions for use.

Handling

Use only clean and properly disinfected and sterilized, or ready to use accessories.

3.1.6 Anterior Placement

Dispense Esthet•X® HD Restorative material directly into the cavity preparation/tooth surface from the Compules® Tip using slow, steady pressure. Excessive force is not necessary. Alternatively, material may be expressed onto a clean pad from the Compules® Tip or Easy•Twist syringe and carried to the preparation with suitable placement instrument. Adapt, contour and shape with appropriate composite instruments. Material may be placed and light cured in increments up to 2 mm (see Curing, Step 3.1.8). **Technique Tip:** For esthetic blending of shades, shade layers may be individually cured to form a foundation for subsequent shade placement.

3.1.7 Posterior Placement

Dispense Esthet•X® HD Restorative material directly into the cavity preparation from the Compules® Tip using slow, steady pressure. Excessive force is not necessary. Alternatively, material may be expressed onto a clean pad from the Compules® Tip or Easy•Twist syringe and carried to the preparation with suitable placement instrument.

- 7.1 Class I and V: Esthet•X® HD Restorative material may be placed and adapted in 2 mm increments, light curing each increment (see Curing, Step 3.1.8).
- 7.2 Class II: Beginning in the proximal box, firmly adapt a 2 mm increment of Esthet•X® HD Restorative material to cavity walls and stabilize matrix band by stepping a clean nonserrated condenser (plugger) from middle to edges, reducing incorporated air and eliminating voids. Light cure (see Curing, Step 3.1.8). The remainder of the preparation may be filled in up to 2 mm increments, each followed by light cure (see Curing, Step 3.1.8). Prior to light curing, contour and shape the final increment with the operator's choice of clean carving and burnishing instruments. **Technique Tips:** To minimize finishing time, contour marginal ridge first before occlusal anatomy. Carve margins and anatomy to final form. Esthet•X® HD Restorative material resists slumping, allowing carving of the majority of anatomical form prior to visible light curing.

3.1.8 Curing

Light cure each area of the restoration surface with a suitable visible light curing unit designed to cure materials containing camphorquinone (CQ) initiator, i.e. spectral output containing 470 nm. Minimum light output must be at least 550 mW/cm² exposure for at least 20 seconds. Some advanced performance curing units have been shown to cure 2 mm increments of most shades of Esthet•X® HD Restorative material in 10 seconds. Refer to curing light manufacturer's recommendations for compatibility and curing recommendations. The Esthet•X® HD Restorative material should be additionally exposed to the curing unit through the proximal, lingual, and buccal enamel walls following matrix removal for the recommended time.

To reduce risk of inadequate polymerization due to insufficient curing.

- Check compatibility of curing light.
- Check curing cycle.
- Check curing output before each procedure.

3.1.9 Finishing and Polishing

- 9.1 Begin finishing immediately after curing. Gross excess may be removed, and general outline form established with Prisma® Finishing Burs or other carbide or diamond finishing instruments. Additional finishing is recommended by the use of Enhance® Finishing System. See manufacturer's complete instructions for use.
- 9.2 To achieve a very high luster on Esthet•X® HD Restorative material, it is necessary to complete the polishing. Enhance® PoGo® One Step Diamond Micro-Polisher System and/or Prisma® Gloss™ Composite Polishing Pastes are recommended. See manufacturer's complete instructions for use. Once polished with Enhance® PoGo® polisher and/or both Prisma® Gloss™ pastes, Esthet•X® HD Restorative material will maintain its high surface gloss due to its micro matrix technology.

3.2 Fabrication of Indirect Inlay/Onlay

3.2.1 Shade Selection

Final desired shade selection should be accomplished prior to tooth preparation. Before selecting the shade, teeth should be clean, hydrated and free of extrinsic material or stain. The TruMatch® shade guide included in the Introductory Kit and the Complete Organizer System has been created to simulate the result of layering specific shades. Select the shade tab that simulates best the final desired result. Following the recipe on the shade tabs will create the desired shade. Mixing layered shades and various thicknesses will allow customized shade adjustability. Please refer to the Optimizing Shade Selection section.

3.2.2 Cavity Preparation

Cavity design requirements are essentially a conventional preparation. Rounding of internal angles and refinement of the cavo-surface margin for enhancement of enamel bonding and finishing procedures is recommended. No residual amalgam or other base material should be left in the internal forms of the preparation which would interfere with light transmission and the hardening of the luting cement. Refer to adhesive and/or luting cement manufacturer's instructions for base/liner/pulp protection requirements.

3.2.3 Impression, Master Cast Fabrication

- 3.1 Make an accurate impression of the preparation following manufacturer's instructions for impression material chosen. A material that allows two pours is recommended. Alternatively, two impressions may be made.
- 3.2 Prepare and place provisional restoration on the prepared tooth. Use non-eugenol containing material(s).
- 3.3 Two casts will be fabricated. One working cast will be indexed and separated, per usual crown and bridge procedures. The second master cast will be used for final contour adjustment. Following impression material manufacturer's instructions for disinfection and impression casting, use of die stone is recommended for cast fabrication. Plating of the impression is not recommended. Additionally, an opposing cast should be prepared and articulated.

3.2.4 Restoration Fabrication

- 4.1 Block out any preparation undercuts present. Apply separating medium and/or die spacer to separated die. In most cases, the restoration may be fabricated in no more than 3 increments, each up to 4 mm depth.
- 4.2 Apply first increment/shade, creating restoration body. Adapt to die, staying short of margins. Recommended light curing is accomplished by placing in Triad® Light Curing Unit (Dentsply Trubyte) for 2 minutes. Apply second, body layer, allowing cut-back for occlusal anatomy and final proximal contour. Repeat Triad® Light Curing. Prior to final "enamel" shade placement, custom staining with appropriate compatible stains may be accomplished per manufacturer's instructions. Apply final layer, slightly overfilling and covering all margins. Replace die into articulated model. Establish all external, proximal and occlusal contacts and anatomy. Slight lubrication of adjacent and opposing model teeth is recommended. Cure briefly (10 seconds) with hand-held visible curing light unit to fix contours. Remove die with restoration, place in Triad Unit for final 2 minute curing.
- 4.3 When removing restoration from die, it may be necessary to scrape die stone away from restoration margins, to prevent accidental chipping of restoration. Clean any residual die stone from restoration. Carefully trim away visible flash beyond preparation margins with acrylic bur.
- 4.4 Gently sandblast the internal surfaces with 50 µm alumina abrasive. Remove visibly undercut material.
- 4.5 Seat restoration onto preparation of uncut master cast, making adjustments as needed. Check for marginal integrity and overall fit and contour. Additional increments may be added if needed by roughening surface and applying compatible adhesive per manufacturer's instructions, followed by placement and light curing as outlined above.

3.2.5 Finishing and Polishing (laboratory)

Complete as outlined in the above Finishing and Polishing Section, Step 3.1.9 for Direct Restorations.

3.2.6 Pre-Cementation Surface Preparation

- 6.1 Sandblast the internal surface of the restoration using 50 µm aluminum oxide at a max pressure of 2.5 bar (35 PSI).
- 6.2 Clean and disinfect the restoration as outlined in Hygiene section below.

**Contamination – To reduce risk of infection**

- Restorations should be polished and sandblasted internally before cleaning and disinfection.
- Restorations should be cleaned and disinfected immediately prior to delivery.
- Follow instructions in Hygiene section below before delivery.

3.2.7 Cementation

Recommended cementation technique is bonding using an adhesive esthetic resin cement. Follow adhesive and cement manufacturer's instructions for prepared tooth and restoration surface pretreatments.

3.2.8 Adjusting, Finishing and Polishing (clinically)

Following cementation, make all necessary occlusal adjustments and polish any clinically adjusted surfaces as outlined in the above Finishing and Polishing Section, Step 3.1.9 for Direct Restorations. All shades of Esthet•X® HD Restorative are radio-opaque, with a 1 mm radio-opacity equivalent to 1.9 mm radio-opacity of aluminum. Aluminum has a radio-opacity equivalent to that of dentin. Thus 1 mm of material having a radio-opacity equivalent to 1 mm of aluminum has a radio-opacity equivalent to that of dentin.

4. HYGIENE, PROCESSING AND DISPOSAL**Cross-contamination – To reduce risk of infection.**

- Do not reuse single use products. Dispose of in accordance with local regulations.
- Do not reuse Compules® Tips. Properly dispose of used Compules® Tips.
- Syringe cannot be reprocessed. Dispose of contaminated syringe in accordance with local regulations.
- During treatment clinicians with patient contact should not handle the syringe without a protective barrier.
- To prevent the syringes from exposure to spatter or spray of body fluids or contaminated hands it is mandatory that the syringes are handled with clean/disinfected gloves, outside of the patient treatment operator.
- Do not reuse syringes if contaminated.
- Protective barriers are designed for single use and must be disposed of after each use in accordance with local regulations. The barrier is not a replacement for applicable cleaning, disinfection, and sterilization.
- Incidental contact of the syringe with water, soap or a water-based hospital-level disinfection solution will not damage syringe body. Do not allow any solution contact with contained material. Discard composite material that has been in contact with any fluid or non-sterile instrument. Repeated liquid contact may damage label. Dry the syringe with a lint-free single-use cloth.
- NOTE:** Vigorous wiping can destroy the label. Wipe syringe gently.
- Reprocess reusable products as described below.
- When not in use, all supplies and instruments should be in covered storage, such as drawers and cabinets, and away from potential contamination.
- Dispensing materials with clean/disinfected gloves or hands or suitable transfer forceps in a separate room, bringing only what will be used into the operator room is strongly recommended.

**4.1 Indirectly Fabricated Inlays and Onlays****Cross-contamination – To reduce risk of infection**

- Restorations must be polished before cleaning and disinfection.
- Restorations must be processed according to the instructions below prior to delivery.
- Locally applicable legal regulations and the hygiene standards applicable for a dental practice must be observed.
- Use only the validated processing procedures specified.
- The responsibility for the microbial processing of the Restoration lies with the user.
- The equipment and devices must be properly maintained and serviced at regular intervals.
- The fabricator (laboratory or in-office dental technician) of the restoration must inform the dentist of the need to process the device before using or inserting it in the patient's mouth.

Process		✓ Acceptable
Cleaning	Mechanical (Ultrasonic)	✓
	Manual	✓
	Automated (Washer/Disinfector)	
Disinfection	Automated (Washer/Disinfector)	
	Manual	✓
Sterilization	Immersion	✓
	Steam Autoclave	
Single Use	Dispose of after Use	
Reusable/Not Reprocessable	Dispose of if contaminated	

Instructions for Cleaning and Disinfecting Indirectly Fabricated Restorations	
Warnings	<ul style="list-style-type: none"> These instructions are for use ONLY for an indirectly fabricated restoration, prior to intraoral delivery. High-level disinfection procedures have not been validated. Intermediate-level disinfection is appropriate for indirect restorations. Do not autoclave in steam autoclave. Do not clean/disinfect in automated washer/disinfector. Follow proper infection prevention activities, such as proper hand washing and donning new puncture and chemical resistant gloves at pertinent steps. Do not use hard wire brushes to clean the restoration as damage can occur. Use only a disinfecting solution which is approved for its efficacy, EPA registered (and/or Health Canada approved) and use in accordance with the IFU of the disinfecting solution manufacturer.
Limitations on Reprocessing	<ul style="list-style-type: none"> Repeated reprocessing has minimum effect on restorations. End of life is normally determined by clinical factors. Restorations are typically processed once prior to patient delivery. Steam autoclaving, automated washer/disinfection, chemical vapor sterilization, and dry heat sterilization methods have not been tested or validated for efficacy and are not recommended for use.

Initial Treatment at the Point of Use	<ul style="list-style-type: none"> Do not allow residue or any form of contamination to dry on the device. Use a new, clean pair of examination gloves. Remove excess soil/gross debris by wiping with a disposable cloth/cotton-tipped applicator saturated with 70% isopropyl alcohol (IPA). Remove all visible soil, ensuring fluid penetrates all crevices. Use fresh wipes to rub fluid into the crevices. Discard used wipes. Additional wipes may be used. Transport to processing area and follow validated processing instructions below.
Preparation before cleaning	Polish external surfaces and prepare internal surfaces for cementation prior to processing.
Cleaning and Disinfection: Automated	Do not use automated washer/disinfectors for processing indirectly fabricated composite restorations. Material may become distorted or compromised.
Cleaning and Disinfection: Mechanical	<ul style="list-style-type: none"> The restoration may be concurrently cleaned and disinfected if using the recommended solution 70% isopropyl alcohol (IPA). Using a cloth and a cotton-tipped applicator saturated with 70% IPA, wipe all the surfaces of the restoration, ensuring to wipe all grooves and edges. Immerse the restoration in an ultrasonic bath in 70% IPA solution for at least 1 minutes unless otherwise specified by the solution manufacturer. Rinse under tap water of at least drinking quality for a minimum of 30 seconds. Allow to air dry until completely dry. Visually inspect for visible soil and repeat the cleaning and disinfection procedure if necessary.
Cleaning and Intermediate-Level Disinfection: Manual	<ul style="list-style-type: none"> The restoration may be concurrently cleaned and disinfected if using the recommended solution 70% isopropyl alcohol (IPA). Using a cloth and a cotton-tipped applicator saturated with 70% IPA, wipe all the surfaces of the restoration, ensuring to wipe all grooves and edges. Immerse the restoration in 70% IPA solution for at least 1 minutes unless otherwise specified by the solution manufacturer. Rinse under tap water of at least drinking quality for a minimum of 30 seconds. Allow to air dry until completely dry. Visually inspect for visible soil and repeat the cleaning and disinfection procedure if necessary.
Drying	Wipe the devices dry with a clean, lint-free cloth. Allow the components to fully air dry before delivery.
Maintenance, Inspection and Testing	Immediately before delivery, visually inspect to ensure that all contamination has been removed. If the device is discolored, damaged, or distorted it should be discarded.
Packaging	No particular requirements.
Sterilization	<ul style="list-style-type: none"> Sterilization is not allowed. No methods have been validated. Do not subject restorations to steam autoclaving or liquid chemical sterilant immersion. Restoration damage may occur.
Storage	<ul style="list-style-type: none"> Store processed restorations at room temperature, away from moisture or excessive humidity ou d'une humidité excessive. To prevent contamination, store processed restorations in covered storage such as a drawer or cabinet until use.
Additional Information	The instructions provided above have been validated by the manufacturer of the medical device as being capable of preparing a medical device for use. It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials, and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process.
Manufacturer Contact	Within the United States, call Dentsply Sirona at 1-302-422-4511. For areas outside the United States, contact your local Dentsply Sirona representative.

4.2 TruMatch® Shade Guide

Processing Instructions for the TruMatch® Shade Guide, Shade Tabs, and Holder	
Warnings	<ul style="list-style-type: none"> High-level disinfection procedures have not been validated. Intermediate-level disinfection is appropriate for the Holder. For individual Shade Tabs, proceed to Sterilization following cleaning and any optional/intermediate disinfection process. Steam autoclave sterilization is appropriate and recommended for the individual shade tabs. Do not autoclave the Holder in steam autoclave. Do not clean/disinfect the Holder automated washer/disinfector. Follow proper infection prevention activities, such as proper hand washing and donning new puncture and chemical resistant gloves at pertinent steps. Do not use hard wire brushes to clean the device as scratching and mechanical degradation can occur. Use only a disinfecting solution which is approved for its efficacy, EPA registered (and/or Health Canada approved) and use in accordance with the IFU of the disinfecting solution manufacturer. Do not allow the shade tabs to exceed 137°C. Do not use phenol-based glutaraldehyde. Always use a pH neutral instrument cleaning solution if authorized. If the integrity of the packaging pouch has been compromised, the device should be processed again before use. The shade tabs should be cleaned and sterilized before each use. For the countries that require a three-step process, in which disinfection is required prior to sterilization, the shade tabs must be cleaned, disinfected, and sterilized before use.
Limitations on Reprocessing	<ul style="list-style-type: none"> Repeated reprocessing has minimal effect on these instruments. End of life is normally determined by wear and damage due to use. Use of the device beyond its useful life may cause damage to equipment and increases the risk of patient cross-contamination. Cold liquid sterilization, chemical vapor sterilization, and dry heat sterilization methods have not been tested or validated for efficacy and are not recommended for use.

Initial Treatment at the Point of Use	<ul style="list-style-type: none"> Immediately after use, remove exposed individual shade tabs from the exposed holder for processing. Do not allow residue or any form of contamination to dry on the device. Remove excess soil/gross debris by wiping, brushing, and/or rinsing under water. Transport to processing area and follow validated processing instructions below. It is recommended that the device be processed as soon as is reasonably practical following use. Start reprocessing within 1 hour after use.
Preparation before cleaning	Always disassemble the individual tabs from the holder before processing.
Cleaning: Mechanical	<ul style="list-style-type: none"> The chemical manufacturer's instructions for the cleaning and disinfection solution should be followed, observing concentration rates, contact times, fill level, degassing processes, etc. The recommended solution is Resurge® Instrument Cleaning Solution mixed to 0.5 oz/gal (14.8 ml/3.8 l), at 35 °C-40 °C/95 °F-104 °F or equivalent enzymatic cleaning solution. Immerse tabs and holder in an ultrasonic bath containing the prepared solution for at least 20 minutes unless otherwise specified by the time specified by the solution manufacturer. Rinse under tap water of at least drinking quality for a minimum of 30 seconds. Allow to air dry. Visually inspect for visible soil, and repeat the cleaning and disinfection procedure if necessary.
Cleaning: Manual	<ul style="list-style-type: none"> Prepare instrument cleaning solution per the manufacturer's recommendation. The recommended solution is Resurge® Instrument Cleaning Solution mixed to 0.5 oz/gal, at 35-40°C Using the cleaning solution per the label directions, brush with a soft nylon brush, paying particular attention to seams and crevices. Rinse under tepid tap water of at least drinking quality (35 °C-40 °C) for a minimum of 30 seconds. Allow to air dry. Visually inspect for visible soil, and repeat the cleaning procedure if necessary. Note: European markets may use FD 370 Cleaner, Dürr Dental AG, Germany.
Disinfection: Manual	<p>Disinfection by Wiping</p> <ul style="list-style-type: none"> To disinfect the device, thoroughly wipe all surfaces, with a new single-use wipe in combination with an alcohol-based, tuberculocidal, quaternary ammonium solution (41.58% IPA, 0.24% ammonium compounds). Use additional wipes as needed to ensure the entire outer surface, including seams and crevices, remain visibly wet for the contact time recommended by the solution manufacturer. Pay special attention to all seams and crevices while wiping. Rinse the delivery device under tepid running water (40 °C, potable) for 30 seconds to remove disinfectant solution. Dry the device using a dry, lint-free disposable cloth or disposable wipe.
Cleaning and disinfection: Automated	<ul style="list-style-type: none"> Note: The following instructions apply only to shade tabs; Do not use a thermal disinfectant for the holder. Holder may become distorted or compromised. Note: Remove the individual shade tabs for automated cleaning and disinfection. A thermal disinfectant that complies with ISO 15883 may be used for automatic cleaning and disinfection. Use only properly maintained, calibrated, and approved washer disinfectors according to ISO 15883. Follow the manufacturer's instructions to select the solutions and programs. Remove excess soil with a disposable cloth/paper wipe. In case of gross contamination of shade tabs, remove it using a soft brush and running water. Place shade tabs in the washer-disinfectant allowing water and detergent to enter into and drain out of device openings. Run washer-disinfection program with A0 value ≥ 3000 (e.g. 5 min at ≥ 90 °C/194 °F) (e.g. Miele Vario TD) using appropriate detergents (e.g. neodisher® MediClean and neodisher® Z, both from Dr. Weigert, Germany or equivalent). After the cycle is complete, visually inspect the shade tabs for degradation due to extended use and reprocessing. If the shade tabs are discolored, cracked, worn, distorted, etc., discard and do not use. Proceed to Sterilization following Automated Cleaning and Disinfection.
Drying	The drying instructions are incorporated in the cleaning and/or disinfection section above.
Maintenance, Inspection and Testing	Visually inspect the shade tabs and holder for degradation due to extended use and reprocessing. If the device is discolored, cracked, worn, distorted, etc., discard and do not use.
Packaging	<ul style="list-style-type: none"> Ensure that the shade tab(s) are completely dry before packaging for sterilization. A lint-free, disposable cloth may be used to dry the shade tabs. Uses of FDA-cleared, ISO 11607 compliant paper/plastic steam sterilization pouches (e.g., AssurePlus® Sterilization Pouches) are recommended.
Sterilization/ Steam autoclaving	<ul style="list-style-type: none"> The shade tabs must be sterilized utilizing steam sterilization. The following Pre-vac sterilization cycles may be used for shade tabs: <ul style="list-style-type: none"> 132 °C (270 °F) for 4 minutes with a minimum drying time of 20 minutes. 134 °C (273 °F) for 3 minutes with a minimum drying time of 18 minutes (recommended only for non-US markets). 135 °C (275 °F) for 3 minutes with a minimum drying time of 16 minutes. The following Gravity Displacement cycles may be used: <ul style="list-style-type: none"> 132 °C (270 °F) for 15 minutes with a minimum drying time of 30 minutes. 135 °C (275 °F) for 10 minutes with a minimum drying time of 30 minutes. <p>An alternate method of sterilization is to place non-bagged instruments into the steam sterilizer and run at one of the above listed cycles</p>

Storage	<ul style="list-style-type: none"> Do not remove from the pouch until ready for use to prevent contamination. Instruments sterilized non-bagged should be used immediately. Store at room temperature, away from moisture or excessive humidity. To prevent contamination, store processed device in covered storage such as a drawer or cabinet until use. Prior to use, inspect the pouch. If the integrity of the pouch has been compromised, the device should be processed again before use.
Additional Information	<ul style="list-style-type: none"> Prior to use, inspect the device. Discard any device which has become damaged, worn, or distorted. Reassemble for use (replace individual shade tabs in holder) as desired. <p>The instructions provided above have been validated by the manufacturer of the medical device as being capable of preparing a medical device for reuse. It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials, and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process.</p>
Manufacturer Contact	Within the United States, call Dentsply Sirona at 1-302-422-4511. For areas outside the United States, contact your local Dentsply Sirona representative.

4.3 Compules® Tips Gun

For complete processing instructions and validated parameters, please refer to the Instructions for Use of the Compules® Tips Gun, which is available on our webpage at <https://www.dentsplysirona.com/ifu>. If requested, we will send you a free printed copy of the Instructions for Use in the language you require within 7 days. In the USA, call 1-800-532-2855. Outside the USA, use the order form provided on <https://www.dentsplysirona.com/en/service-contact/ifu.html#orderform-ifu> for this purpose.

4.4 Disposal

- Based on the information in the SDS, the waste generated by the disposal of the product does not present characteristics of physical danger, for the person or the environment; It can therefore be managed, according to local regulations, as a non-hazardous waste.
- Contaminated containers and accessories that have been in contact with the patient must be cleaned and disinfected as outlined in Hygiene above prior to disposal or disposed of as hazardous waste at risk of biological contamination.
- Always wear gloves when handling packaging and accessories.

5. LOT NUMBER, EXPIRATION DATE AND CORRESPONDENCE

- Do not use after expiration date. Do not use if the lot number and/or expiration date are missing or illegible. ISO standard uses: "YYYY-MM-DD."
- The following numbers should be quoted in all correspondences:
 - Reorder number
 - Lot number
 - Expiration date
- Any serious incident in relation to the product should be reported to the manufacturer and the competent authority according to local regulations.
- A summary of the safety and clinical performance (SSCP) for this product can be found (upon activation) at <https://ec.europa.eu/tools/eudamed> by searching for the Basic UDI-DI number: ++D002DENTALCOMPOSITES68

© 2025 Dentsply Sirona. All Rights Reserved. **530600WEB** (R 1/13/25)

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