

1 Identification of the substance/preparation and of the company/undertaking

- **Product details**
- **Trade name:** CARIES DETECTOR
- **Article number:** 043
- **Application of the substance / the preparation** Caries detection device
- **Manufacturer:**
Kuraray Europe GmbH
Building F821, Hoechst Industrial Park
65926 Frankfurt am Main
Germany
Tel: +49 (0)69 305 35 840
Fax: +49 (0)69 305 35 640
E-mail: dental@kuraray.de
- **Supplier:**
J & S Davis
Unit 6, Argyle Way Trading Estate,
Fulton Close, Stevenage,
Hertfordshire, SG 1 2AF,
Great Britain
Tel: +44- (0)1438 735800
Fax: +44- (0)1438 735801
E-mail: jsdsales@js-davis.co.uk
- **Further information obtainable from:** Manufacturer or Supplier
- **Information in case of emergency:**
National Poisons information Service (London Centre)
Guy's & St. Thomas' Hospital Trust
Medical Toxicology Unit
Avonley Road
LONDON
SE14 5ER
Emergency telephone: (local) 0870 600 6266 or +44 20 7771 5394
Administrative telephone: +44 20 7 771 53 10
Facsimile number: +44 20 7 771 53 09

2 Hazards identification

- **Hazard description:** Not applicable.
- **Information concerning particular hazards for human and environment:**
The product does not have to be labelled due to the calculation procedure of the "General Classification guideline for preparations of the EU" in the latest valid version.
- **Classification system:**
The classification is according to the latest editions of the EU-lists, and extended by company and literature data.
- **GHS label elements**
Safety data sheet available on request.

3 Composition/information on ingredients

- **Chemical characterization**
- **Description:** Mixture of substances listed below with nonhazardous additions.

- **Dangerous components:**

CAS: 57-55-6 EINECS: 200-338-0	propane-1,2-diol	> 90%
-----------------------------------	------------------	-------

(Contd. on page 2)

Printing date 27.08.2009

Revision: 24.08.2009

Trade name: CARIES DETECTOR

· **Other ingredients:** Dyes

(Contd. of page 1)

4 First aid measures

- **After inhalation:** Supply fresh air; consult doctor in case of complaints.
- **After skin contact:** Generally the product does not irritate the skin.
- **After eye contact:** Rinse opened eye for several minutes under running water.
- **After swallowing:** If symptoms persist consult doctor.

5 Fire-fighting measures

- **Suitable extinguishing agents:** Use fire extinguishing methods suitable to surrounding conditions.
- **Protective equipment:** No special measures required.

6 Accidental release measures

- **Person-related safety precautions:** Not required.
- **Measures for environmental protection:**
Dilute with plenty of water.
Do not allow to enter sewers/ surface or ground water.
- **Measures for cleaning/collecting:**
Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).
- **Additional information:** No dangerous substances are released.

7 Handling and storage

- **Handling:**
- **Information for safe handling:** No special measures required.
- **Information about fire - and explosion protection:** No special measures required.
- **Storage:**
- **Requirements to be met by storerooms and receptacles:** Store in a cool location (2-25°C).
- **Information about storage in one common storage facility:** Not required.
- **Further information about storage conditions:** None.

8 Exposure controls/personal protection

- **Additional information about design of technical facilities:** No further data; see item 7.

· **Ingredients with limit values that require monitoring at the workplace:**

57-55-6 propane-1,2-diol

WEL	Long-term value: 474* 10** mg/m ³ , 150* ppm *total vapour and particulates **particulates
-----	--

- **Additional information:** The lists valid during the making were used as basis.

· **Personal protective equipment:**

· **General protective and hygienic measures:**

The usual precautionary measures are to be adhered to when handling chemicals.

· **Respiratory protection:** Not required.

· **Protection of hands:**

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.

Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.

(Contd. on page 3)

Printing date 27.08.2009

Revision: 24.08.2009

Trade name: CARIES DETECTOR

(Contd. of page 2)

Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation

- **Material of gloves**

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

- **Penetration time of glove material**

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

- **Eye protection:** Goggles recommended during refilling

9 Physical and chemical properties

- **General Information**

Form:	Liquid
Colour:	Red
Odour:	Odourless

- **Change in condition**

Melting point/Melting range:	-68 °C (propane-1,2-diol)
Boiling point/Boiling range:	185 °C (propane-1,2-diol)

- **Flash point:**

101 °C (propane-1,2-diol)

- **Ignition temperature:**

371.0 °C (propane-1,2-diol)

- **Danger of explosion:**

Product does not present an explosion hazard.

- **Explosion limits:**

Lower:	2.6 Vol % (propane-1,2-diol)
Upper:	12.6 Vol % (propane-1,2-diol)

- **Vapour pressure:**

0.1 hPa (20°C, propane-1,2-diol)

- **Density:**

Not determined.

- **Vapour density**

Not applicable

- **Solubility in / Miscibility with water:**

Fully miscible.

- **pH-value:**

5.6

- **Segregation coefficient (n-octanol/water):** Not applicable

10 Stability and reactivity

- **Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.
- **Dangerous reactions** No dangerous reactions known.
- **Dangerous decomposition products:** Carbon monoxide and carbon dioxide

11 Toxicological information

- **Acute toxicity:**
- **Primary irritant effect:**
- **on the skin:** No irritant effect.
- **on the eye:** No irritating effect.
- **Sensitization:** No sensitizing effects known.

(Contd. on page 4)

Printing date 27.08.2009

Revision: 24.08.2009

Trade name: CARIES DETECTOR

(Contd. of page 3)

· Additional toxicological information:

The product is not subject to classification according to the calculation method of the General EU Classification Guidelines for Preparations as issued in the latest version.
When used and handled according to specifications, the product does not have any harmful effects to our experience and the information provided to us.

12 Ecological information**· General notes:**

Water hazard class 1 (German Regulation) (Self-assessment): slightly hazardous for water
Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.

13 Disposal considerations**· Product:**

· Recommendation Smaller quantities can be disposed of with household waste.

· Uncleaned packaging:

· Recommendation: Disposal must be made according to official regulations.

· Recommended cleansing agents: Water, if necessary together with cleansing agents.

14 Transport information**· Land transport ADR/RID (cross-border)**

· ADR/RID class: -

· Maritime transport IMDG:

· IMDG Class: -

· Marine pollutant: No

· Air transport ICAO-TI and IATA-DGR:

· ICAO/IATA Class: -

· UN "Model Regulation": -

15 Regulatory information**· Labelling according to EU guidelines:**

Observe the general safety regulations when handling chemicals.

The product is not subject to identification regulations under EU Directives and the Ordinance on Hazardous Materials (German GefStoffV).

· Special labelling of certain preparations:

Safety data sheet available for professional user on request.

16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

· Version: 3**· Abbreviations and acronyms:**

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)

RID: Règlement international concernant le transport des marchandises dangereuses par chemin de fer (Regulations Concerning the International Transport of Dangerous Goods by Rail)

IMDG: International Maritime Code for Dangerous Goods

(Contd. on page 5)

Safety Data Sheet
according to 1907/2006/EC, Article 31

Printing date 27.08.2009

Revision: 24.08.2009

Trade name: CARIES DETECTOR

(Contd. of page 4)

IATA: International Air Transport Association

IATA-DGR: Dangerous Goods Regulations by the "International Air Transport Association" (IATA)

ICAO: International Civil Aviation Organization

ICAO-TI: Technical Instructions by the "International Civil Aviation Organization" (ICAO)

GHS: Globally Harmonized System of Classification and Labelling of Chemicals

-GB-