

CONVENTIONAL, ISO COLOUR-CODED,  
ISO TAPER AND MULTI TAPER GUTTA-PERCHA POINTS.



## Application:

- For permanent hermetic sealing of root canals.

## Why choose R&S gutta percha points?

- + Excellent radiodensity**  
Highly contrasted for increased visibility in x-rays.
- + Extremely rigid and highly flexible**  
Rigid so they will not bend but flexible enough to use in curved canals.
- + Production quality**  
All of R&S's points are in compliance with and surpass Standard ISO 6877 and ADA specifications, and they are individually inspected by laser.
- + Uniformity**  
The hand-rolling method allows for:
  - Taper and uniform consistency.
  - Smooth, rounded points.
  - Perfect for preparing canals.



## Advice from our dental experts

- >** To ensure the gutta-percha points will be perfectly clean before inserting them into the canal, I recommend dipping them in sodium hypochlorite 1% solution for one minute and then drying them on a sterile compress.

## Technical specifications and instructions for use

ISO STANDARD	6877 - DENTISTRY - ROOT-CANAL OBTURATING POINTS			
Type	CONVENTIONAL	ISO COLOUR-CODED	ISO TAPER	ISO MULTI-TAPER
Use	With all types of sequences, such as accessory tips or a master tip	With all conventional instruments	With all instruments with constant tapers, such as PROFILE®, CMA®, HERO 642®, HERO SHAPER® and R-endo®	With all instruments with variable tapers, such as PROTaper®, Go Taper Universal & Flex
Shelf Life (in years)			4	
Taper (%)	Variable	2	4 and 6	Variable
ISO (Ø)				
Length (mm)	20	28		30
Sterilisation			Sterilised by gamma rays	
Packaging		Box of 120 points		Box of 60 points

## Additional products:

Dams with or without latex



Sodium hypochlorite



Nerve broaches



Root canal pluggers



Easyaccess Niti



K-files, K-reamers and H-files



Endo EZ burs



Gutta percha points, Sodium hypochlorite, Nerve broaches, Easyaccess Niti, Root canal pluggers and Endo EZ burs: Class IIa medical devices. CE Marking 0120. Certifying body: SGS UK LIMITED.  
 K-files, K-reamers and H-files: class I medical devices. CE marking.  
 Photos are examples only.  
 Manufactured by CFPM.

GACD - 25 rue Bleue 75009 PARIS - France 435 051 776 R.C.S PARIS - SASU.  
 Tel.: +33 (0)1 48 01 27 14 - Fax: +33 (0)1 45 23 51 89.  
 export@rs-dental-products.com - www.rs-dental-products.com.

