

CFPM

www.rs-dental-products.com

### Description

**R&S Diamond Burs** are available in a wide selection of shapes and sizes, and in five color-coded grits: s1 per coarse, coarse, medium, fine, and super fine. The burs are made by the application of a homogeneous, plated diamond surface to a one-piece stainless steel blank. They are available in packages of five in a plastic box with a clear sliding top for easy identification and selection. The burs are guaranteed for 3-5 crown preparations. Three shapes (football, wide round-end taper, narrow round-end taper) in a coarse grit of the **R&S Diamond Burs** were evaluated by 29 consultants in over 380 uses.

This product received a **91% clinical rating**.

**91%**  
overall  
rating




**R&S**

### Consultants' Comments

- "Good bulk reduction with the coarse grit."
- "Can prepare several teeth with each bur."
- "Nice assortment of shapes."
- "The catalog is well laid out for easy ordering."
- "I'd like to have these in a short-shank version."

"THE BUR WORKS WELL IN MY ELECTRIC HANDPIECE."

Percentage of Consultants Who Would:



### Product Features

Consultants consistently gave **R&S Diamond Burs** high marks for packaging, fit to handpiece, shank length, concentricity, and resistance to clogging. The cutting efficiency and durability of the burs were features that consultants appreciated. Several noted that the burs lasted longer and cut better than other burs they had used, making them very cost-effective. Eighty-two percent of consultants would switch to **R&S Diamond Burs** and 90% would recommend them to colleagues.