

W&H Sterilizers

Cycle Type(s)



How to choose the right class of sterilizer!

It is important to ensure that the chosen cycle type is capable of correctly processing the intended load.



N

Type – non-vacuum sterilization suitable for ‘Naked’, solid, unwrapped instruments - not suitable for hollow or porous load (i.e handpieces)

S

Type – may be vacuum or non-vacuum sterilization, only suitable for load types ‘Specified’ by the manufacturer

B

Type – vacuum sterilization; can process the same load types as a ‘Big, Best or Bag’ hospital sterilizer