

Fast-Fill Fast-Pack pro



— Cordless Obturation System —


Fast-Fill




Useable cartridge tip, economical design
GP pellet refill into device
Easy maintenance



Motorized gutta-percha
Push & stop, better flow controlling



Digital displaying of rest
gutta percha



Real-time temperature displaying
Set from 100~200°C



3 ranges pushing speed, low, mid, high



Powerful battery



| Obturation Fill Handpiece

Fast-Pack pro



| Obturation Pack Handpiece

Fast-Fill Fast-Pack pro

