

MIXPAC™

T-Mixer™ Mixing Tips



A new generation of mixing tips based on the latest state-of-the-art mixing technology. Its compact design with a shorter geometry delivers less material waste with consistent, homogenous mixing results compared to the conventional Helical mixing tips.

Mixing tips for 18 mL, 25 mL, 30 mL, 50 mL and 75 mL Systems.




T-Mixer, 6.5 mm

- Impression taking (tray material)



Outlet	Ratio	Description	Inner Ø (mm)	Elements	Length (mm)	Waste Vol. (mL)	Order No.	Mixing Tips per bag
D	1:1, 2:1	teal, natural elements	6.5	12	67.6	1.66	186992	50




T-Mixer, 5.4 mm

- Impression taking
- Bite registrations



Outlet	Ratio	Description	Inner Ø (mm)	Elements	Length (mm)	Waste Vol. (mL)	Order No.	Mixing Tips per bag
D	1:1, 2:1	pink, natural elements	5.4	12	59.4	1.06	186988	50




T-Mixer, 4.2 mm

- Impression taking (wash material)
- Core build-ups
- Relining




Outlet	Ratio	Description	Inner Ø (mm)	Elements	Length (mm)	Waste Vol. (mL)	Order No.	Mixing Tips per bag
D	1:1, 2:1	yellow, natural elements	4.2	12	50.1	0.63	186991	50




Intra Oral Tips

for mixing tips with D outlet

For intraoral applications in confined areas to ensure for precise and safe dispensing of materials.



Ratio	Description	Order No.	Tips per bag
1:1, 2:1	yellow, for 4.2 mm diameter mixing tips	186985	100
	natural, for 5.4 mm and 6.5 mm diameter mixing tips	186986	100




T-Mixer™, 3.2 mm

- Temporary crown & bridge



Outlet	Ratio	Description	Inner Ø (mm)	Elements	Length (mm)	Waste Vol. (mL)	Order No.	Mixing Tips per bag
S	1:1, 2:1	light blue, natural elements	3.2	12	49.8	0.37	186989	50




T-Mixer, 3.2 mm

- Temporary crown & bridge



Outlet	Ratio	Description	Inner Ø (mm)	Elements	Length (mm)	Waste Vol. (mL)	Order No.	Mixing Tips per bag
S	4:1, 10:1	light blue, orange elements	3.2	16	57.6	0.46	186990	50



T-Mixer Colibri™ Plus Breakable


- Impression of root canal
- Impression taking (wash material)
- Core build-ups
- Cementation of root post



Ratio	Description	Outer Ø (mm)	Inner Ø (mm)	Gauge	Elements	Ø (mm)	Length (mm)*	Waste Vol. (mL)	Order No.	Mixing Tips per bag
1:1, 2:1	mint, natural elements	1.1	0.8	19	12	4.2	63.8+15	0.66	186987	25

* length mixing tip + length cannula

All of the MIXPAC Candy Colors™ are available for mixing tips in any diameter, length, type or application upon request. The Candy Colors are Trademarks representing genuine Swiss quality.



medmix Switzerland AG | Ruetistrasse 7 | 9469 Haag, Switzerland
P +41 81 414 70 00 | www.medmix.swiss

Legal Notice: This brochure is a general presentation. The information contained in this publication is believed to be accurate and reliable but is not to be construed as implying any warranty or guarantee of performance. We waive any liability and indemnity for effects resulting from its application. Please, contact us for a description of the warranties and guarantees offered with our products. Directions for use and safety will be given separately. All information herein is subject to change without notice. MIXPAC, T-Mixer, and Colibri are trademarks of medmix Switzerland AG. medmix is a trademark of medmix AG. © medmix Switzerland AG 2024



DENTAL