



Toffeln®

**FOOTWEAR & SURGICAL
PRODUCTS CATALOGUE 2023/2024**

HEALTHCARE INNOVATION FOR YOUR WELLBEING

For your convenience, we have put all Toffeln products into one catalogue for you.

We are known and loved for our footwear, and we are currently the number one comfort brand with an extensive footwear range designed specifically for healthcare professionals. Our footwear products are presently being worn by around 90% of UK operating theatre staff and by thousands of healthcare professionals in over 18 countries worldwide.

Did you know that Toffeln also offers an extensive range of surgical products designed for your well-being and, like our footwear, are all guaranteed certified? Better still, many products are available through NHS supply chain contracts and are also used by most Private Healthcare providers. We back all Toffeln products with stock in our Bristol warehouse for fast delivery, and after 45 years, we know and understand the pressures on healthcare, and we are always here to support you.



Charlie Leflaive
Director



Our price promise

We want everyone to be able to wear Toffeln shoes, so we work hard to deliver at the most reasonable price possible. If you find equivalent footwear at a lower price, we will better that price (reasonable proof required).

Our guarantee

True to our mission, we guarantee the quality of every pair of shoes that we dispatch. If they fail to reach the high quality you would expect, we will replace them free of charge. And for further assurance, we are an ISO 9001 quality registered company.

3 Steps to easy ordering

- 1 Choose your product
- 2 Use the product code, followed by the size required
- 3 Email your order through to sales@toffeln.com

Alternatively, we are **PEPPOL ready** for fully integrated ordering, so please get in touch with us to arrange integration.

Need help?

Get in touch **0345 500 4433** sales@toffeln.com
Check out even more choices toffeln.com

The Toffeln Difference

1-6

Footwear

7-66

Surgical products

67-106

Key product features

 CE Marked	 Shock Absorbing	 UKCA
 Slip Resistant Sole	 E-tech Sole Comfort	 Anti-Static Sole
 Leather Upper	 Wash Temperature	 12 Month Guarantee
 Fully Autoclavable	 Safety Toe-Cap	 Action Breathe Upper
 Recycled	 Latex Free	 Fluid protection
 LiquidTex	 Sterile	 Single use
 OEKO-TEX	 BarrierTex	 EcoTex
 ENI3795	 SoftTex	

TOFFELN POINTS OF DIFFERENCE

Personal service from people who care

At Toffeln, we take pride in providing the latest healthcare products and a personal service that stands out from the crowd from a team who genuinely cares about healthcare and your well-being.

Our highly trained customer care professionals are on-hand to go the extra mile and assist you with any questions you may have about our diverse range of footwear and surgical products.

We also understand the importance of offering hands-on experience with our products, which is why our dedicated Key Account Managers can provide a product demonstration, a free trial or a product sample.

Reliable & Trustworthy

We understand that experience matters, and as a friendly family-run business with over 45 years in the industry, we have gained extensive knowledge that our customers have come to trust and rely on.

Our expertise enables us to design footwear and surgical products that are explicitly tailored to meet the healthcare industry's unique demands. We are proud to say that all products within this catalogue exceed the quality control standards you'd expect and deserve.

Additionally, our long-standing relationship with the NHS means that most of our products are already available on NHS Supply Chain for your convenience and peace of mind.



Look out for the logo and see how many products are already available on NHS Supply Chain.

Meet our friendly Customer Care Team



Prioritising Efficiency & Sustainability

At Toffeln, we understand the importance of offering a reliable delivery service that is efficient and sustainable. We take great pride in our extensive stock holding facility in Bristol, UK, which is manned by our experienced and capable team of Warehouse Operatives. This facility enables us to deliver the right products to the right places on a next day service.

Our distribution partner has an impressive on-time delivery rate of 98% and shares our carbon reduction goals by committing to a 100% electric vehicle fleet by 2030. Our surgical gowns and drapes manufacturing facility also produces 35% of their electricity from solar energy. Better still, their proximity to the UK dramatically reduces our carbon footprint and lead times, allowing us to respond promptly to any increase in demand.

Meet our experienced Warehouse Operations Team



Customer Service Promise



We are here to help
8am-5pm Monday – Friday



Order by **3pm** for **same day despatch**



All Toffeln products undergo rigorous testing and carry a **full quality guarantee**.

“Our extensive stock holding facility and our experienced team, enables us to deliver the right products to the right places on a next day service.”

Corrie Leflaive
Operations Director



TOFFELN GIVING BACK

At Toffeln we believe in giving back. We support a wide range of Charities engaged in humanitarian work. In fact, we have made large donations to charitable causes. In addition, our New Lease of Light programme donates reconditioned headlights to worthy causes around the world along with footwear and medical equipment.



Toffeln is dedicated to supporting charities including:



TOFFELN SUPPORTING THE ROYAL COLLEGE OF SURGEONS OF EDINBURGH

We are pleased to be supporting the Royal College of Surgeons of Edinburgh to launch a new fellowship that will support trainee surgeons to gain international experience, share knowledge and best practice in surgical care. By sharing expertise and bringing back new skills, the fellows will contribute to the continuous improvement of surgical care in the NHS and UK Healthcare, ultimately leading to better patient outcomes.



“As a company committed to improving healthcare, we believe that supporting our trainee surgeons will make a real difference. We are thrilled about the possibilities that our sponsorship will bring to the surgical community and are excited to welcome the first RCSEd Toffeln Travelling Fellows, who will undoubtedly make a significant contribution to the NHS and UK Healthcare”

Ben Turner
Managing Director



SUSTAINABLE COMMITMENTS

At Toffeln, we are passionate about reducing our carbon footprint and working towards a more sustainable future for us, our customers, and the planet. As a responsible company, we care about the future and believe we must support the NHS Net Zero Plus program and reduce our influenceable emissions to help combat climate change.

Step to the Future

As a company, we are committed to becoming **carbon neutral by 2026**, ahead of both the NHS 2045 Net Zero date and the 2036-2039 80% reduction target date. To achieve this goal, we have established a carbon savings committee with members from every department within Toffeln. We aim to achieve a 20% reduction in carbon emissions each year until we hit carbon-neutral status in 2026.

CARBON
NEUTRAL



Key Initiatives

We are excited to share some of our key sustainability initiatives, that will help us on our journey to becoming carbon neutral:

- 1 Developing new space-saving packaging solutions from recycled FSC rated material.
- 2 Transitioning towards producing products from sustainable sources.
- 3 Reducing single-use plastic consumption across all our product ranges.
- 4 Switching to LED lighting to reduce carbon emissions by 55%.
- 5 Introducing PIR (Passive Infra-Red) detection sensors to reduce energy consumption from lighting when areas are not in use.
- 6 Using electric vehicles to deliver our products where possible.
- 7 Installing electric vehicle charging points at our work premises to encourage our staff to invest in sustainable methods of transport.
- 8 By selecting sustainable manufacturing facilities who demonstrate green credentials and who are in closer proximity to the UK.

We believe that these initiatives will help us reduce our carbon footprint and make a positive contribution to a more sustainable future for us, the planet and its future.



We have implemented sustainable practices, including the installation of electric vehicle charging points for staff, the use of electric vehicles for UK distribution, and by selecting sustainable manufacturing facilities – Our Surgical Gowns & Drapes Manufacturer generates 35% of their electricity from solar panels, and they are looking to increase this number to 50% in 2024.



We are developing new ways to manufacture our products sustainably and already offer hospital clothing made from **75% recycled material** (see pages 80-81) and have recently developed a **Carbon Neutral Clog** (see page 19-20). We want to extend this range in 2024 by creating a carbon-neutral surgical gown.

Introducing Toffeln Green T

We understand that the journey towards offering sustainable products is a long one. This is why in 2022, we took our first steps on this journey and developed a sub-brand of Toffeln that embodies our sustainable values by offering a collection of products that have a reduced carbon footprint and are manufactured with Renewable and/or Recycled sources. We called this sub-brand the Green T Collection.



Discover the first product of our **Green T Collection - Green T EziKlog**

See page 19 for more details



“Through the changes implemented so far, we have successfully reduced our Electricity Carbon Footprint by 40%, which shows encouraging progress”

Jim Leflaive Systems Manager





FOOTWEAR

At Toffeln we've dedicated our lives to enhancing foot wellbeing, helping healthcare professionals to have healthier, supremely well supported feet.

We know how harshly poorly designed, ill-fitting and inappropriate footwear can impact tender toes and hard-pressed arches. That's why we work with the leading biomechanics researchers to design and develop the most advanced range of fully compliant footwear.

We combine the functional and the fashionable to create ergonomically designed, super lightweight, beautifully fitting footwear, to provide ultimate comfort for hour after hour of standing, lifting and walking.

And driven by our mission to make every minute of those long working shifts seem shorter, at Toffeln we're passionate about...

Loving your feet

Footwear Products



9-12

Anthony Stoker

Meet Toffeln'
Footwear Designer



13-18

Science of comfort



19-20

Green T EziKlog

The Worlds first
carbon neutral clog



21-24

EziKlog V2.0



25-34

SmartSole Collection



35-38

Aktiv Collection



39-40

FlexiKlog



41-44

UltraLite



45

ClassicKlog



46

ErgoBoot



47

KlimaFlex



48

SteriKlog



49-58

Safety Footwear



59-60

Accessories



61-62

AktivCool Socks



63-64

Footwear Index

Meet Toffeln's Footwear Designer

At Toffeln, we are passionate about delivering an exceptional customer experience. And, even though comfort and support are our top priority, we also understand that design plays a pivotal role in shaping this experience.

That's why, in 2015, we started collaborating with the renowned footwear designer Anthony Stoker. Since then, we have nurtured a great partnership with Anthony, and together, we have embarked on some awe-inspiring projects that have led to the evolution of some of our best selling collections - from our famous SmartSole Range to this year's latest addition, the Green T EziKlog.

In this section, we offer a glimpse into the creative mind of Anthony as he shares sources of inspiration and gives you an insight into the design process, from his initial sketches to the final product.



“When I first began working with Toffeln a few years ago, I was amazed by the fantastic position they were in, in being able to serve such a unique and caring audience. It’s been an honour to design shoes that offer the comfort and support you require to carry out your role so that you can provide your patients with the comfort and support they need.”

Anthony Stoker
Footwear Designer



Developing SmartSole

At Toffeln, we understand the unique position that our healthcare professionals are in and that their specific needs for workplace footwear would not be met with just one solution. With this in mind, we embarked on the remarkable journey of SmartSole, with a vision to develop a range of footwear that prioritised proper standing posture and incorporated interchangeable insoles that catered to the specific contours of each individual's footbed.

The result? A revolutionary solution that dramatically reduces the risk of musculoskeletal problems in healthcare workers by delivering unparalleled comfort and customisable support tailored to the individuals' preferences.

Anthony's Inspiration

The seeds of the **SmartSole** range began with wanting to develop a style - a language of sorts - that people working within healthcare would come to recognise as something designed specifically for them - as something of their own. A range that they would feel confident in, supported by, proud to wear and to recognise them on their colleagues' feet and feel an affinity and a sense of trust. A visual language that would be recognised as comfortable and supportive. Not just physically but psychologically too. It's about building that trust and honouring it.

I've delivered comfort levels for other brands that have seen significant success. Together with the Footwear Development Team at Toffeln, we've provided a footwear range that we can be proud to see across the healthcare industry.

SmartSole Shoe Development



SmartSole Collection

See page 25 for more details

Developing Green T EziKlog

Our dedication to sustainability led us to create footwear that prioritises comfort and functionality for the end-users and respects the environment during the manufacturing process. As a starting point for this project, we took our best selling **EziKlog** design and looked into ways to reduce our carbon emissions during manufacturing.

After two years of diligent research and development, we proudly present our latest addition to the Toffeln footwear range and the **World's First Carbon Neutral Clog**. We have achieved a remarkable milestone in our commitment to a greener future by utilising sustainable materials from renewable sources - in the case of **Green T EziKlog**, we replaced part of the traditional EVA content with Sugar Cane Waste and used recycled rubber to create a slip-resistant sole. With this innovation, healthcare professionals can now step forward with confidence, knowing that their choice of footwear aligns with their eco-conscious values.

Anthony's Inspiration

Sometimes it's about evolving what we have, to be better, rather than starting anew. And when it came to developing a new sustainable **EziKlog**, 'better' started with looking at nature's sustainable substitutes to traditional EVA.

We found that by incorporating sugar cane, we could reduced our usual EVA compound usage by half. Pairing this with a recycled rubber sole then elevated our **EziKlog** offering to new levels by attaining a carbon-neutral status.

Better still, the reference to sugar cane also provides a tasty inspiration for some **fruity-coloured clogs!**

Green T EziKlog Development



EziKlog Green T

See page 19-20 for more details

The Science of Comfort & Wellbeing

The Difference with Toffeln

There are many brands of footwear promoting comfort. The difference with Toffeln footwear is that it is designed specially for the healthcare industry, for doctors and nurses... people like you.

Better by Design

All shoes are not equal - sports footwear is designed for intensive, limited duration activities, the popular high street brands are not designed for healthcare environments or the length of shift for which they are worn; hard working feet need supportive footwear, so soft is not always the best solution. That is why Toffeln footwear has gone through extensive testing in healthcare environments in order to create a unique product that really provides the support and all day comfort that healthcare professionals deserve. They are also marked with recognised standards and approvals as proof that they are fit for purpose which you will not find on most other brands.



Toffeln Biomechanics Working for You

We have over 40 years of expertise in healthcare footwear. Toffeln work closely with industry experts and scientists in foot health and have commissioned the first in depth studies on the impact of operating theatre work on feet and the musculoskeletal system.

Working alongside experts at the University of Salford we have carried out biomechanical studies to understand the effects of long hours of standing and walking in Operating Theatres, Hospitals, and related Healthcare environments.



University of
Salford
MANCHESTER




Toffeln Footwear Manufactured with Love

Every pair of Toffeln shoes is manufactured with painstaking care to ensure they give you lasting comfort and reliability.

To ensure you enjoy the best service we carry stock of more than 70,000 pairs for fast despatch to you.

Easy Tips for Footcare

Get into the habit of carrying out this simple footcare routine for 21 days and feel the difference...




3 Cut toenails regularly. Always cut straight across and never at an angle or down the edges, which can cause ingrown toenails



4 Wear good socks with a high natural fibre content, or technical smart fibres and change daily. Compression socks can help to reduce fluid in legs and feet and reduce fatigue



5 Dry thoroughly, paying particular attention to drying between the toes. This will help to prevent fungal infections



6 Remove hard skin. Avoid the build-up of thick and painful hard skin by using a foot file or pumice stone to gently get rid of any dead skin cells, especially on the heels



7 Select well-fitting footwear designed specifically for healthcare use. These usually have a cushioned sole, closed front, are quiet to wear and will be marked SRC for exceptional slip resistance



8 Check feet daily for any blisters, cuts, spots, or swellings and don't just 'put up' with foot pain. Report any serious concerns to a Podiatrist or Doctor



1 As part of your daily routine wash your feet in mild soapy water to gently remove impurities



2 Moisturise with a good urea-based moisturiser - this improves hydration and prevents dry, cracked skin

See our New range of **Bamboo Compression Socks** on page 61-62



Are you looking after your "highly valued teams" feet?

To ensure good levels of staff wellbeing and to improve infection control footwear should be regularly replaced. New footwear provides maximum levels of support for hard working feet and is more hygienic.

"We suggest replacing footwear twice a year, it's a small price to pay for comfort and safety."

Matt Oxbrough
Product Development & Quality Assurance



Plantar Fasciitis

If you've ever experienced a persistent stabbing pain on your heel – especially when you get out of bed in the morning – then you know all about plantar fasciitis.

This common orthopaedic complaint can cause nagging discomfort that makes walking almost unbearable. While many runners battle this condition while exercising, for many that are on their feet all day at work, this can turn your job into a nightmare.

What causes plantar fasciitis?


Plantar fasciitis is caused by straining the part of your foot that connects your heel bone to your toes (plantar fascia). It is not always clear why this happens. You may be more likely to get plantar fasciitis if: you are 40 to 60 years of age, do a lot of walking running or standing, wear shoes with poor cushioning or support or are overweight.

Possible Solutions to the list of ways of easing PF:

- Wearing wide, comfortable shoes with a low and soft heel
- Using insoles in your shoes
- Comfortable arch support
- Resting your foot raised when you can
- Avoid putting pressure on your feet if you're able to (e.g. try swimming rather than running)

“I suffer from plantar fasciitis on my 9 hour+ shifts. With Toffeln SmartSoles, I have comfort throughout my shift and am pain free when I get home.”

Catherine Nurse



What is plantar fasciitis?

Plantar fasciitis is a common cause of heel pain, which involves inflammation to the thick band of tissue which runs across the bottom of your foot. This band of tissue is used to connect your heel bone to your toes, is also known as the plantar fascia. Plantar fasciitis can often cause a sharp pain that normally occurs as soon as you wake up.

The pain commonly decreases through the day, but can return after walking, standing or running for a long period of time. Despite its insidious nature, it's one of the most prevalent orthopaedic ailments, with an estimated 1 in 10 people suffering from it in their lifetime.

What are the best shoes for plantar fasciitis?

The **Toffeln SmartSole Range** receives praise from sufferers of this condition with its cupped insole, contoured footbed and customisable arch support that stops your feet from flattening out and limits pronation. It offers extra support around the mid-foot to prevent impact on the heel and evenly distributes your weight across the whole foot to help ease the pain for those with plantar fasciitis.



Toffeln SmartSole Clog



See page 29 for more details

Toffeln SmartSole Trainer



See page 33 for more details

Do you have staff who suffer from plantar fasciitis, foot pain, or tired and aching feet?

We suggest you try the SmartSole range or one of the Toffeln insoles. Thousands of happy healthcare workers have proved this for themselves.

If you suffer from any other types of foot conditions, please contact us and we'll be glad to advise.

Green T EziKlog Toffeln

MEET THE WORLD'S FIRST CARBON NEUTRAL CLOG

Using the same comfort technology trusted worldwide by countless doctors, nurses and other healthcare professionals, we've taken our best-selling EziKlog and made it carbon neutral.

MADE FROM SUGAR CANE

Using a renewable green EVA compound, the new Green T EziKlog is made from **51% sugar cane waste**.



12 Month Guarantee



RECYCLED RUBBER SOLE

It's **recycled rubber** sole is anti-static, slip-resistant and meets all the regulations required by healthcare professionals.

NEW COMPACT PACKAGING

Our newly designed, packaging solution transports up to **40% less air**. Helping you do your bit with a reduced carbon footprint.



0899GU
Guava
Sizes 3-12



0899BB
Blueberry
Sizes 3-12



Scan me to find out more about the Green T EziKlog by Toffeln



“The Footwear Development Team takes great pride in their accomplishment of attaining carbon neutral status with the Green T EziKlog and consider it a significant step towards a sustainable future for Toffeln and Healthcare.”

Roy Leflaive
Footwear Development



UK	Euro
3	36
4	37
5	38
6	39
6.5	40
7	41
8	42
9	43
10	44
10.5	45
11	46
12	47

Features

- 1 Enhanced Comfort Footbed**
Built-in anti-fatigue technology for day-long comfort.
- 2 Softer Feel E-Tech Material**
Updated E-Tech material to give an improved feel.
- 3 GripSafe Extreme Sole**
Cutting-edge slip-resistant technology helps reduce accidents at work offering up to 50% more slip resistance than its predecessor.
- 4 Improved Anti-Static System**
Shock absorbing anti-static device specifically incorporated to meet static sensitive situations.
- 5 Contemporary Styling**
Our best seller restyled by Anthony Stoker.
- 6 Adjustable Heel Strap**
Heel strap can be worn in forward or back position.



With side vents 0865

0865W
White
Sizes 3-12



0865N
Navy
Sizes 3-12



0865LB
Light Blue
Sizes 3-12



12 Month Guarantee



0865PR
Purple
Sizes 3-9



0865BK
Black
Sizes 3-12



UK	Euro
3	36
4	37
5	38
6	39
6.5	40
7	41
8	42
9	43
10	44
10.5	45
11	46
12	47

Taking Comfort To A New Level

We took the legendary EziKlog and made it even better. The new, soft feel E-Tech material combines with an enhanced footbed to deliver even greater shift-long comfort. Whilst the GripSafe Extreme pattern offers improved slip resistance, making the **EziKlog**, the safest, most stylish clog we've ever made.



GripSafe Extreme SRC rated sole for high slip resistance to reduce slip accidents at work, it is also puncture resistant to reduce risk of sharps injury.



EziKlog[™] Toffeln



UK	Euro
3	36
4	37
5	38
6	39
6.5	40
7	41
8	42
9	43
10	44
10.5	45
11	46
12	47



12 Month Guarantee



With side vents O865

O865BA
Navy Blue & Aqua
Sizes 3-12



O865MBN
Mid Blue & Navy
Sizes 3-12



O865BR
Raspberry & Navy
Sizes 3-8



Without side vents O875

O875W
White
Sizes 3-12



O875N
Navy
Sizes 3-12



O802
Insoles
Sizes 3-12



Comforttech

Underfoot cushioning with built in anti-fatigue technology, for day long comfort



Contemporary styling

Our best seller restyled by a leading designer



Gripsafe extreme sole

Highest level of slip resistance to keep you safe



Our SmartSole Collection

The ultimate foot pain relief

Our feet are unique - the way they look, the way they feel, even the way we walk. So, when it comes to our shoes we don't take a 'one size fits all' approach. Two years in the making, it's time to discover the Toffeln SmartSole range.

Our advanced range of SmartSoles and customizable insoles will make your working day far more comfortable, no matter the shape or size of your feet.

SmartSole Insoles Available individually



0301 SOFT

0302 MEDIUM

0303 FIRM

A perfect fit with our free insoles

Because everyone's feet are different, we include 3 types of insoles free with each pair of SmartSole shoes.

UK	Euro
3	36
4	37
5	38
6	39
6.5	40
7	41
8	42
9	43
10	44
10.5	45
11	46
12	47

“I have sensitive feet and started having pain due to plantar fasciitis...very difficult to find comfortable shoes. I bought other clogs from other brands and had to send them back, but these are lovely! I bought 3 pairs, different colours, as I work in different places, and have no regrets. Best thing I ever did, don't feel pain any more as they have elevated heels. And the fact that you can choose the insole you want to use is a plus.”

Christina

Target Foot Pain Fast with SmartSole



Toffeln[®] SmartSole Breeze

A new generation of breathable, washable footwear, **SmartSole Breeze** provides slipper-like comfort with gentle support to reduce fatigue. All day cool comfort, to make your working day a Breeze.



12 Month Guarantee



370M
Mauve
Sizes 3-8



UK	Euro
3	36
4	37
5	38
6	39
6.5	40
7	41
8	42
9	43
10	44
10.5	45
11	46
12	47



370CB
Cameo Blue
Sizes 3-8



370BK
Black
Sizes 3-12



370N
Navy
Sizes 3-12



Features

- 1 Guaranteed Certified**
Toffeln Breeze are certified to EN20347:2012, OB, E, SRC. This is the reason that it makes sense to wear Toffeln rather than conventional high street brands or sports shoes as they are tested and approved safe and comfortable for use at work.
- 2 GripSafe Extreme SRC rated sole**
Slips and falls cost the average Health Trust £90,000 per year. The GripSafe Extreme sole is independently tested and approved to the high standards of SRC slip resistance.
- 3 StretchMesh Breathable Uppers**
The comfy stretch uppers mould comfortably to the feet, allowing room for natural swelling that takes place during the day and accommodating painful or swollen joints. The uppers are highly breathable, keeping feet cool and comfortable, the mesh is tested to EN ISO 20244:2011.
- 4 Washable**
SmartSole can be washed to 40°C (allow to dry in a well-ventilated area do not dry with direct heat). Remove insole prior to washing.
- 5 Lightweight**
Have you ever considered how many miles a day you walk? The average healthcare worker takes over 10,000 steps a day. Wearing lighter footwear in a healthcare setting can reduce fatigue substantially.
- 6 Raised Heel**
As a result of our research the raised heel is at the perfect pitch to distribute pressure more evenly within the shoe for improved comfort and less fatigue.



See page 62 for trainer socks

Toffeln SmartSole Clog

Target foot pain fast

Just when you thought Toffeln clogs couldn't be more comfortable, we've raised the bar again with glove-like fitting and a choice of three insoles that come free with every pair.

Features

- 1 Stretch Lining**
Smooth lining eliminates rubbing and discomfort whilst providing a comfortable fit.
- 2 GripSafe Extreme**
Slips and falls cost the average Health Trust £90,000 per year. The GripSafe Extreme sole is independently tested and approved to the high standards of SRC slip resistance.
- 3 Breathable Uppers**
The uppers are extremely breathable, washable, antibacterial, eco-friendly and highly resistant to stains.
- 4 Washable**
SmartSole can be washed to 40°C (allow to dry in a well-ventilated area do not dry with direct heat). Remove insole prior to washing.
- 5 Lightweight**
Minimises the stress and strain on feet and legs throughout the day.
- 6 Anti-Static Sole**
Specifically incorporated to meet static sensitive situations.
- 7 Raised Heel**
Helps distribute pressure within the shoe for improved comfort and less fatigue.



0340BK
Black
Sizes 3-12



0340N
Navy
Sizes 3-12



0340W
White
Sizes 3-12



12 Month Guarantee



All **SmartSole** footwear includes **THREE** pairs of insoles for your perfect fit. Also available to buy separately in UK sizes 3-12:
0301 Soft \ 0302 Medium
0303 Firm

Now available in larger sizes

UK	Euro
3	36
4	37
5	38
6	39
6.5	40
7	41
8	42
9	43
10	44
10.5	45
11	46
12	47

Toffeln[®] SmartSole Slip-On

Now you can have comfort you only dreamed of in a full slip-on shoe that has soft feel lining and has a choice of three insoles that come free with every pair.

Features

- 1 Stretch Lining**
Smooth lining eliminates rubbing and discomfort whilst providing a comfortable fit.
- 2 GripSafe Extreme**
Slips and falls cost the average Health Trust £90,000 per year. The GripSafe Extreme sole is independently tested and approved to the high standards of SRC slip resistance.
- 3 Breathable Uppers**
The uppers are extremely breathable, washable, antibacterial, eco-friendly and highly resistant to stains
- 4 Washable**
SmartSole can be washed to 40°C (allow to dry in a well-ventilated area do not dry with direct heat). Remove insole prior to washing.
- 5 Lightweight**
Minimises the stress and strain on feet and legs throughout the day.
- 6 Anti-Static Sole**
Specifically incorporated to meet static sensitive situations.
- 7 Raised Heel**
Helps distribute pressure within the shoe for improved comfort and less fatigue.

Now available
in larger sizes

12

12 Month Guarantee

COLLABORATIVE
Toffeln
COMFORT
TECH
DESIGN INNOVATION

UK
CA



0350W
White
Sizes 3-12



0350N
Navy
Sizes 3-12



0350BK
Black
Sizes 3-12

Fully lined
and easy
slip-on

UK	Euro
3	36
4	37
5	38
6	39
6.5	40
7	41
8	42
9	43
10	44
10.5	45
11	46
12	47

-
-
-
-
-
-

All SmartSole footwear includes **THREE** pairs of insoles for your perfect fit.
Also available to buy separately in UK sizes 3-12:
0301 Soft \ 0302 Medium
0303 Firm



TOFFELN SmartSole Trainer

Peak performers

At last a trainer that provides you with unparalleled comfort and complies with all the needed CE certification, has a wipe-clean upper and slip-resistant soles. Every pair comes with a choice of 3 free insoles.



GripSafe Extreme SRC rated sole, for high slip resistance to reduce slip accidents at work.



Features

- 1 Guaranteed Certified**
Toffeln Breeze are certified to EN20347:2012, OB, E, SRC. This is the reason that it makes sense to wear Toffeln rather than conventional high street brands or sports shoes as they are tested and approved safe and comfortable for use at work.
- 2 GripSafe Extreme SRC rated sole**
Slips and falls cost the average Health Trust £90,000 per year. The GripSafe Extreme sole is independently tested and approved to the high standards of SRC slip resistance.
- 3 Washable**
SmartSole can be washed to 40°C (allow to dry in a well-ventilated area do not dry with direct heat). Remove insole prior to washing.
- 4 Lightweight**
Have you ever considered how many miles a day you walk? The average healthcare worker takes over 10,000 steps a day. Wearing lighter footwear in a healthcare setting can reduce fatigue substantially.
- 5 Breathable Uppers**
The uppers are liquid repellent but highly breathable, keeping feet cool and comfortable.
- 6 Raised Heel**
As a result of our research the raised heel is at the perfect pitch to distribute pressure more evenly within the shoe for improved comfort and less fatigue.
- 7 Anti-Static Sole**
Specifically incorporated to meet static sensitive situations.



12 Month Guarantee



UK	Euro
3	36
4	37
5	38
6	39
6.5	40
7	41
8	42
9	43
10	44
10.5	45
11	46
12	47



0361BKW
Black & White
Sizes 3-12



0361BKB
Black
Sizes 3-12



0361W
White
Sizes 3-12

AktivKlog[™] Toffeln

This super-lightweight range takes our advanced clog heritage and raises the bar, with a higher closed-in heel, enhanced vents, Grip-Safe[®] sole intelligence and removable insole. This refreshingly ergonomic design was created in collaboration with biomechanics researchers from the University of Salford.



12 Month Guarantee



0230OR
Orange
Sizes 2-13



0230W
White
Sizes 2-13



0230MB
Metallic Blue
Sizes 2-13



0230N
Navy
Sizes 2-13



0230PR
Purple
Sizes 2-8



0230BK
Black
Sizes 2-13



0230G
Green
Sizes 2-13



0201
Insole
Sizes 2-13



Features

- 1 Closed-In Heel**
Designed to hold the heel securely, with enhanced arch support. The heel is raised for greater comfort.
- 2 Insole Included**
Walk with even greater comfort with the specially designed ComfortTech insole.
- 3 Adjustable Heel Strap**
Heel strap can be worn in forward or back position.
- 4 Anti-Static System**
Shock absorbing anti-static device specifically incorporated to meet any static sensitive situation.
- 5 Side Vents**
Allows your feet to breathe naturally, keeping them cool throughout the day.
- 6 Grip-Safe[®] Sole**
This unique sole not only gives excellent slip resistance, but also offers a much greater penetration resistance than other lightweight footwear available.

UK	Euro
2	35
3	36
4	37
5	38
6	39
6.5	40
7	41
8	42
9	43
10	44
10.5	45
11	46
12	47
13	48

AktivFlex[™] Toffeln

At last a trainer that provides you with unparalleled comfort and complies with all the needed CE certification, has a wipe-clean upper and slip-resistant soles. Every pair comes with a choice of 3 free insoles.

Features

- 1 Elasticated Gusset**
Slip on with ease and enjoy a secure fit.
- 2 Side Ventilation**
Allows your feet to breathe naturally, keeping them cool throughout the day.
- 3 Insole Included**
Walk with even greater comfort with the specially designed ComfortTech insole.
- 4 Stylish Rubber Laces**
Enjoy the laced up look without having any time-consuming laces to tie.
- 5 Grip-Safe[®] Sole**
This unique sole not only gives excellent slip resistance, but also offers a much greater penetration resistance than other lightweight footwear available.



“AktivFlex offer me the perfect combination of comfort and durability, and I would happily recommend them.”

Andrew Clements
Warehouse Team Member



12 Month Guarantee



Now available
in larger sizes



220MF
Metallic Fuchsia
Sizes 3-8



220MB
Metallic Blue
Sizes 3-12



220W
White
Sizes 3-12



220BK
Black
Sizes 3-12

UK	Euro
3	36
4	37
5	38
6	39
6.5	40
7	41
8	42
9	43
10	44
10.5	45
11	46
12	47

Toffeln FlexiKlog

Built in comfort

FlexiKlog is designed for comfort. It has a built-in arch and a slightly raised heel to help keep your body in the correct position, which helps you stay more comfortable all day long. There's also plenty of room for your toes to spread out, thanks to the high toe box.



UK	Euro
3	36
4	37
5	38
6	39
6.5	40
7	41
8	42
9	43
10	44
10.5	45
11	46
12	47

Features

- 1 Soft Insole**
Soft insole cushioned for ultimate comfort.
- 2 Leather Upper**
PU coated easy-to-clean leather uppers.
- 3 Anti-Static Sole**
Shock absorbing sole specifically incorporated to meet static sensitive situations.
- 4 High Toe Box**
The high toe box gives room for your toes to spread.
- 5 Ergonomically Designed**
The fully flexible sole is ergonomically designed to provide day-long comfort.



0723N
Navy
Sizes 3-12



0723W
White
Sizes 3-12



0723BK
Black
Sizes 3-12



12 Month Guarantee



0723BUP
Butterfly
Sizes 3-8



0724BFL
Blue Flower
Sizes 3-8



UltraLite Washable Clog



12 Month Guarantee



Features

- Adjustable Heel Strap**
Can be worn in forward or backward position.
- Extraordinary Comfort**
Extraordinary levels of comfort for those long shifts.
- Reduce Fatigue**
Contoured insole and raised heel distributes pressure evenly to reduce fatigue.
- Nursing Focussed Design**
Specially designed for Nursing based on our scientific research with Salford University.
- Washable**
Breathable wipe clean upper washable to 40C.
- Anti-Static Sole**
Shock absorbing sole specifically incorporated to meet static sensitive situations.
- Grip-Safe Sole**
SRC certificated slip resistance.
- Fully Certificated**
UKCA CE EN ISO 20347:2012 OB E A SRC.

O616LDP
Black Large Dog
Sizes 3-8



O618W
White (perforated)
Sizes 3-12



O618N
Navy (perforated)
Sizes 3-12



UK	Euro
3	36
4	37
5	38
6	39
6.5	40
7	41
8	42
9	43
10	44
10.5	45
11	46
12	47

O619W
White
Sizes 3-12



O619N
Navy
Sizes 3-12



O619BK
Black
Sizes 3-12



“The UltraLite Clog has remained a top choice among my customers, seamlessly blending essential healthcare prerequisites with a fun and vibrant print design.”

Harvey Leflaive
Key Account Manager



UltraLite Washable Comfort Trainer

For all day comfort and support



Features

- 1 Extraordinary Comfort**
Extraordinary levels of comfort for those long shifts.
- 2 Reduce Fatigue**
Contoured insole and raised heel distributes pressure evenly to reduce fatigue.
- 3 Nursing Focussed Design**
Specially designed for Nursing based on our scientific research with Salford University.
- 4 Washable**
Breathable wipe clean upper washable to 40C.
- 5 Anti-Static Sole**
Shock absorbing sole specifically incorporated to meet static sensitive situations.
- 6 Grip-Safe Sole**
SRC certificated slip resistance.
- 7 Fully Certificated**
UKCA CE EN ISO 20347:2012 OB E A SRC.

UK	Euro
3	36
4	37
5	38
6	39
6.5	40
7	41
8	42
9	43
10	44
10.5	45
11	46
12	47



12 Month Guarantee



0663BK
Black
Sizes 3-12



0663W
White
Sizes 3-12



Bamboo Trainer Socks
Page 62





Toffeln ClassicKlog

If you love the firmness and feel of natural wood supporting you underfoot, these clogs are ideal. With wipe-clean leather uppers, they're easy to keep clean too.

020W
White
Sizes 3-12



Features

- 1 Wooden Footbed**
Enjoy the firm wooden footbed that gives your body the correct anatomical support.
- 2 Leather Upper**
PU coated easy-to-clean leather uppers.



UK	Euro
3	36
4	37
5	38
6	39
6.5	40
7	41
8	42
9	43
10	44
10.5	45
11	46
12	47



Toffeln ErgoBoot

These white unisex boots are strong, light and have non-marking soles. They're also easy to clean, as they're washable to 50°C.

Features

- 1 Non-Marking Soles**
- 2 Anti-Static System**
Flexible PU anti-static soles.
- 3 Offers Full Protection**
Offers the best protection from fluids.



UK	Euro
3	36
4	37
5	38
6	39
6.5	40
7	41
8	42
9	43
10	44
10.5	45
11	46
12	47
13	48
14	49
15	50



12 Month
Guarantee

12 Month
Guarantee

Toffeln KlimaFlex

Made from PU material, the KlimaFlex range gives you a firm, comfortable footbed. You can also easily keep them hygienically clean, using a hot wash up to 85°C.

Features

- 1 Heel Strap**
Heel strap can be worn in forward or back positions.
- 2 Side Vents (0165 & 0168)**
Allows your feet to breathe naturally, keeping them cool throughout the day.
- 3 Lightweight PU Material**
Helps cushion your walk.
- 4 Nature-Form Footbed**
Provides long lasting comfort.

Up to
85°C
Hot Wash



With heel strap and side vents 0168



0168N
Navy
Sizes 2-15



0168W
White
Sizes 2-15

With side vents 0165



0165N
Navy
Sizes 2-15



0165MB
Mid Blue
Sizes 2-12

Without side vents 0160



0160W
White
Sizes 4-12



12 Month
Guarantee

UK	Euro
3	36
4	37
5	38
6	39
6.5	40
7	41
8	42
9	43
10	44
10.5	45
11	46
12	47
13	48
14	49
15	50

Toffeln SteriKlog

If you need the ultimate in hygiene, these long-lasting, ergonomically-shaped shoes are just the job. Ideal for visitors or short-stay personnel, they're made from strong, thermoplastic material that can be autoclaved at 134°C.

Features

- 1 Thermoplastic Material**
Long-lasting, strong thermoplastic material.



UK	Euro
3-4	35-36
4.5-5	37-38
6-6.5	39-40
7-7.5	41-42
8-9	43-44
10-11	45-46



058G
Green
Sizes 3-11



12 Month
Guarantee



Toffeln & V12 Safety Footwear

1

Did you know the average cost of slip accidents to Healthcare Trusts is £90,000. Provision of approved slip resistant footwear massively reduces this risk and associated costs.

2

An estimated one third of reported accidents result from moving and handling. The procedure guidance states: "Before undertaking any manual handling procedure, it is important to be wearing non-slip regulation footwear".

3

All V12 Safety Footwear is approved to EN 20345 and carries a slip rating.

4

All products on pages 51-58 have a safety toe cap protection to EN ISO 20345.



Intelligent Grip System

Three years in the making, V12 have developed the SRC approved Intelligent Grip System (IGS) to provide grip when it really counts in the Healthcare Industry. Following a rigorous testing process, the HSL awarded the IGS a highly prestigious 4-star.

Stability Traction System

Inspired by tractor tyre technology the SRC compliant V12 STS is a revolution in sole design. Made from highly durable PU, the STS provides superb stability and traction on loose surfaces. With firefighter standard ladder grips, a re-designed heel strike pad and non-clogging channels for exceptional water dispersal, the STS sets a new benchmark in rugged sole technology.





Cover the ground more comfortably

The EziProtekta range also comes with a **FREE** comfort insole.



Toffeln EziProtekta

Protecting the protectors

With a strong and light-as-air safety toe cap, EziProtekta keeps your pinkies protected if you drop something on your foot. They're also super-light and have a Grip-Safe® slip-resistant sole to help prevent falls.

Features

- 1 Closed-In Heel**
Specially designed to hold the heel securely.
- 2 EziProtekta Insole**
A complimentary insole is provided with every EziProtekta but its use is optional.
- 3 Adjustable Heel Strap**
Heel strap can be worn in forward or back positions.
- 4 Composite Toe Cap To EN20345:2011 SB SRC**
Tested to 200 joule impact and approved to EN20345:2011 SB SRC, this toe cap offers protection against a high level of impact and compression.
- 5 Grip-Safe® Sole**
This unique sole not only gives excellent slip resistance, but also offers a much greater penetration resistance than other lightweight footwear available.

UK	Euro
3	36
4	37
5	38
6	39
6.5	40
7	41
8	42
9	43
10	44
10.5	45
11	46
12	47



12 Month Guarantee



855W
White
Sizes 3-12



855BK
Black
Sizes 3-12



803
Insoles
Sizes 3-12



Womens Safety Footwear Range

Designed by women for women

This range of footwear is built on a last or shape specially designed for womens feet for a more comfortable fit.



1/2 of all major workplace injuries are caused by slips

Features

- 1 Intelligent Grip System (IGS) heat resistant sole
- 2 Electrostatic Dissipative (ESD)
Womens fit
- 3 Suitable for Vegans
Lightweight composite safety toecap
- 4 Flexible anti-penetration composite midsole
- 5 EN ISO 20345:2011 SIP HI HRO ESD SRC



V1925
Octane Womens IGS
Black metal free hiker
Sizes 2-8



V1945
Torque Womens IGS
Grey SIP Metal Free Hiker
Sizes 2-8



UK	Euro
2	35
3	36
4	37
5	38
6	39
6.5	40
7	41
8	42

V1955
Velocity Womens IGS
Cobalt SIP metal free safety trainer.
Sizes 2-8



V1915
Boost Womens IGS
Black S3 metal free trainer
Sizes 2-8



V1935
Turbo Womens
Grey SIP Metal free Trainer
Sizes 2-8



Mens Safety Footwear Range

Toffeln Approved safety distributor for Healthcare sector.

1/2 of all major workplace injuries are caused by slips

UK	Euro
2	35
3	36
4	37
5	38
6	39
7	41
8	42
9	43
10	44
11	46
12	47
13	48



VT201
Agile
Black SIP Safety Trainer
Sizes 3-13

- Lightweight Rubber/EVA sole
- Sport mesh and microfibre upper
- Wide fitting steel safety toecap
- Flexible anti-penetration composite midsole
- EN ISO 20345:2011 SIP SRC and protective midsole.



VT203
Speed
Blue SIP Safety Hiker
Sizes 7-11

- Lightweight Rubber/EVA sole
- Nubuck leather and woven fabric upper
- Wide fitting steel safety toecap
- Flexible anti-penetration composite midsole
- EN ISO 20345:2011 SIP SRC



VT204
React
Sand SIP Safety Hiker
Sizes 3-13

- Lightweight Rubber/EVA sole
- Nubuck leather and woven fabric upper
- Wide fitting steel safety toecap
- Flexible anti-penetration composite midsole
- EN ISO 20345:2011 SIP SRC



V1940
Torque ICS
Grey SIP Metal free Trainer
Sizes 6-13

- Intelligent Grip System (IGS) heat resistant sole
- Electrostatic Dissipative (ESD)
- Suitable for Vegans
- Lightweight composite safety toecap
- Flexible anti-penetration composite midsole
- EN ISO 20345:2011 SIP HI HRO ESD SRC



V1501.01
Caiman IGS
Graphite S3 Waterproof hiker
Sizes 3-13

- Intelligent Grip System (IGS) heat resistant sole
- Waterproof and breathable construction
- Suitable for Vegans
- Lightweight composite safety toecap
- Flexible anti-penetration composite midsole
- EN ISO 20345:2011 S3 HRO WR SRC



V6400.01
Otter STS
Black metal free, Derby Boot
Sizes 3-13

- Stability Traction System (STS) rugged sole
- Metal-free construction
- Padded moisture-wicking lining
- Lightweight composite safety toecap
- Flexible anti-penetration composite midsole
- EN ISO 20345:2011 S3 SRC



V6411.01
Oxen STS
Black metal free, S3 Derby Shoe
Sizes 3-13

- Stability Traction System (STS) rugged sole
- Metal-free construction
- Padded moisture-wicking lining
- Lightweight composite safety toecap
- Flexible anti-penetration composite midsole
- EN ISO 20345:2011 S3 SRC



VC100
Diplomat
Black S1 Derby Shoe
Size 6-12

- Rubber heat resistant sole
- Classic leather brogue style
- Soft moisture-wicking lining
- Wide fitting steel safety toecap
- Removable comfort insole
- EN ISO 20345:2011 S1 HRO SRA



V1750
Defender STS
Black S3 Waterproof Scuff Cap Boot
Size 5-13

- Stability Traction System (STS) rugged sole
- Waterproof and breathable construction
- Poron XRD side-impact ankle protection
- Lightweight composite safety toecap
- Flexible anti-penetration composite midsole
- EN ISO 20345:2011 S3 AN WR SRC



V2120

Lynx

Honey S3 Waterproof Hiker

Size 6-13

- Waterproof and breathable V-Dri lining
- This product is carbon neutral
- Intelligent Grip System (IGS) heat-resistant sole
- Lightweight composite safety toecap
- Flexible anti-penetration composite midsole
- EN ISO 20345:2011 S3 HRO WR SRC
- Metal-free construction



VT2130

Lynx

Brown S3 Waterproof IGS Hike

Size 6-13

- Waterproof and breathable V-Dri lining
- This product is carbon neutral
- Intelligent Grip System (IGS) heat-resistant sole
- Lightweight composite safety toecap
- Flexible anti-penetration composite midsole
- EN ISO 20345:2011 S3 HRO WR SRC
- Metal-free construction



V1920

Octane IGS

Black metal free hiker

Sizes 2-8

- Intelligent Grip System (IGS) heat resistant sole
- Women's fit
- Electrostatic Dissipative (ESD)
- Suitable for Vegans
- Lightweight composite safety toecap
- Flexible anti-penetration composite midsole
- EN ISO 20345:2011 S3 HI HRO ESD SRC



V1930

Turbo

Grey SIP Metal free Trainer

Sizes 6-13

- Intelligent Grip System (IGS) heat resistant sole
- Electrostatic Dissipative (ESD)
- Suitable for Vegans
- Lightweight composite safety toecap
- Flexible anti-penetration composite midsole
- EN ISO 20345:2011 SIP HI HRO ESD SRC



V1950

Velocity

Cobalt SIP metal free safety trainer

Sizes 6-13

- Intelligent Grip System (IGS) heat resistant sole
- Electrostatic Dissipative (ESD)
- Womens fit
- Suitable for Vegans
- Lightweight composite safety toecap
- Flexible anti-penetration composite midsole
- EN ISO 20345:2011 SIP HI HRO ESD SRC



V1910

Boost

Black S3 metal free trainer

Sizes 6-13

- Intelligent Grip System (IGS) heat resistant sole
- Electrostatic Dissipative (ESD)
- Womens fit
- Suitable for Vegans
- Lightweight composite safety toecap
- Flexible anti-penetration composite midsole
- EN ISO 20345:2011 S3 HI HRO ESD SRC

SafetyLite[™] Toffeln

You'd hardly know these are safety shoes; they're light, supremely comfortable and sport high tech, ultra-healthy action-breathe uppers. The wide safety toe cap means no pinching, and the slip-resistant soles could save your day.

Features

- 1 High Grade Materials**
Have been carefully selected to give you an enhanced level of comfort.
- 2 Slip-Resistant Soles**
You'll appreciate the extra protection that the slip-resistant soles give you on slippery surfaces.
- 3 Ergonomically Designed**
Foot shape provides room for your feet to spread naturally, and a wider safety toe-cap prevents toes from being pinched.

UK	Euro
3	36
4	37
5	38
6	39
6.5	40
7	41
8	42
9	43
10	44
10.5	45
11	46
12	47



12 Month Guarantee

Slip on
04165BK
Black
Sizes 3-12



04165W
White
Sizes 3-12



Lace up
04195BK
Black
Sizes 3-12



Toffeln SterileClean

Following the many requests for a simple footwear disinfecting system, we decided to research what infection control professionals advised. This research led us to produce our SterileClean™ system.


This system is simple and safe to use on all Toffeln products coming out of contaminated areas. The whole process is generally completed in less than one minute.




1 Connect cleaning pad onto handle.




2 Spray solution onto cleaning pad.



3 Wipe over footwear.



4 Discard cleaning pad.



TPH1 SterileClean

TCP50 Sponge



TSS1 Disinfectant 500ml

Toffeln EziBag



Shoes have a habit of walking when other theatre staff pick them up accidentally. So this nameable bag will make sure your Toffeln shoes stay yours.

EBO1 EziBag

Toffeln Odour Eliminator



Keep footwear smelling good for longer with this handy spray. It blitzes unpleasant odours using a fast-acting, extra-strength treatment that's antibacterial and anti-fungal.

5500 OdourEliminator

Toffeln PowerLite Insoles

For those who suffer with foot problems, these ergonomic inserts alleviate aches all day long. The Bi-Layer design has a firm but flexible support shell, built-in foot arch support, stability-enhancing heel cradle and ultra-thin VCT cushioning for complete comfort. Prevents and reduces aches associated with; plantar fasciitis, mild pronation, heel or arch pain and sore/aching feet.



PowerLite Womens 804

Code	UK	Euro
8040305	3-5	36-38
8040607	6-7	39-41
8040809	8-9	42-43

PowerLite Mens 805

Code	UK	Euro
8050607	6-7	39-41
80508	8	42
80509	9	43
80510105	10-10.5	44-45
8051112	11-12	46-47


Toffeln Insoles

Comfortable cushioning

Our range of insoles offer extra support when you need it.



EziKlog 0802 Available in sizes 3-12



EziProtekta 803 Available in sizes 3-12



AktivKlog Insole 0201 Available in sizes 2-13



AktivLite Insole 0201 Available in sizes 2-13

AktivCool[™] Bamboo Compression Socks

Soft and breathable bamboo knee-high graduated compression socks for nurses that reduces swelling and fatigue in the legs.

The **AktivCool** nurses' compression socks are made from bamboo which is known for its absorbency and amazing softness which is often compared to silk. The bamboo material is more breathable than cotton and has channels that absorb moisture, and wick it away from the skin, keeping your feet cool and dry. It is also naturally antibacterial and is completely skin friendly and free from any harmful substances according to the International Oeko-Tex 100 Standard.

Bamboo is a fast-growing renewable resource that is 100% biodegradable.

UK	Euro
2-3	34-36
4-6	37-39
6.5-8	40-42
9-10.5	43-45
11-13	46-48



Do not
tumble dry

Knee length Bamboo Socks

Material: 50% bamboo, 44% polyamide, 6% elastane (class 1, 15-21 mmHg)



ACSPOS
Pink Orange Stripe
Sizes 2-3, 4-6, 6.5-8



ACSBK
Black
Sizes 2-3, 4-6, 6.5-8, 9-10.5, 11-13



ACSBFL
Blue Flower
Sizes 4-6, 6.5-8



ACSLGS
Light Grey Stripe
Sizes 4-6, 6.5-8, 9-10.5



ACSB
Blue Stripe
Sizes 4-6, 6.5-8, 9-10.5, 11-13



ACSBUPD
Burgundy & Pink Dots
Sizes 4-6, 6.5-8



ACSBOS
Blue Orange Stripe
Sizes 4-6, 6.5-8, 9-10.5

Bamboo Trainer Socks

Material: 72 % Bamboo viscose, 26% Polyester, 2% Elastane



ATSW
White
Sizes 3-6.5, 7-10.5

Pack of 3 pairs



ATSBK
Black
Sizes 3-6.5, 7-10.5

Pack of 3 pairs

<p>GreenT EziKlog pg 19 Guava 0899GU</p>  <p>Blueberry 0899BB</p> 	<p>EziKlog V2.0 pg 23 Mid Blue & Navy 0865MBN</p>  <p>Raspberry & Navy 0865BR</p> 	<p>SmartSole Breeze pg 27 Black 370BK</p>  <p>Navy 370N</p> 	<p>SmartSole Trainer pg 33 Black & White 0361BKV</p>  <p>Black 0361BKB</p> 	<p>AktivKlog pg 35 Navy 0230N</p>  <p>Black 0230BK</p> 	<p>FlexiKlog pg 39 White 0723W</p>  <p>Black 0723BK</p> 	<p>UltraLite Washable Clog pg 42 Navy 0619N</p>  <p>Black 0619BK</p> 
<p>EziKlog V2.0 pg 21 White 0865W</p>  <p>Navy 0865N</p> 	<p>White 0875W</p>  <p>Navy 0875N</p> 	<p>SmartSole Clog pg 29 Black 0340BK</p>  <p>Navy 0340N</p> 	<p>White 0361W</p>  <p>AktivKlog pg 35 Metallic Blue 0230MB</p> 	<p>Insole 0201</p>  <p>AktivFlex pg 37 Metallic Fuchsia 220MF</p> 	<p>Butterfly 0723BUP</p>  <p>Blue Flower 0723BFL</p> 	<p>UltraLite Washable Comfort Trainer pg 43 Black 0663BK</p>  <p>White 0663W</p> 
<p>Light Blue 0865LB</p> 	<p>Insoles 0802</p> 	<p>White 0340W</p> 	<p>Purple 0230PR</p> 	<p>Metallic Blue 220MB</p> 	<p>UltraLite Washable Clog pg 41 Large Dog 0616LDP</p> 	<p>ClassicKlog pg 45 White 020W</p> 
<p>Purple 0865PR</p> 	<p>SmartSole Insoles pg 25 0301 Soft 0302 Medium 0303 Firm</p> 	<p>SmartSole Slip-On pg 31 White 0350W</p> 	<p>Green 0230G</p> 	<p>White 220W</p> 	<p>White (perforated) 0618W</p> 	<p>ErgoBoot pg 46 Half 0900W</p> 
<p>Black 0865BK</p> 	<p>SmartSole Breeze pg 27 Maive 370M</p> 	<p>Navy 0350N</p> 	<p>Orange 0230OR</p> 	<p>Black 220BK</p> 	<p>Navy (perforated) 0618N</p> 	<p>Full 0153W</p> 
<p>Navy Blue & Aqua 0865BA</p> 	<p>Cameo Blue 370CB</p> 	<p>Black 0350BK</p> 	<p>White 0230W</p> 	<p>FlexiKlog pg 39 Navy 0723N</p> 	<p>White 0619W</p> 	

<p>KlimaFlex pg 47 Navy 0168N</p>  <p>White 0168W</p> 	<p>V12 Safety Womens pg 53 Velocity Womens IGS V1955</p>  <p>Boost Womens IGS V1915</p> 	<p>V12 Safety Mens pg 56 Torque ICS V1940</p>  <p>Caiman ICS V1501.01</p> 	<p>V12 Safety Mens pg 57 Turbo V1930</p>  <p>Octane V1920</p> 	<p>Accessories pg 59 SterileClean TPHI</p>  <p>Sponge TCP50</p> 	<p>AktivCool Bamboo Trainer Socks pg 62 White ATSW</p>  <p>Black ATSBK</p> 	<p>pg 62 Blue Orange Stripe ACSBOS</p> 
<p>White 0165W</p>  <p>Mid Blue 0165MB</p> 	<p>Octane Womens IGS V1925</p>  <p>Torque Womens IGS V1945</p> 	<p>Otter STS V6400.01</p>  <p>Oxen STS V6411.01</p> 	<p>Velocity V1950</p>  <p>Boost V1910</p> 	<p>Disinfectant 500ml TSSI</p> 	<p>AktivCool Bamboo Compression Socks pg 61 Pink Orange Stripe ACSPOS</p> 	<p>Light Grey Stripe ACSLGS</p> 
<p>Navy 0165N</p> 	<p>Turbo Womens IGS V1935</p> 	<p>Diplomat VC100</p> 	<p>SafetyLite pg 58 Black 04165BK</p> 	<p>EziBag EB01</p> 	<p>Blue Stripe ACSBS</p> 	<p>Blue Flower ACSBFL</p> 
<p>SteriKlog pg 48 Green 058G</p> 	<p>V12 Safety Mens pg 55 Agile VT201</p> 	<p>Defender STS V1750</p> 	<p>White 04165W</p> 		<p>Black ACSBK</p> 	<p>Burgundy & Pink Dots ACSBUPD</p> 
<p>EziProtekta pg 51 White 855W</p> 	<p>VT203 220MB</p> 	<p>pg 57 Lynx VT2130</p> 	<p>Black 04195BK</p> 	<p>Odour Eliminator 5500</p> 	<p>Black ACSBK</p> 	
<p>Black 855BK</p> 	<p>VT204 220W</p> 	<p>Lynx V2120</p> 	<p>Accessories pg 60 Insoles 803</p> 			

SURGICAL PRODUCTS

Tried and trusted, Toffeln has been supplying hospitals for over 40 years.

Working alongside surgeons, clinicians, theatre managers and procurement professionals, we provide innovative and cost-effective solutions that enhance surgical capability and improve patient outcomes.

Delivered from our UK stock and guaranteed certified, we supply a well-loved range of theatre equipment and consumables.

“We are excited to share that our surgical product line has experienced substantial growth this year. With numerous new additions available through the NHS

Supply Chain, you can be confident that each product surpasses the healthcare industry’s quality standards”

Harry Turner Sales Director



Surgical Products



67-69

UltraView Headlight


Meet the world's lightest and brightest Cordless Headlight



70-72

Postpartum Balloon

Discover the NEW Postpartum Balloon System



101-105

ErgoTrics

Patient moving with compressed air.



73-74

Instrument Cassettes



75-81

Disposable Gowns & Clothing



82-88

Sterile Drapes



89-92

Face & Eye Protection



93-95

Theatre Safety



96-98

Laparoscopic Products



99-104

Patient Moving & Positioning

UltraView Headlights

The world's brightest headlight. The Toffeln surgical headlight is a match for any challenge – and any challenger.

It has an output of up to 275,000 lux and a beam size of 100mm up to 200mm. Even better; The HotSwap technology means its rechargeable batteries can be swapped while in use, so there's no downtime.

Because it has no leads, this surgical headlight helps to eliminate surgeons neck pain caused by the weight of the leads in a conventional headlight and eliminates the safety issues caused by trailing cables.

“At Toffeln, our team is dedicated to offering expert advice on tailoring the ideal headlight setup for your unique surgical specialisation. We also take pride in offering an exceptional aftercare service that sets us apart from the rest.”

Taylor Whiteside
Headlight Technician



UltraView Extreme

The UltraView Extreme cable-free super-bright headlight is a real help to surgeons. Weighing only 312g and with a comfortable padded headband, its beam can be adjusted easily from wide-reaching, down to a very precise 90 millimetres.

Features

- 1 Li-Po Cell Battery**
Two rechargeable batteries are embedded within the headband. Don't worry about which battery is charged - the batteries talk to each other. With our HotSwap technology, a battery can be changed from either side without losing light. Batteries last up to 8 hours.
- 2 Adjustable Headstrap**
Fully adjustable by twisting
- 3 Adjustable Headband**
Fully adjustable by twisting
- 4 Padded Headband**
Ergonomically designed for comfort throughout extended periods of use.
- 5 Adjustable Light Positioning**
Easy to direct light beam.
- 6 Focusable LED Light**
Zoom focus high power LED light.

Bezel Options

We offer three different types of bezel options giving you the flexibility to select the headlight that suites you.



SHL100

Brightness: 100,000 lux
Spot Range: 10-100mm
(At 25cm working distance)
Specialities:
• Dental
• Oral & Maxillofacial
• ENT
• Veterinary



SHL200

Brightness: 200,000 lux
Spot Range: 70-90mm
(At 25cm working distance)
Specialities:
• Plastic, Cosmetic & Recon
• Oral & Maxillofacial
• Vascular
• Obs & Gynae
• ENT
• Veterinary



SHL250

Brightness: 275,000 lux
Spot Range: 100-200mm
(At 25cm working distance)
Specialities:
• Cardio & Cardiothoracic
• Plastic, Cosmetic & Recon
• Spinal
• Trauma
• Orthopaedic
• Neuro
• General & Colorectal
• Urology



UltraView Accessories



BC200

Four Bay Charger

Charging time: 2.5 hours

For convenient fast charging



RB200

Rechargeable Batteries

(2 included)

Type: LiPo

Capacity: 2260mAh

Weight: 48g



Protective Case (included)



Magnification Lenses



SHLGPI

GoPro Mount



SHLFV005

High Impact Clear Faceshield

Mounts on Swivel

Introducing The ObCare Postpartum Balloon

New

We are proud to introduce the ObCare Postpartum Balloon, the latest addition to our surgical range and exclusively available through Toffeln in the UK. As a trusted supplier to the medical sector, we understand the importance of providing healthcare professionals with safe, effective, and reliable medical solutions.

The ObCare Postpartum Balloon System is an innovative solution for the management of postpartum haemorrhage. Backed by medical research and meeting the BS EN ISO 20696:2018 requirements for single-use sterile urethral catheters, the ObCare Postpartum Balloon systems ensures you a safe and reliable solution.

If you have any questions or would like to request a demonstration or free trial, please don't hesitate to contact our dedicated customer care team or your Key Account Manager, who will be happy to assist.




ObCare[™] Postpartum Balloon

The disposable sterile system comes complete with two syringes, a bag spike with check valve, and a dual lumen catheter. For the treatment of postpartum haemorrhage, the deflated balloon can be positioned in the uterine cavity either trans-vaginally or via the caesarean section incision. Once the correct position has been confirmed, the balloon can be promptly inflated with saline solution to provide the necessary pressure to the uterine wall. The dual lumen catheter allows for easy monitoring of the haemorrhage during the procedure until the bleed has stopped.

Features


- 1 The systems come complete with:**
 - Bag spike with check valve
 - Dual lumen catheter
 - Two syringes
- 2 Made from 100% medical grade silicone**
- 3 Durable design**
Can hold up to 4ltrs of saline solution
- 4 No protrusion at the distal end**
- 5 Safe & Reliable Supply**

PPB0206
ObCare Postpartum Balloon
Box of 1



Pear-Shaped Form

The ObCare Postpartum Balloon has a unique pear-shaped form designed to fit the uterine anatomy perfectly and providing the patient with a secure and ergonomic fit.



Reduced Risk Of Harm

Unlike other systems, the blood drainage channel is positioned inwards, allowing the passage of blood to flow more effectively and, with no protrusion at the distal end, the risk of further harm to the uterus is significantly reduced.



Safe & Reliable

The ObCare Postpartum Balloon System by Toffeln UK is available to dispatch within 24 hours of ordering, enabling you to confidently manage the risk of postpartum haemorrhage safely, efficiently, and reliably.



Scan me to find out more about the ObCare Postpartum Balloon



“We’re really excited about this new addition to the Toffeln Surgical range. We are always on hand to answer any questions, demonstrate the ObCare or arrange a free trial.”

Braden Leflaive
Key Account Manager



Instrument Cassettes



Introducing the latest addition to our surgical range: SteriLockBox, instrument protection cassettes by Toffeln. Designed specifically for healthcare, this exciting new range offers a cost-effective solution that minimises the risk of damage to delicate surgical instruments, ensuring a safe and secure environment for storage, transportation and decontamination by Sterile Services.

Our instrument cassettes come in a range of sizes, specifically designed to cater to the unique needs of your department. You can browse our complete range of sizes and options in this section or on our website. As always, if you have any questions or require further information, please don't hesitate to contact our dedicated customer care team or your assigned key account manager.

“Ask about our bespoke customisation service.”

Ryan Storey
Key Account Manager



Benefits

- Robust protection during transportation of delicate instruments
- No need to remove instruments during decontamination process
- Teflon coated silicone reduces adhesion to allow for optimal cleaning
- Easily customised to individual instruments or sets
- Reduction in repair and replacement costs
- Inexpensive and long lasting
- Barcode panel to allow for quick identification.
- Wide range of sizes

Features

- 1** 305 surgical grade stainless steel
- 2** Electro-Plated and Polished
- 3** No sharp edges
- 4** Teflon coated silicone to aid washing
- 5** V-Design strip corrects tension



SteriLockBox

Our instrument cassettes come in a range of sizes, specifically designed to cater to the unique needs of your department.



IPC5940 Micro
Size 59 x 40 x 15mm
Ideal for little clips & adaptors.



IPC7459 Mini
Size 74 x 59 x 24mm
Ideal for micro instruments, larger clips & adaptors.




IPC18045 Uno
Size 180 x 45 x 25mm
Available with a standard cut to hold a single instrument, or uncut for optimal customisation.



IPC18050 Trio
Size 180 x 50 x 24mm
IPC20050 Trio Plus
Size 200 x 50 x 24mm
Available with standard cuts to hold 3 instruments, or uncut for optimal customization.



IPC18060 Quintet
Size 180 x 60 x 22mm
IPC20080 Quintet Plus
Size 200 x 80 x 22 mm
Available with standard cuts to hold 5 instruments, or uncut for optimal customization.



IPC185134 Dect
Size 180 x 134 x 25mm
Available with standard cuts to hold 10 instruments, or uncut for optimal customization.

Disposable Gowns & Clothing

Why Toffeln Gowns?

Put simply, you can trust Toffeln because we carefully control not only the manufacture of the gowns but also the production of our technical fabrics. The seams of our garments are ultrasonically welded for the highest levels of seam integrity. All processes follow the requirements of ISO 13485:2016 and our sterile products are packed by our manufacturing partner in a state of the art facility, class 100,000 cleanroom environment. Each Gown is latex free and exceeds the standards of EN13795. We carry stock of our gowns in our Bristol warehouse for fast despatch and have a large production capacity to fulfil NHS and Private Healthcare contracts.

“Thanks to our extensive stock-holding facility located in Bristol and our strong collaborative relationship with our manufacturing facility, we can promise a reliable source of high quality surgical gowns to the NHS and UK Healthcare.”

Oliver Turner Key Account Manager

Our gowns are generously cut to allow freedom of movement

Material Guide



EcoTex

Manufactured from a mix of 100% Recycled SS spunbonded and SMS spunbonded melt-blown fabric. It is soft and comfortable to wear and is embossed with a texture to provide no-see through modesty properties. It is a favourite for patient gowns, scrubs and warming jackets.



SoftTex

Non-woven Spunlace fabric using pulp and polyester fibres. Spunlace is popular for its softness and cotton-like feel. It has fluid repellent qualities.



LiquidTex

This is a four-layer SMMS fabric made of an inner and outer layer of spunbonded non-woven fabric, sandwiched between this are two layers of melt-blown non-woven. This fabric provides an excellent bacterial and fluid repellent barrier, whilst being cool and breathable to wear.



BarrierTex

SMMS with PE coating this is the Ultimate protection. The coating is impervious but allows a degree of breathability where the highest levels of fluid protection are required.



EN13795-1: 2019

This is the harmonised performance standard for single-use and reusable surgical gowns, intended to minimise the transmission of infective agents between patients and clinical staff during medical procedures. The standard specifies mutually agreed test methods including resistance of the fabric to both wet and dry microbial penetration, wet and dry bursting and tensile strength, liquid penetration, linting and microbial cleanliness. All Toffeln Surgical Gowns meet or exceed EN13795.



OEKO-TEX

Our gowns are made from textiles that have been tested for harmful substances, from yarns to the finished product.



NHS Supply Chain

Our strong working relationship with the NHS means that many of our surgical gowns and clothing are available direct from NHS Supply Chain. Look out for this symbol and NPC code.

Latex Free	Fluid protection	LiquidTex
Sterile	OEKO-TEX	BarrierTex
CE Marked	EcoTex	EN13795
Single use	SoftTex	Recycled

Toffeln SterileGowns

Gown Usage

Department	TG002/ TG010	TG011	TG011R/ TG003	TG012
ENT		✓	✓	
General Surgery		✓	✓	✓
Laparoscopic		✓	✓	
Patient Examination	✓			
Ophthalmic	✓	✓		
Major General			✓	✓
Cardiovascular			✓	✓
Orthopaedic		✓	✓	✓
Urology			✓	✓
Maxillofacial		✓	✓	
Neurosurgical			✓	✓
Gynaecology				✓

TG003R Comfort Spunlace, Reinforced, High Performance Gown

Sizes: L, XL, XXL
Cases of 50 per size

Features:

- Soft Comfort-Spunlace for outstanding comfort and breathability
- Double sewn seams for strength
- Knitted cuffs
- Supplied with two hand towels in a sterile wrap



Scan me
to find out
more about
our Gowns



NHS
Supply Chain
Product



TG011
Standard Surgical Gown
Sizes: S, M, L, XL, XXL, Extra-Long: XL-L, XXL-L

NHS NPC: BWK85081 (S) / BWK85089 (M)
BWK85085 (L) / BWK85082 (XL)*

Cases of 50 per size

Features:

- Manufactured from soft SMSS comfort fabric, with fluid protection
- Supplied with two hand towels in a sterile wrap
- Ultrasonic welded seams
- Knitted cuffs

NHS
Supply Chain
Product



TG011R
Reinforced Surgical Gown, High Performance
Sizes: S, M, L, XL, XXL, Extra-Long: XL-L, XXL-L

NHS NPC: BWK85080 (S) / BWK85078 (M)
BWK85077 (L) / BWK85083 (XL) / BWK85075 (XL-Long)
BWK85079 (XXL) / BWK85076 (XXL-Long)

Cases of 50 per size

Features:

- Reinforced in frontal areas and forearms for higher fluid protection
- Manufactured from soft SMSS comfort fabric, with fluid protection
- Supplied with two hand towels in a sterile wrap
- Ultrasonic welded seams
- Knitted cuffs

NHS
Supply Chain
Product



TG012
Impervious Surgical Gown, High Performance
Sizes: L, XL, Extra-Long: XL-L, XXL-L

NHS NPC: BWK85084 (L) / BWK85087 (XL)
BWK85088 (XL-L) / BWK85086 (XXL-L)*

Cases of 50 per size

Features:

- Impervious laminated material for high performance fluid protection
- Supplied with two hand towels in a sterile wrap
- Ultrasonic welded seams
- Knitted cuffs

Non-Sterile Gowns Toffeln



- LATEX FREE
- OT
- LT
- CE
- ♻️

TG010 Non-Sterile Isolation Gown

Sizes: L, XL
Cases of 80 per size

- Features:**
- Manufactured from SMSS comfort fabric, with fluid protection
 - Ultrasonic welded seams
 - Knitted cuffs



- ET
- LATEX FREE
- ♻️
- ♻️

TG020 Non-Sterile Patient Gown

Sizes: L, XL
NHS NPC: BWK85029 (L) / BWK85030 (XL)
Cases 100 per size

- Features:**
- Darker blue fabric for modesty
 - Velcro neck fastener and waist ties
 - Manufactured from a mix of 75% recycled SS spunbonded and SMS spunbonded melt-blown fabric.



- ♻️
- ♻️

TG002 Sleeved Fluid Protection Gown

Sizes: One Size
Cases of 100

- Features:**
- Economical Fluid repellent plastic gown
 - Hydro-static pressure head tested
 - Microbial cleanliness tested
 - Thumb loop cuff
 - Lightweight & flexible
 - 16 Micron Polyethylene

Warming Jackets Toffeln 75% Recycled



- ET
- LATEX FREE
- OT
- CE
- ♻️
- ♻️

TG015 Disposable Short Sleeve Warming Jacket

Sizes: M, L, XL, XXL
NHS NPC: BWK2206 (M) / BWK2209 (L) / BWK2203 (XL) / BWK2202 (XXL)
Cases of 50 per size

- Features:**
- Extra layer of comfort
 - Compatible with infection control measures
 - Manufactured from recycled low lint, comfort fabric
 - Ultrasonic welded seams
 - Knitted cuffs
 - Popper stud front



- ET
- LATEX FREE
- OT
- CE
- ♻️
- ♻️

TG014 Disposable Long Sleeve Warming Jacket

Sizes: M, L, XL, XXL
NHS NPC: BWK2208 (M) / BWK2207 (L) / BWK2205 (XL) / BWK2204 (XXL)
Cases of 50 per size

- Features:**
- Extra layer of comfort
 - Compatible with infection control measures
 - Manufactured from recycled low lint, comfort fabric
 - Ultrasonic welded seams
 - Knitted cuffs
 - Popper stud front

Call us for free samples of Gowns & Clothing

Disposable Scrubs 75% Recycled



ET

LATEX FREE

OT

CE

2

♻️

TG016 Disposable Scrub Suit Set

Sizes: S, M, L, XL, XXL

NHS NPC: BWO169 (S) / BWO176 (M)
BWO177 (L) / BWO162 (XL)

Cases of 75 per size

Features:

- Includes jackets and trousers
- Darker Blue Fabric for modesty
- Compatible with infection control measures
- Manufactured from recycled low lint, comfort fabric
- Stitched seams
- Short sleeves

TG017 Disposable Scrub Suit Tunic

Sizes: S, M, L, XL, XXL

NHS NPC: BWO175 (S) / BWO170 (M)
BWO166 (L) / BWO163 (XL) / BWO172 (XXL)

Cases of 125 per size

TG018 Disposable Scrub Suit Trousers

Sizes: S, M, L, XL, XXL

NHS NPC: BWO179 (S) / BWO165 (M)
BWO164 (L) / BWO171 (XL) / BWO168 (XXL)

Cases of 125 per size



Sterile Drapes

Toffeln now offer a wide range of equipment and patient drapes for your convenience. From C-Arm, Mini C Arm, and Microscope covers, to light handle covers, u-split, trolley and table drapes, to full operation drape sets.

Let Toffeln make it simple, call us to arrange for a free audit of your sterile drape requirements and we can offer a cost-effective single source supply that's based on your specific usage and requirements.

“Our close partnership with our new FlexiDrape manufacturer means we can offer a cost effective and reliable supply to the NHS and UK Healthcare. Call us today for a Free Drapes Audit.”

Harvey Leflaive
Key Account Manager



Equipment Drapes

Toffeln

Toffeln provides an ever increasing range of Equipment Drapes. The covers are latex free and manufactured to exacting standards of quality, offering a competitive and easy to use option.

Contact us for a free audit of all your Equipment Cover needs, whatever the make or size we can offer you a one stop shop for all your equipment covers, you can simply call off at competitive rates.

We can offer covers for many brands including **Siemens, Philips, Ziehm, Toshiba, GE, OEC, C-Bow, Shimadzu, Hologic.**



C Arm Covers

Product Number	NHS NPC	Description	Upper Cover Intensifier (W x L)	Lower Cover X-Ray Tank (W x L)	Units per box	Units per case
TC55130C	VEB1199	C-Arm Cover Siemens Siremobil - Set of 3 parts	75x90cm	75x90cm	8	32
TC55150C	VEB1200	C-Arm Cover Siemens Arcadis - Set of 3 parts	75x90cm	75x90cm	8	32
TC55220C	VEB1210	C-Arm Cover Philips BV29				

Mini C Arm Covers

Product Number	NHS NPC	Description	Bag Size (W x L)	Units per box	Units per case
TC55466C	VEB1209	Mini C-Arm Cover with Band	138x160cm	10	40
TC55476C	VEB1206	Mini C-Arm Cover with Band	138x220cm	8	32



QwikCover

Toffeln



TC45310C QwikCover Video Camera Drape – Closed System

Box of 100

NHS NPC: VEB1202

Features:

- Covers all types and brands of video cameras
- Easy to use
- Secures tightly to camera head
- Latex free
- Adhesive free
- Tear resistant polyethylene ring



TC45301C QwikCover Video Camera Ring Drape

Box of 100

NHS NPC: VEB1204

Features:

- Covers all types and brands of video cameras
- Easy to use
- Secures tightly to camera head
- Latex free
- Tear resistant polyethylene ring



TC35030C QwikCover Ultrasound Probe Cover

Toffeln Probe Cover ensures the protection of ultrasound transducers and avoids cross contamination.

18 x 30cm

Box of 200

NHS NPC: VEB1207

Features:

- Offers additional sterile protection during procedures
- Available in different lengths for complete coverage from transducer array to connector
- Easy application thanks to telescopic folding
- Protecting expensive imaging equipment during operation



TC45350C QwikCover Folded Camera Cover

Box of 100

NHS NPC: VEB1205

Features:

- Centralized hanger hall
- Telescopically folded
- Single layer thickness = 40 micron
- Eight-layer folding allows shorter cover length
- Easier and faster to pass camera
- Compatible with all camera models



LightHandleCovers



NHS
Supply Chain
Product

TC65404C ProCover

(Single Pack) Box 480

NHS NPC: VEB1203

These high-quality light handle covers are guaranteed not to split or tear. The ProCover offers an ideal solution that eliminates the costs of stocking reusable handles.



NHS
Supply Chain
Product

TC65460 UniCover

Box 60

NHS NPC: VEB1208

TC65460C UniCover

Carton quantity 240

NHS NPC: VEB1201

This patented light handle cover has been cleverly designed with twin sets of teeth to give maximum grip to the handle. It fits well on all manufacturers original light handles and all universal adapters, stopping the need for expensive OEM brands of light adaptors.

CamCover



TC65470 Square CamCover

Box of 120

Foam ring light handle cover with clear lens. This cover fits the KLS Martin SurgiCam video camera system. Covers for other brands of camera available on request





TC65471 Round CamCover

Box of 120

Foam ring light handle cover with clear lens. For use with video handle covers.

FlexiDrape



Trolley & Table Drapes

TSD0020019 Small

Box of 125

NHS NPC: VEB85072

Dimensions:
W100 x L150cm

Reinforced Area:
W50 x L150cm

TSD0020001 Medium

Box of 70

NHS NPC: VEB85071

Dimensions:
W150 x L200cm

Reinforced Area:
W100 x L200cm

TSD0020010 Large

Box of 70

NHS NPC: VEB85081

Dimensions:
W150 x L240cm

Reinforced Area:
W75 x L240cm

NHS
Supply Chain
Product

U-Split Drapes

TSD0470310 Small

Box of 80

NHS NPC: VEB85078

Dimensions:
W120 x L130cm

Reinforced Area:
W75 x L75cm

Split: W7 x L70cm

Adhesive Strip: W5cm

TSD0470311 Medium

Box of 40

NHS NPC: VEB85074

Dimensions:
W200 x L260cm

Reinforced Area:
W70 x L100cm

Split: W7 x L102cm

Adhesive Strip: W5cm

TSD0470312 Large

Box of 40

NHS NPC: VEB85079

Dimensions:
W228 x L260cm

Reinforced Area:
W75 x L100cm

Split: W20 x L100cm


Adhesive Strip: W5cm

NHS
Supply Chain
Product



“Our close working relationship with our manufacturing facility means that we are able to create bespoke FlexiDrapes to your specific requirements. Ask your key Account Manager for more information.”

Danielle Leflaive
Sales Support



General Adhesive Drapes



TSD0420965
Small

Box of 200

NHS NPC: VEB85075


Dimensions:
W75 x L75cm

TSD0420447
Medium

Box of 150

NHS NPC: VEB85073

Dimensions:
W100 x L85cm



OP Sheet




TSD0420301

Box of 80

NHS NPC: VEB85077

Dimensions:
W150 x L240cm

Adhesive Strip: L100cm



U-Split Suction & Diathermy




TSU4004808


Box of 250

NHS NPC: VEB85076


Dimensions:
W35 x L42cm




4x hand towels



4x Adhesive Drapes




1x Instrument Table Drape




1x Op-tape


Length (L)



1x Back Table Drape



2x Adhesive Side Drapes



2x Adhesive Side Drapes



Adhesive Strip



Reinforced Area

Universal Drape Set 1



TSP0200492

Box of 13

NHS NPC: VEB85080

Set comes complete with:

- 1x Op-tape 10 x 55cm
- 4x hand towel 40 x 40cm
- 2x Adhesive drape 100 x 85cm
- 1x Adhesive drape 183 x 183cm
- 1x Adhesive drape 250 x 150cm
- 1x Instrument table drape 200 x 150cm

Universal Drape Set 2



TSP0200491

Box of 15

NHS NPC: VEB85070

Set comes complete with:

- 2x Adhesive drape 100 x 75cm
- 2x Adhesive drape 170 x 170cm
- 1x Back Table Drape 250 x 150cm

Face & Eye Protection

Face and Eye Protection is highly personal, which is why Toffeln can be trusted to supply comfortable, well-made products. Our face visors and FlexiSpecs offer premium comfort with crystal clear vision. Our masks are designed for best-in-class protection, that can be worn with ease. The EyeSafe Tabs are a great new addition to the range helping to improve patient outcomes. So play safe every day with these tried and tested products.

New



FMH50E
Hypoallergenic IIR Medical Face Mask with Earloops

Boxes of 50

Features:

- Hypoallergenic
- Type IIR EN14683:2019
- Bacterial Filtration > 98%
- Triple layer spunbonded filtration

Face-Shields & Masks ^{Toffeln}



SMV050
Type II R Tie Back Face Mask with Anti-Fog Visor

Boxes of 50

Features:

- UKCA, EN14683:2019
- Optically clear visor with anti-fog coating
- Anti-fog adjustable nose strip to prevent fogging
- Triple layer spunbonded filtration



FM001
Disposable Face Masks, 3 ply, Type II R Fluid Resistant, Ear-Loop

Boxes of 50

Features:

- Type IIR EN14683:2019
- Bacterial Filtration > 98%
- Triple layer spunbonded filtration



FMIITB
Disposable Face Masks, 3 ply, Type II R Fluid Resistant, Tie Back

Boxes of 50

Features:

- Type IIR EN14683:2019
- Bacterial Filtration > 98%
- Triple layer spunbonded filtration



FV001
Premium Comfort Face Shield

Boxes of 100

Features:

- ISO: 13485:2016
- Ergonomic foam headband
- Can be used with eyeglasses and masks
- Elastic headband for stability
- Double sided protective film to protect from scratches



Toffeln FlexiSpecs



FSF0050 FlexiSpecs Premium Eye Shield Frame



Boxes of 50
Assorted Colours



Features:

- ISO: 13485:2016
- Antibacterial and anti-fog
- Ergonomic patented S-Flex Frame design for increased comfort and flexibility.
- Double-side protective film with easy peel tabs
- Can be fitted over prescription eyewear and surgical masks



FSL0100

FlexiSpecs Premium Eye Shield Lens

Boxes of 100

Features:

- Easy to replace, super clear lens
- Anti-fog

Toffeln EyeSafe



ESA0050 EyeSafe Adult Eye Safety Tabs



ESC0050 EyeSafe Child Eye Safety Tabs

Boxes of 50



Features:

- Protects the eyes during anaesthesia
- Easy peel does not harm skin, eyebrows, or eyelashes, when removing post operation
- Eyes do not dry out even in long operations
- Medical grade material, suitable for sensitive skin
- Reduces risk of corneal abrasion



Advantage for Patients

- Does not hurt the skin, eyebrows and eyelashes while removing (post operation)
- Protects eyes from external splashes during the operation
- Eyes do not dry out even during long operations
- Medical grade material avoids any skin irritation (even on sensitive skin)
- Reduces corneal abrasion risk

Advantages for Healthcare Professionals

- The transparent structure makes for easy pupil dilation control and is visible at any time without removing the EyeSafe tab
- The coloured, non-sticky handle allows for easy and fast application and removal. This reduces the operational risk and saves on application time

Theatre Safety

Specialists in theatre safety, Toffeln tackles the most common risk in hospitals. Everyday thousands of medical staff rely on TripStop to reduce the risk of falling over leads and cables, use our SharpSecure system to avoid injury from sharps, or our SpillStop absorbent sheets to effortlessly soak up fluids.



SharpSecure

Needle and blade disposal system

Toffeln SharpSecure needle counters provide the safest and easiest way of disposing of needles and blades during procedures and includes features such as a blade remover and magnetic and foam areas. Across the UK, hospitals and operating theatres are fast becoming aware that they need a safer alternative than basic sharps pads. SharpSecure provides a safe and effective solution, whilst providing a decreased risk of injury for staff.

TC65659C
SharpSecure

Boxes of 100



TripStop



Trip prevention pipe and cable cover

Slips and trips are the biggest cause of serious injuries to health care workers. This brightly coloured, single use, cable and pipe cover helps to prevent accidents.



TS40R
TripStop Cable Cover
40m Roll

NHS NPC: VWL308

Features

- Fixes with low tack adhesive on outer edges only
- No sticky residue left behind
- Quality guaranteed Toffeln product
- Quick and easy to use
- Fitted guillotine for easy dispensing



TSDISP
TripStop Dispenser

Features

- Pre-drilled holes to attach to a wall if necessary
- Fitted guillotine for easy dispensing



New




SpillStop[™]

Super Absorbent Medical Matting

Ideal for use around scrub sinks or for operations that involve higher levels of body fluids. Helps to reduce risks of slips related to contaminated floors and cuts down cross contamination.

Features:

- Multi-layer matting with non-slip backing
- Triple-layer structure for high absorbency up to 3.5 litres per metre sq
- Micro-perforated surface, allows ultra-fast absorption
- Sealed back and edges prevents leakage of fluids through the edge or back of the mat protecting the floor
- Available as cut sheets or rolls




The upper layer is a **micro-perforated sheet** which allows instant liquid penetration and fast dry; holes are reverse-cone-shaped to avoid backflow of fluids to the surface. The central layer is added with **super absorbent polymer (SAP)** capable of retaining large volumes of liquids.

Product Number	NHS NPC	Description	Size (W x L)	Box Qty
SSROLL500090	TBR85001	Spill Stop Medical Absorbent Roll	50m x 90cm	1 roll
SSMAT15090	TBR85003	Spill Stop Medical Absorbent Mats	150cm x 90cm	35 pieces
SSMAT10090	TBR85000	Spill Stop Medical Absorbent Mats	100cm x 90cm	50 pieces
SSMAT5090	TBR85002	Spill Stop Medical Absorbent Mats	50cm x 90cm	100 pieces

Laparoscopic Equipment

When there is no margin for error you can rely on Toffeln. Our Laparoscopic equipment is trusted by surgeons for its precision and reliability. EndoZenith is ergonomically designed for maximum comfort in use, helping to reduce muscle strain and fatigue. Our ViewMax scope warmer activates in seconds to warm laparoscopes and endoscopes to body temperature, preventing lens condensation. Thanks to its convenient popper fastening no tapes are required, speeding up procedures whilst reducing costs. Trust Toffeln every time.





LA75YHGR-R
EndoZenith Yohanne Grasper with Ratchet

Box of 6

Features:

- Smooth action shaft wheel-fail
- Very easy to use – fingers stay positioned
- Easy action, positive engaging mechanism
- Very ergonomically designed for maximum comfort during use
- Helps reduce muscle strain and fatigue



LA75YHGR
EndoZenith Yohanne Grasper

Box of 6

Features:

- Smooth action shaft wheel-fail
- Safe on/off ratchet mechanism
- Very easy to use – fingers stay positioned
- Easy action, positive engaging mechanism
- Very ergonomically designed for maximum comfort during use
- Helps reduce muscle strain and fatigue



LA75CVSC
EndoZenith Curved Scissors

Box of 6

Features:

- Smooth action shaft wheel-fail
- Safe on/off ratchet mechanism
- Very easy to use – fingers stay positioned
- Easy action, positive engaging mechanism
- Very ergonomically designed for maximum comfort during use
- Helps reduce muscle strain and fatigue



 **TC65LSWM**
ViewMax Sterile Liquid Scope Warmer

Box of 30



Features:

- Activates in seconds to warm laparoscopes and endoscopes to body temperature to prevent fogging
- Simple to use popper fastening



 **TC65701**
Anti Fog Solution

Box of 200



NHS NPC: FBE24632

Features:

- Helps prevent and remove condensation on the distal lens of the laparoscope, for improved visibility in procedures
- Sterile pouch with 6ml easy squeeze bottle and adhesive backed sponge pad
- Non-abrasive sponge, radio opaque



 **TC65SCWM**
Electric Scope Warmer

Box of 25



TC65SCWB
Replacement Rechargeable Battery

Qty 1



Features:

- Prevents and removes condensation from the distal lens of the laparoscope, for improved visibility in procedures
- 1x rechargeable battery unit and transformer supplied with your first order
- Ideal for longer procedures
- Adjustable temperature control

New

Patient Moving & Positioning

We are delighted to unveil our brand-new ErgoTrics range! This new addition to Toffeln's Patient Positioning & Moving Range empowers healthcare professionals to effortlessly position and transfer patients with the power of compressed air. Exclusively available through Toffeln in the UK, this new range ensures the safety and well-being of both patients and healthcare professionals.

In addition to this new range, we continue to offer our trusted patient gel positioners, which provide superior pressure relief to patients, whilst being easy to clean and x-ray translucent.

For more information on any of our Patient Positioning & Moving products, please contact your Toffeln Key Account Manager today.

Product demonstrations & free trials are available on request.

DISCOVER THE POWER OF COMPRESSED AIR

“Our new ErgoTrics product line will ensure the patient’s safety and dramatically reduce the risk of injury to healthcare professionals when manoeuvring patients under anaesthetic.”

Danielle Leflaive
Sales Support



Er^{toffeln}goTrics

Do you suffer from back pain as a result of moving your patients?

With the use of compressed air instead of manual force, **ErgoTrics** products provide better support for patients without adding ergonomic burden and potential injury to healthcare professionals. This new technology ensures that patients are kept comfortable and correctly positioned during medical procedures and improves the quality of care provided.

Features

- Optimal patient positioning capabilities
- The use of compressed air means less chance of injury for the patient and the healthcare professional.
- Fully x-ray radiolucent
- Reusable depending on the infection risk



Optimal Positioning

The ErgoTrics cushions are designed as a mirror image of the human skeleton. In this way, they support the ribcage and pelvis, while leaving the abdomen free from pressure.



Better Ergonomics

Using compressed air instead of manual labour is an ergonomic breakthrough for hospital staff. 90% of existing users say that "minimal to no force is required to position patients into prone position".




X-Ray Transmission

The ErgoTrics solutions are fully X-ray radiolucent, which facilitates X-ray applications.



Optimal Hygiene

COVID-19 has heightened the world's awareness of attention to hygiene. ErgoTrics offers disposable, single patient use and reusable inflatable solutions depending on the infection risk.



Scan me to see ErgoTrics in action.

“At Toffeln, we are committed to providing innovative products that improve the well-being of healthcare staff and their patients. With the addition of ErgoTrics, we are confident that we can better serve the medical industry, improve patient outcomes and reduce the risk of injury to healthcare professionals.”

Harry Turner
Sales Director at Toffeln



PP-IBO
ErgoTrics Inflatable Board IBO

Box of 1
The IBO is reusable for 250 applications and is provided with a pressure sensor. The sensor counts automatically how often the IBO has been used and displays the usages with the LED lights. A handy feature to ensure safe use, also in case of regular shift changes.
Turning patients from supine to prone position remains a difficult task, both technically and ergonomically, especially with patients weighting over 80kg (175lb). A patient under general anaesthesia should not be tilted by one person alone. The Inflatable Board (IBO) is the ideal device for the safe and controlled tilting of patients.

- Features:**
- Makes process of proning (tilting) easier and faster
 - Better ergonomics for healthcare workers by using compressed air
 - Controlled and safe for patient (without torsion and specific pressure points)
 - Tilting with less people.

Inflate using Hummingbird Air Pistol (see page 103)

Topfel
ErgoTrics

Free demonstrations available on request




PP-IXSR / PP-IXSRT (with tape)
ErgoTrics Inflatable Support IXS Rectangular Cushion

Shape: Rectangle **Sizes:** L, XL **Qty:** 5

PP-IXSS / PP-IXSST (with tape)
ErgoTrics Inflatable Support IXS Sickle Cushion

Shape: Sickle **Sizes:** L, XL **Qty:** 5
Used to provides optimal access to the surgical or medical procedure site while maintaining body function and structural integrity.

- Features:**
- Adjustable (inflate/deflate)
 - 100% hygienic
 - X-ray permeable
 - Different shapes and sizes available



PP-IPS
ErgoTrics Inflatable Prone Support (IPS)

Size: L (up to 150kg) / XL (up to 180kg)
Box of 5

Used to avoid venous congestion and pressure on nerves, by mirroring the human skeleton and variations in size, the IPS supports the rib cage and pelvis whilst allowing the soft tissue and abdomen to hang freely. This reduces the risk of bedsores, Intra-Abdominal Pressure (IAP) and prevents venous congestion. The IPS also reduces the risk of injury to healthcare professionals.

- Benefits:**
- Safer positioning on air
 - Good hygiene
 - Better Ergonomics
 - Supports patient weight of up to 180kg



PP-OPEROLL
ErgoTrics Operoll

Size: L (up to 150kg) / XL (up to 180kg)
Box of 2

Used to provides optimal access to the surgical or medical procedure site while maintaining body function and structural integrity. The Operoll is placed under the patient on the operating table (or any other supporting surface). When the patient needs (re)positioning, the Operoll is partially or fully inflated using compressed air and will lift the patient into the desired position.

- Benefits:**
- Adjustable (inflate/deflate)
 - Better ergonomics for healthcare professional by using compressed air
 - 100% hygienic
 - X-ray permeable



“Better support
for patients
without the
ergonomic
burden on
healthcare
professionals”



The Hummingbird is compatible with all ErgoTrics products

PP-HUMMAP
ErgoTrics Hummingbird

Box of 1

A specially designed reusable air pistol for applications in the medical world. Ensures the safe use of medical compressed air in the hospital. The capacity of 4-6 bar can safely be used in a dosed manner to inflate the positioning devices from the ErgoTrics range.

Features:

- Compatible with all ErgoTrics devices
- Safe and easy to use
- Used to inflate and deflate the ErgoTrics positioning cushions and to adjust the pressure

CE
Topfelni
GelPositioners

Features

- Superior pressure relief
- Help prevent decubitus ulcers and nerve injury
- Ideal for sensitive skin
- Easy to clean
- High strength and elasticity
- Reusable
- Can be warmed or cooled
- Will not harden over time
- X Ray translucent
- Anti-bactericidal

For other shapes and forms, please contact your Key Account Manager



AP003
Closed Head Ring

Size: 20 x 7.5 x 4.5cm



AP013
Open Head Ring

Size: 20 x 7.5 x 4.5cm



AP002
Closed Head Ring

Size: 14 x 6 x 3.5cm



AP013/C
Open Head Ring With Cut

Size: 20 x 7.5 x 4.5cm



AP031
Heel Support Pad

Size: 20 x 11 x 7cm



AP071
Armboard Pad

Size: 40 x 12.5 x 2cm



AP072
Armboard Pad

Size: 60 x 12.5 x 2cm



AP043
Chest Roll (Adult)

Size: 50 x 15 x 10cm



AP022
Prone Head Rest

Size: 26 x 23 x 14.5cm



AP032
Heel Cups Without Strap

Size: 17 x 11 x 5cm



Toffeln[®]

Unit 2, Bridge Road, Kingswood,
Bristol BS15 4FW, United Kingdom

0345 500 4433 · sales@toffeln.com

toffeln.com