

MICRO-TOUCH® AFFINITY™

Performance Series

Tough neoprene exam gloves, delivering elastic comfort and robust tear resistance

- **Heightened defenses:** MICRO-TOUCH® AFFINITY™ medical exam gloves' tough neoprene formulation offers protection against snags and tears
- **Specialized protection:** These medical exam gloves are also tested against cytotoxic/chemotherapy drugs*, ensuring hand protection in high-risk environments
- **Enhanced grip:** These neoprene gloves' textured fingers grant them stronger grip performance in both wet and dry conditions
- **Elevated comfort:** The elastic design offers users excellent everyday comfort
- **Reduced allergy risk:** Neoprene material prevents Type 1 (latex) allergic reactions




*see product packaging for details

Recommended For

- Tested for use with chemotherapy drugs
- Protection from Type I latex allergy in HCW's or patients
- Enhanced grip in wet environments
- Sustained, high-risk examination procedures

FEATURED TECHNOLOGIES



MICRO-TOUCH® Everyday, Everywhere Protection

MICRO-TOUCH® premium medical examination gloves are scientifically designed to complement every need of medical professionals. Sophisticated design features include Ansell ergonomic technologies that ensure comfort during extended procedures, wet grip innovations that prevent harmful mishandling, and micro-thin films for maximum sensitivity. Caregivers who are fully prepared and comfortably protected by versatile MICRO-TOUCH styles have the peace of mind to focus on what matters most: the health and protection of themselves and their patients.

MICRO-TOUCH® AFFINITY™

Performance Series

PRODUCT INFORMATION

Material	Neoprene (Polychloroprene)	Not Made From Natural Rubber Latex	Yes
Color	Green	Cuff Length	Standard
Powder Content	Powder-Free	External Glove Surface	Textured Fingers
Freedom from Holes (Inspection level I)	1.5 AQL	Finger Thickness (mm/mil)	0.16 / 6.3
Palm Thickness (mm/mil)	0.13 / 5.1	Allergy Prevention	Latex (Type I)
Available Sizes	XS (5.5 - 6), S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10)	Tested For Use With Chemotherapy Drugs	Yes, in accordance with ASTM D6978 and US FDA cleared
Sterile	No	Antistatic	Not Tested
Glove Length (mm/inches)	240 / 9.5	Product Segmentation	Performance Series

ORDERING INFORMATION

Size	S (6.5 - 7)	M (7.5 - 8)	L (8.5 - 9)	XL (9.5 - 10)
Product Code	553352	553353	553354	553355

PACKAGING AND STORAGE

Packaging	100 gloves per dispenser; 10 dispensers/1000 gloves per case
Shelf Life	5 Years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.

STANDARDS AND CERTIFICATIONS

ASTM D6977, CA Prop 65, CE Number, Medical Device Regulation (EU) 2017/745 Class I, EN 455 Part 1, EN 455 Part 3, ISO 13485, ISO 10993-10:2002, UKCA

Vulcanization Chemical Accelerators

Zinc Dibutyldithiocarbamate (ZDBC)

Only a very small number of users may be sensitive to this ingredient(s) and hence may develop irritant and/or allergic contact reactions.

Contact your Ansell representative for ordering or more information.

North America
US Tel: 800 952 9916
CA Tel: 1-844-494-7854
Email: insidesalesus@ansell.com

Latin America & Caribbean
Tel: +52(442) 296 20 50
Email: cslac@ansell.com
Brazil Tel: +55-11-3356-3100
Brazil Email: luvas.medicas@ansell.com

Europe, Middle East & Africa
Tel: +32 (0) 2 528 74 00
Email: info@ansell.eu

Asia Pacific
Tel: +603 8310 6688
Email: apac.medical@ansell.com

Australia & New Zealand
Tel: +61 3 9270 7270
Email: protection@ap.ansell.com