


IMS™ Infinity Series™ Cassettes

Improve your Dental Practice



Infection Prevention and Instrument Management

Uniquely sophisticated, ergonomic design with revolutionary new hole pattern for optimal use in thermodisinfectors.

The IMS™ Infinity Series™ Cassettes are made of long-lasting, low-maintenance stainless steel. Durable, coloured silicone rails securely hold instruments throughout the processing cycle. Available in different sizes (1 DIN, 1/2 DIN, 1/4 DIN and container cassette), configurations (Ortho, Bone Spreading, Surgical, Diagnostics and Prophylaxis, etc.) and colour options, and can be specialised to meet your needs.

Ergonomic latch with positive locking mechanism is safe and easy-to-use, allows one-handed opening and visual confirmation if the cassette is locked /unlocked.

Revolutionary hole pattern design provides more access to instruments during cleaning and sterilization and increases compatibility with today's cleaning equipment, including automated washers.

CHOOSE THE SIZE.

Available in different sizes: 1 DIN, 1/2 DIN, 1/4 DIN and container cassette.



High quality, electropolished stainless steel protects against corrosion and provides more durability than plastic.

Innovative color-coded silicone rail system, available in 11 color options, significantly reduces instrument contact, allowing for more water flow while protecting the instruments during reprocessing.

CHOOSE THE COLOR.

Available in 11 colors for easy color coding by procedure or operator.

Red		0
Grey		1
White		2
Orange		3
Purple		4
Yellow		5
Lavendar		6
Black		7
Blue		8
Green		9
Ocean Blue		11



Smooth round corners and slotted edges increase drainage and reduce drying time.

Beveled edge allows for a more ergonomic grasp.



EFFICIENCY

IMS™ Cassettes eliminate time-consuming steps to streamline reprocessing, **which can save facilities an hour or more each day.**



ORGANIZATION

IMS™ Cassettes help facilities better organize their reprocessing flow, freeing up to seven linear feet of counter space, **so they can maximize the space in their sterilization area.**



SAFETY & COMPLIANCE

IMS™ Cassettes minimize how often contaminated instruments are handled, **improving compliance during transportation and reprocessing while also reducing the risk of a costly sharps injury.**





INSTRUMENT PROTECTION

IMS™ Cassettes are made of durable, stainless steel to protect the integrity of your instruments. **They also keep your instruments held in place and spaced out for an optimal environment for proper reprocessing.**






STANDARDIZATION

IMS™ Cassettes help facilities create a system for managing their procedural setups. **Each setup includes all of the instruments that are needed to help with onboarding and standardizing care** throughout a single facility and even across multiple locations.



Free Call: 00800 48 37 43 39 | Free Fax: 00800 48 37 43 40

E-Mail: info@hu-friedy.eu | Website: HuFriedyGroup.eu

   HuFriedyGroupEurope

Hu-Friedy Mfg. Co. LLC • European Headquarters • Lyoner Str. 9 • 60528 Frankfurt am Main, Germany • HuFriedyGroup.eu
All company and product names are trademarks of Hu-Friedy Mfg. Co. LLC, its affiliates or related companies, unless otherwise noted.
©2024 Hu-Friedy Mfg. Co., LLC All rights reserved. IMS-GB-IME/0224