

access

ENDO MOTOR

Wireless endodontic motor with integrated apex locator

- Facilitates root canal preparation
- Real-time feedback technology and dynamic torque control
- Automatic stop, slow down or reverse rotation when the file approaches or reaches the apical foramen
- Equipped with **27** pre-registered systems and **10** customisable programmes

Autostart mode

The motor starts when entering the canal

ATR mode

When the torque limit is reached the motor changes the rotation direction



ENDO MOTOR

Brushless motor !

- **360°** rotation of the contra-angle
- Precise continuous rotation and alternate movement
- Clear, easy-to-use **LED screen** and convenient **wireless** charging!

- **Speed:** 100 - 1800 rpm
- Torque range: 0,4 - 5 Ncm
- Weight of the handpiece with the contra-angle: **158g**
- Frequency range: 112 - 205 KHz

Motor
Garantie 2 ans

Motor
2 year warranty

Contra-angle
1 year warranty

Supplied with :

Instructions for use, quick guide, charging base, handpiece, contra-angle, measuring cable, lip hook, file holder, mains adaptor, protective silicone cap, spray nozzle, touch probe.

ENDO MOTOR

3166410

