

access

STERILE R

GO-TAPER FLEX

Continuous rotation



250 - 350 rpm


NiTi alloy and heat treatment

for controlled memory to prevent
ledge formation and perforations

Try GO-TAPER FLEX



Regain
FLEXIBILITY AND ENDURANCE !



Variable taper

STANDARD AND/OR
DIFFICULT CASES



PRE - OP




POST- OP



Clinical case by
Dr. Hugo Setbon
Endodontics specialist


GO-TAPER FLEX



Flexibility



Cutting power



Resistance




Single use

DURABLE FILES FOR RELIABLE TREATMENTS !

x6

FOR TREATMENT



GO-TAPER FLEX	LENGTH			
	19 mm	21 mm	25 mm	31 mm
ASSORTED KIT A0/A1/A2/B1/B2/B3	–	3148376	3148384	3148392
A0 – 4% – ISO 19	3148368	–	–	–
A1 – 2% – ISO 18	–	3148369	3148377	3148385
A2 – 4% – ISO 20	–	3148370	3148378	3148386
B1 – 7% – ISO 20	–	3148371	3148379	3148387
B2 – 8% – ISO 25	–	3148372	3148380	3148388
B3 – 9% – ISO 30	–	3148373	3148381	3148389
B4 – 6% – ISO 40	–	3148374	3148382	3148390
B5 – 5% – ISO 50	–	3148375	3148383	3148391

access

A COMPLETE RANGE OF STERILE FILES

	STANDARD CANALS STRAIGHT/WIDE	DIFFICULT CANALS NARROW/CURVED
Pre-enlargement		
	PREFILE	✓
	GO-GLIDER	✓
	EASY-FILE FLEX <i>Glider flexibility +++</i>	✓
Treatment		
	GO-TAPER FLEX <i>flexibility +++</i>	✓
	GO-TAPER BLUE <i>flexibility +++</i>	✓
	GO-TAPER X	✓
	GO-TWO	✓
	GO-TAPER UNIVERSAL	✓
	REVERSO BLUE <i>flexibility +++</i>	✓
	EASY-FILE FLEX <i>flexibility +++</i>	✓
Retreatment		
	GO-TAPER RETREATMENT	✓
	REVERSO SILVER	✓

access


STERILE R

PREFILE


FILES FOR CATHETERISATION



Try **PREFILE**



Offer yourself **MECHANISED
CATHETERISATION !**



STANDARD AND/OR
DIFFICULT CASES



-  **FLEXIBILITY**
-  **CYCLIC FATIGUE
RESISTANCE**
-  **FOAM GUIDE TIPS**

PROTOCOL with **PREFILE**




+ preparation




PREFILE


CATHETERISATION



Flexibility



Cutting power




Resistance



Single use

ENSURE THE PERMEABILITY OF THE ROOT CANAL BEFORE ROTARY INSTRUMENTATION !



x6

PREFILE	LENGTH		
	21 mm	25 mm	31 mm
ASSORTED KIT	3148396	3148400	3148404
2% – ISO 013 STE	3148393	3148397	3148401
2% – ISO 016 STE	3148394	3148398	3148402
2% – ISO 019 STE	3148395	3148399	3148403

access

A COMPLETE RANGE OF STERILE FILES

		STANDARD CANALS STRAIGHT/WIDE	DIFFICULT CANALS NARROW/CURVED
Pre-enlargement			
	PREFILE		✓
	GO-GLIDER		✓
	EASY-FILE FLEX <i>Glaser flexibility +++</i>		✓
Treatment			
	GO-TAPER FLEX <i>flexibility +++</i>	✓	✓
	GO-TAPER BLUE <i>flexibility +++</i>	✓	✓
	GO-TAPER X	✓	✓
	GO-TWO	✓	✓
	GO-TAPER UNIVERSAL	✓	✓
	REVERSO BLUE <i>flexibility +++</i>	✓	✓
	EASY-FILE FLEX <i>flexibility +++</i>	✓	✓
Retreatment			
	GO-TAPER RETREATMENT	✓	
	REVERSO SILVER	✓	

access

STERILE R


EASY-FILE FLEX *Glider*

Pre-enlarger




NiTi alloy and heat treatment

Try EASY-FILE FLEX Glider



Secure YOUR ENDODONTIC PROTOCOL !



Variable taper

STANDARD CASES

- FLEXIBILITY
- ENDURANCE
- FOAM GUIDE TIPS



POST - OP


Clinical case by
Dr. Hugo Setbon
Endodontics specialist

with **EASY-FILE FLEX Glider**
then EASY-FILE FLEX
for the preparation


PRE - OP




EASY-FILE
FLEX *Glider*



Flexibility



Cutting power




Resistance



Single use

TO ANTICIPATE A CONTROLLED ROOT CANAL PREPARATION !

PRE-ENLARGER



x6

EASY-FILE FLEX Glider	LENGTH		
	21 mm	25 mm	31 mm
2% – ISO 015	3148405	3148406	3148407

access

A COMPLETE RANGE OF STERILE FILES

		STANDARD CANALS STRAIGHT/WIDE	DIFFICULT CANALS NARROW/CURVED
Pre-enlargement			
	PREFILE		✓
	GO-GLIDER		✓
	EASY-FILE FLEX <i>Glider</i> flexibility +++		✓
Treatment			
	GO-TAPER FLEX flexibility +++	✓	✓
	GO-TAPER BLUE flexibility +++	✓	✓
	GO-TAPER X	✓	✓
	GO-TWO	✓	✓
	GO-TAPER UNIVERSAL	✓	✓
	REVERSO BLUE flexibility +++	✓	✓
	EASY-FILE FLEX flexibility +++	✓	✓
Retreatment			
	GO-TAPER RETREATMENT	✓	
	REVERSO SILVER	✓	

access

STERILE R

EASY-FILE FLEX

Alternate movement

Clockwise and counter-clockwise



250-350 rpm

NiTi alloy and heat treatment

for controlled memory to prevent
ledge formation and perforations

Try EASY-FILE FLEX



Simplify YOUR ENDODONTIC PROTOCOL
WITH A **SINGLE INSTRUMENT !**



Variable taper


STANDARD AND/OR
DIFFICULT CASES



PRE - OP



Master cone




POST - OP



Clinical case by
Dr. Hugo Setbon
Endodontics specialist


EASY-FILE FLEX



Flexibility



Cutting power




Resistance



Single use

ALL THE PREPARATION WITH 1 SINGLE INSTRUMENT !

FOR TREATMENT



EASY-FILE FLEX		LENGTH		
		21 mm	25 mm	31 mm
ASSORTED KIT	x4	3166615	3166620	3166625
S20 - 7% - ISO 20	x6	3166611	3166616	3166621
R25 - 7% - ISO 25	x6	3166612	3166617	3166622
M35 - 6% - ISO 35	x6	3166613	3166618	3166623
L45 - 5% - ISO 45	x6	3166614	3166619	3166624

access

A COMPLETE RANGE OF STERILE FILES

		STANDARD CANALS STRAIGHT/WIDE	DIFFICULT CANALS NARROW/CURVED
Pre-enlargement			
	CAT-FILE		✓
	GO-GLIDER		✓
	EASY-FILE <i>Glider flexibility +++</i>		✓
Treatment			
	GO-TAPER FLEXI <i>flexibility +++</i>	✓	✓
	GO-TAPER BLUE <i>flexibility +++</i>	✓	✓
	GO-TAPER X	✓	✓
	GO-TWO	✓	✓
	GO-TAPER UNIVERSAL	✓	✓
	REVERSO BLUE <i>flexibility +++</i>	✓	✓
	EASY-FILE FLEX <i>flexibility +++</i>	✓	✓
Retreatment			
	GO-TAPER RETREATMENT	✓	
	REVERSO SILVER	✓	

access

STERILE R

GO-GLIDER

Pre-enlarger




250 rpm


NiTi alloy

for optimal **fracture resistance**

Try GO-GLIDER



Regain FLEXIBILITY AND ENDURANCE !



Variable taper

STANDARD CASES



- FLEXIBILITY**
- ENDURANCE**
- FOAM GUIDE TIPS**


PROTOCOL with **GO-GLIDER**



+ preparation




GO-GLIDER



Flexibility



Cutting power




Resistance



Single use

TO ANTICIPATE A CONTROLLED ROOT CANAL PREPARATION !

PRE-ENLARGER



x6

GO-GLIDER	LENGTH		
	21 mm	25 mm	31 mm
2% – ISO 016	3166663	3166664	3166665

access

A COMPLETE RANGE OF STERILE FILES

	STANDARD CANALS STRAIGHT/WIDE	DIFFICULT CANALS NARROW/CURVED
Pre-enlargement		
	CAT-FILE	✓
	GO-GLIDER	✓
	EASY-FILE <i>Glider flexibility +++</i>	✓
Treatment		
	GO-TAPER FLEXI <i>flexibility +++</i>	✓
	GO-TAPER BLUE <i>flexibility +++</i>	✓
	GO-TAPER X	✓
	GO-TWO	✓
	GO-TAPER UNIVERSAL	✓
	REVERSO BLUE <i>flexibility +++</i>	✓
	EASY-FILE FLEX <i>flexibility +++</i>	✓
Retreatment		
	GO-TAPER RETREATMENT	✓
	REVERSO SILVER	✓

access

STERILE R

GO-TAPER UNIVERSAL


Continuous rotation



NiTi alloy

for optimal **fracture resistance**

Try GO-TAPER



Optimize **DEBRIS REMOVAL !**



Variable taper

STANDARD CASES



PRE - OP




POST - OP



Clinical case by
Dr. Hugo Setbon
Endodontics specialist


GO-TAPER UNIVERSAL



Flexibility



Cutting power



Resistance




Single use

SAME SEQUENCE FOR ALL STANDARD CANALS !

FOR TREATMENT

x6



GO-TAPER UNIVERSAL	LENGTH			
	19 mm	21 mm	25 mm	31 mm
ASSORTED KIT	–	3166578	3166586	3166594
A0 – 4% – ISO 19	3166570	–	–	–
A1 – 2% – ISO 18	–	3166571	3166579	3166587
A2 – 4% – ISO 20	–	3166572	3166580	3166588
B1 – 7% – ISO 20	–	3166573	3166581	3166589
B2 – 8% – ISO 25	–	3166574	3166582	3166590
B3 – 9% – ISO 30	–	3166575	3166583	3166591
B4 – 6% – ISO 40	–	3166576	3166584	3166592
B5 – 5% – ISO 50	–	3166577	3166585	3166593

access

A COMPLETE RANGE OF STERILE FILES

	STANDARD CANALS STRAIGHT/WIDE	DIFFICULT CANALS NARROW/CURVED
Pre-enlargement		
	CAT-FILE	✓
	GO-GLIDER	✓
	EASY-FILE <i>Glider flexibility +++</i>	✓
Treatment		
	GO-TAPER FLEXI <i>flexibility +++</i>	✓
	GO-TAPER BLUE <i>flexibility +++</i>	✓
	GO-TAPER X	✓
	GO-TWO	✓
	GO-TAPER UNIVERSAL	✓
	REVERSO BLUE <i>flexibility +++</i>	✓
	EASY-FILE FLEX <i>flexibility +++</i>	✓
Retreatment		
	GO-TAPER RETREATMENT	✓
	REVERSO SILVER	✓

access

STERILE R


GO-TAPER RETREATMENT

Continuous rotation



NiTi alloy

Try GO-TAPER



Optimize THE ELIMINATION OF A PREVIOUS
OBTURATION !



PRE - OP



Master cone




POST - OP



Clinical case by
Dr. Hugo Setbon
Endodontics specialist


GO-TAPER RETREATMENT



Flexibility



Cutting power



Resistance




Single use

FOR AN ELIMINATION OF THE ENDODONTIC FILLING !

FOR RETREATMENT

ADAPTED INSTRUMENTS to each third of the root canal



x6

GO-TAPER RETREATMENT	LENGTH	REFERENCE
ASSORTED KIT	Assorted	3166598
C1 – 9% – ISO 30	16 mm	3166595
C2 – 8% – ISO 25	18 mm	3166596
C3 – 7% – ISO 20	22 mm	3166597

access

A COMPLETE RANGE OF STERILE FILES

	STANDARD CANALS STRAIGHT/WIDE	DIFFICULT CANALS NARROW/CURVED
Pre-enlargement		
	CAT-FILE	✓
	GO-GLIDER	✓
	EASY-FILE <i>Glider flexibility +++</i>	✓
Treatment		
	GO-TAPER FLEXI <i>flexibility +++</i>	✓
	GO-TAPER BLUE <i>flexibility +++</i>	✓
	GO-TAPER X	✓
	GO-TWO	✓
	GO-TAPER UNIVERSAL	✓
	REVERSO BLUE <i>flexibility +++</i>	✓
	EASY-FILE FLEX <i>flexibility +++</i>	✓
Retreatment		
	GO-TAPER RETREATMENT	✓
	REVERSO SILVER	✓

access

STERILE R

GO-TAPER BLUE

Continuous rotation




250 - 350 rpm


NiTi alloy and heat treatment

for controlled memory to prevent
ledge formation and perforations

Try GO-TAPER BLUE



Regain FLEXIBILITY AND ENDURANCE !



Variable taper

STANDARD AND/OR
DIFFICULT CASES



PRE - OP




POST - OP



Clinical case by
Dr. Hugo Setbon
Endodontics specialist


GO-TAPER
BLUE



Flexibility



Cutting
power



Resistance




Single
use

NEW THERMAL TREATMENT "BLUE TREATMENT"

x6

FOR TREATMENT



GO-TAPER BLUE	LENGTH			
	19 mm	21 mm	25 mm	31 mm
ASSORTED KIT A0 à B3	–	3166634	3166642	3166650
A0 – 4% – ISO 019	3166626	–	–	–
A1 – 2% – ISO 018	–	3166627	3166635	3166643
A2 – 4% – ISO 020	–	3166628	3166636	3166644
B1 – 7% – ISO 020	–	3166629	3166637	3166645
B2 – 8% – ISO 025	–	3166630	3166638	3166646
B3 – 9% – ISO 030	–	3166631	3166639	3166647
B4 – 6% – ISO 040	–	3166632	3166640	3166648
B5 – 5% – ISO 050	–	3166633	3166641	3166649

access

A COMPLETE RANGE OF STERILE FILES

		STANDARD CANALS STRAIGHT/WIDE	DIFFICULT CANALS NARROW/CURVED
Pre-enlargement			
	CAT-FILE		✓
	GO-GLIDER		✓
	EASY-FILE <i>Glider flexibility +++</i>		✓
Treatment			
	GO-TAPER FLEXI <i>flexibility +++</i>	✓	✓
	GO-TAPER BLUE <i>flexibility +++</i>	✓	✓
	GO-TAPER X	✓	✓
	GO-TWO	✓	✓
	GO-TAPER UNIVERSAL	✓	
	REVERSO BLUE <i>flexibility +++</i>	✓	✓
	EASY-FILE FLEX <i>flexibility +++</i>	✓	✓
Retreatment			
	GO-TAPER RETREATMENT	✓	
	REVERSO SILVER	✓	

access

STERILE R


GO-TAPER X

Continuous rotation 
250-300 rpm


NiTi alloy

for controlled memory to prevent
ledge formation and perforations

Try GO-TAPER X



Optimize DEBRIS REMOVAL !



Variable taper

STANDARD CASES



CLINICAL CASE 1

PRE - OP



POST - OP



CLINICAL CASE 2

PRE - OP




POST - OP



Clinical case by
Dr. Hugo Setbon
Endodontics specialist


GO-TAPER X



Flexibility



Cutting power



Resistance




Single use

FELXIBILITY AND ENDURANCE !

Off-centre rectangular cross-section design

Asymmetrical rotary movement

FOR TREATMENT



x6

GO-TAPER X	LENGTH		
	21 mm	25 mm	31 mm
ASSORTED KIT	3166711	3166717	3166723
4% – ISO 017	3166706	3166712	3166718
6% – ISO 017	3166707	3166713	3166719
7% – ISO 030	3166708	3166714	3166720
6% – ISO 040	3166709	3166715	3166721
6% – ISO 050	3166710	3166716	3166722

access

A COMPLETE RANGE OF STERILE FILES

		STANDARD CANALS STRAIGHT/WIDE	DIFFICULT CANALS NARROW/CURVED
Pre-enlargement			
	CAT-FILE		✓
	GO-GLIDER		✓
	EASY-FILE <i>Glider flexibility +++</i>		✓
Treatment			
	GO-TAPER FLEXI <i>flexibility +++</i>	✓	✓
	GO-TAPER BLUE <i>flexibility +++</i>	✓	✓
	GO-TAPER X	✓	✓
	GO-TWO	✓	✓
	GO-TAPER UNIVERSAL	✓	
	REVERSO BLUE <i>flexibility +++</i>	✓	✓
	EASY-FILE FLEX <i>flexibility +++</i>	✓	✓
Retreatment			
	GO-TAPER RETREATMENT	✓	
	REVERSO SILVER	✓	

access

STERILE R

GO-TWO

Continuous rotation




NiTi alloy

for controlled memory to prevent ledge formation and perforations

Try GO-TWO



Optimize DEBRIS REMOVAL !




Constant taper


STANDARD AND/OR
DIFFICULT CASES



PROTOCOL



BASIC SEQUENCE



Basic sequence :

10/.04 - 15/.05 - 20/.06 - 25/.06

Apical enlargement :

30/.05 - 35/.04 - 40/.04

Retreatment :

15/.05 - 25/.06 - 25/.07


GO-TWO

CHARACTERISTICS


Length 21 mm : working section 16 mm

Length 25 mm : working section 16 mm et 21 mm


Length 31 mm : working section 21 mm



Flexibility



Cutting power



Resistance



Single use

x6

LENGTH OF THE WORKING SECTION : 16 mm					
TAPER	ISO	SPEED	TORQUE	21 mm	25 mm
Assorted	N°10-25	-	-	3166668	3166669
	N°25-40	-	-	3166670	3166671
4%	N°10	250	12.	3166673	3166674
	N°35	250	12.	3166675	3166676
	N°40	250	2	3166666	3166677
5%	N°15	250	12.	3166678	3166679
	N°30	250	12.	3166680	3166681
6%	N°20	250	2	3166682	3166683
	N°25	250	2	3166684	3166685
	N°30	250	1.5	3166667	3166686
	N°35	250	1.5	3166687	3166688
7%	N°40	250	2	-	3166689
	N°25	250	2	3166690	3166691

LENGTH OF THE WORKING SECTION : 21 mm					
TAPER	ISO	SPEED	TORQUE	25 mm	31 mm
Assorted	N°10-25	-	-	3166672	-
4%	N°10	250	12.	3166692	3166693
	N°35	250	12.	-	3166694
	N°40	250	2	-	3166695
5%	N°15	250	12.	3166696	3166697
	N°30	250	12.	-	3166698
6%	N°20	250	2	3166699	3166700
	N°25	250	2	3166701	3166702
	N°30	250	1.5	-	3166703
	N°35	250	1.5	-	3166704
	N°40	250	2	-	3166705

access

A COMPLETE RANGE OF STERILE FILES

		STANDARD CANALS STRAIGHT/WIDE	DIFFICULT CANALS NARROW/CURVED
Pre-enlargement			
	CAT-FILE		✓
	GO-GLIDER		✓
	EASY-FILE <i>Glider flexibility +++</i>		✓
Treatment			
	GO-TAPER FLEXI <i>flexibility +++</i>	✓	✓
	GO-TAPER BLUE <i>flexibility +++</i>	✓	✓
	GO-TAPER X	✓	✓
	GO-TWO	✓	✓
	GO-TAPER UNIVERSAL	✓	
	REVERSO BLUE <i>flexibility +++</i>	✓	✓
	EASY-FILE FLEX <i>flexibility +++</i>	✓	✓
Retreatment			
	GO-TAPER RETREATMENT	✓	
	REVERSO SILVER	✓	

access

STERILE R

REVERSO BLUE

Alternate movement


Clockwise and counter-clockwise




NiTi alloy and heat treatment

for controlled memory to prevent
ledge formation and perforations

Try REVERSO BLUE



Simplify YOUR ENDODONTIC PROTOCOL !




Variable taper

STANDARD AND/OR
DIFFICULT CASES




FLEXIBILITY
**CYCLIC FATIGUE
RESISTANCE**

PRE - OP




POST - OP



Clinical case by
Dr. Hugo Setbon
Endodontics specialist


REVERSO BLUE



Flexibility



Cutting power




Resistance



Single use

ALL THE PREPARATION WITH 1 SINGLE INSTRUMENT !

FOR TREATMENT



x6

REVERSO BLUE	LENGTH		
	21 mm	25 mm	31 mm
ASSORTED KIT	3166654	3166658	3166662
SO25 - 8% - ISO 25	3166651	3166655	3166659
SO40 - 6% - ISO 40	3166652	3166656	3166660
SO50 - 5% - ISO 50	3166653	3166657	3166661

access

A COMPLETE RANGE OF STERILE FILES

	STANDARD CANALS STRAIGHT/WIDE	DIFFICULT CANALS NARROW/CURVED
Pre-enlargement		
	CAT-FILE	✓
	GO-GLIDER	✓
	EASY-FILE <i>Glider flexibility +++</i>	✓
Treatment		
	GO-TAPER FLEXI <i>flexibility +++</i>	✓
	GO-TAPER BLUE <i>flexibility +++</i>	✓
	GO-TAPER X	✓
	GO-TWO	✓
	GO-TAPER UNIVERSAL	✓
	REVERSO BLUE <i>flexibility +++</i>	✓
	EASY-FILE FLEX <i>flexibility +++</i>	✓
Retreatment		
	GO-TAPER RETREATMENT	✓
	REVERSO SILVER	✓

access

STERILE R

REVERSO SILVER

Alternate movement

Clockwise and counter-clockwise



250-350 rpm

NiTi alloy

for controlled memory to prevent
ledge formation and perforations

Try **REVERSO SILVER**



Simplify YOUR **ENDODONTIC PROTOCOL !**




Variable taper

STANDARD CASES



PRE - OP



Disobturation



Master cone



POST - OP




Clinical case by

Dr. Hugo Setbon

Endodontics specialist


REVERSO SILVER



Flexibility



Cutting power




Resistance



Single use

ALL THE PREPARATION WITH 1 SINGLE INSTRUMENT !

FOR TREATMENT
AND RETREATMENT



x6

REVERSO SILVER	LENGTH		
	21 mm	25 mm	31 mm
ASSORTED KIT	3166602	3166606	3166610
SO25 - 8% - ISO 25	3166599	3166603	3166607
SO40 - 6% - ISO 40	3166600	3166604	3166608
SO50 - 5% - ISO 50	3166601	3166605	3166609

access

A COMPLETE RANGE OF STERILE FILES

	STANDARD CANALS STRAIGHT/WIDE	DIFFICULT CANALS NARROW/CURVED
Pre-enlargement		
	CAT-FILE	✓
	GO-GLIDER	✓
	EASY-FILE <i>Glider flexibility +++</i>	✓
Treatment		
	GO-TAPER FLEXI <i>flexibility +++</i>	✓
	GO-TAPER BLUE <i>flexibility +++</i>	✓
	GO-TAPER X	✓
	GO-TWO	✓
	GO-TAPER UNIVERSAL	✓
	REVERSO BLUE <i>flexibility +++</i>	✓
	EASY-FILE FLEX <i>flexibility +++</i>	✓
Retreatment		
	GO-TAPER RETREATMENT	✓
	REVERSO SILVER	✓