



| Designation | Sterilisation pouches Medibase | | | | | | | | | |
|--|---|---|------------|------------|------------|-------------|------------|----------------------------------|----------|------------|
| Reference(s) | SELF SEALING POUCHES | | | | | | | | | |
| | Size | 8,5 x 16,5cm | 7 x 25,5cm | 9 x 25,5cm | 14 x 28cm | 25,5 x 38cm | 5,7 x 13cm | 18 x 33cm | 9 x 13cm | 30 x 45 cm |
| | Ref | 3135469 | 3135470 | 3135471 | 3135472 | 3135473 | 3135474 | 3135475 | 3135476 | 3135477 |
| | HEAT SEALING POUCHES | | | | | | | | | |
| | Size | 7 x 25,5cm | | | 9 x 25,5cm | | | 20 x 33cm AVEC SOUFFLET DE 5,5cm | | |
| Ref | 3135479 | | | 3135480 | | | 3135478 | | | |
| Product picture |  | | | | | | | | | |
| Product description | Sterilization bags recommended for autoclaving or chemiclave. | | | | | | | | | |
| Characteristics and composition | Paper | The basis paper weight is $\pm 60\text{g/m}^2$, fitting the requirements of sterile packaging manufactured paper which list on ISO11607-1, 2 and EN868-5; According to the pouch designs, the printed inscription on the rolls is navy blue (Pantone298U), printing steam and EO indicator. | | | | | | | | |
| | Film | Lucid wathet blue hyperthermal PET/CPP complex film, thickness is $\pm 0,052\text{mm}$, weight is 56 g/m^2 , conforming to the related requirements of EN868-5. | | | | | | | | |
| | Ink | Using ink; the min printing area of indicator should be 100mm^2 , may not fade or become stained. | | | | | | | | |

| | |
|---------------------------|--|
| Packaging | Box of 200 pouches, available in 10 sizes |
| Shelf life | 5 years from manufacturing date |
| Indications of use |  |
| Legal information | <p>Medical Device of class I CE Marking Manufacturer: FOMED Industries Inc</p> |

Version 2.3

Updated 15/11/2022 by RG