


ImpressPLUS Clear

Transparent Vinyl Polysiloxane Material



ImpressPLUS Clear is a transparent medium viscosity vinyl polysiloxane (VPS) bite registration material.

Ideal for where a clear view of the area will optimise the accuracy of the impression. Other indications include as a silicone guide for composite restorations, needle point registrations, standard bite registrations in the end bite position and production of small model segments.



Delivery System

- 2 x 50ml cartridge & tips

Presentation

REF 010/0348

Characteristics

- Outstanding transparency
- Excellent dimensional stability
- Very high tear strength
- Non-sticky

Technical Data

Shore A hardness	70±2
Linear dimensional change (%)	<0.3
Tear strength (N/mm)	>7.0
Colour	Colorless transparent
Strain in compression (%)	3-5
Mixing ratio	1:1
Mixing time (min)	Auto-mix system
Total working time (min)	1:00
Time in mouth (min)	2:00
Total setting time (min)	3:00